

## Goat Anti-rabbit IgG(H&L), Alexa Fluor 405 conjugated

Cat. Number: bs-0295G-AF405

Quantity size: 100ul, 1ml

Concentration: 1mg/ml. Buffer = 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50%

Glycerol.

**Background:** Immunoglobulin G (IgG), is one of the most abundant proteins in serum

with normal levels between 8-17 mg/mL in adult blood. lgG is important for our defence against microorganisms and the molecules are produced by B lymphocytes as a part of our adaptive immune response. The lgG molecule has two separate functions; to bind to the pathogen that elicited

the response and to recruit other cells and molecules to destroy the antigen. The variability of the IgG pool is generated by somatic

recombination and the number of specificities in an individual at a given

time point is estimated to be 1011 variants.

**Specificity:** Goat Polyclonal IgG, affinity purified by Protein G, nonspecific adsorbed.

· Reacts with: Rabbit

· Immunogen: Native rabbit IgG

**Application:** · Immunofluorescence: 1:100-1000

· Flow-Cyt : 1ug/Test

Excitation spectrum: 405nmEmission spectrum: 424nm

· Optimal working dilutions must be determined by the end user.

· Protect from light.

Storage: Shipped at 4°C, Store at -20°C (Avoid repeated freeze/thaw cycles).

**Important Note:** This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.