



## Goat Anti-Rabbit IgG, Co-Gold conjugated

Cat. Number: bs-0	)295G-Gold
-------------------	------------

**Quantity size:** 0.5 ml / 1 ml  $\square$  5 nm  $\square$  10 nm  $\square$  15 nm  $\square$  35 nm

**Concentration:** 0.4mg/ml. Buffer = 20mM TBS(pH8.0) with 1% BSA, 0.03% Proclin300.

**Background:** Immunoglobulin G (IgG), is one of the most abundant proteins in serum

with normal levels between 8-17 mg/mL in adult blood. IgG is important for our defence against microorganisms and the molecules are produced by B lymphocytes as a part of our adaptive immune response. The IgG molecule has two separate functions; to bind to the pathogen that elicited the response and to recruit other cells and molecules to destroy the antigen. The variability of the IgG pool is generated by somatic

recombination and the number of specificities in an individual at a given

time point is estimated to be 1011 variants.

**Specificity:** • Goat Polyclonal IgG, affinity purified by Protein G.

Reacts with: Rabbit

· Immunogen: Rabbit IgG

**Application:** • IEM=1:20-200 IGS=1:20-200 GICA=1:20-200

· Optimal working dilutions must be determined by the end user.

Storage: Store at 2-8 °C for 3-6 months. Avoid repeated freeze/thaw cycles...

**Important Note:** This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.