

**&** 400-901-9800 sales@bioss.com.cn support@bioss.com.cn

## Streptavidin, Cy3 conjugated

**Catalog Number :** bs-0437P-Cy3

Quantity size: 100ul / 1ml

Concentration: 1mg/ml. Buffer = 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

**Background:** Streptavidin is a tetrameric protein composed of identic subunits. Each subunit binds one biotin molecule with a KD of ~1x10-15M. The preparation contains an N- and C-terminal shortened variant (core streptavidin) with improved properties concerning homogeneity, solubility, resistance towards proteolytic degradation and accessibility of the biotin binding pocket as compared to native streptavidin. The high affinity recognition of biotin and biotinylated molecules has made streptavidin one of the most important components in diagnostics and laboratory kits.

Specificity:	• mol wt: 66kDa
	<ul> <li>Purity: ≥98%</li> </ul>
	Cy3 conjugated.
Application:	Excitation spectrum: 552nm
	Emission spectrum: 570nm
	<ul> <li>Optimal working dilutions multiple</li> </ul>

Optimal working dilutions must be determined by the end user.

552nm

Storage: Shipped at 4°C, Store at -20°C (Avoid repeated freeze/thaw cycles).

Important Note: This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.