

Mouse IgG1, Isotype Control

Cat. Number:	bs-30296P
Quantity Size:	100ul, 500ul
Purity:	≥95%, Purified by Protein G
Form:	Liquid
Concentration:	1ug/ul
Buffer:	0.01M PBS, 0.5%BSA, 0.03%Proclin300.
Storage:	Store at 2-8°C.
Application:	Isotype Control
	0.125ug/1X10 ⁶ cells

Isotype controls are negative controls that are designed to measure the amount of non-specific background signal caused by the primary antibody. Isotype Control have no specificity for the target protein, yet retain all of the non-specific characteristics of an antibody.