Recombinant Human ACE2 (Fc Tag)

Cat. No. bs-46002P



Description	
Protein Sequence	Recombinant human ACE2 with a Fc tag and Avi at the C-terminal (Gln18-Ser740).
Source	Mammalian Expression System
Accession	Q9BYF1
Mol wt	The protein has a predicted MW of 111.2kDa. Due to glycosylation, the protein migrates to 115-125KD based on Bis-Tris PAGE result.
Endotoxin	Less than 0.5EU per ug by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE
	>95%as determined by HPLC
Activity assay	Not tested.
Formulation an	d Storago

Formulation and Storage


Storage	The product should be stored at -70°C or -20°C.
Formulation	trehalose is added as protectant before lyophilization.)
Formulation	Lyophilized powder (Lyophilized from 0.22um filtered solution in 20mM PB (pH 7.4). Normally 5%

Background

ACE2 (Angiotensin I Converting Enzyme 2) is a Protein Coding gene. Diseases associated with ACE2 include Severe Acute Respiratory Syndrome and Neurogenic Hypertension. The protein encoded by this gene belongs to the angiotensin-converting enzyme family of dipeptidyl carboxydipeptidases and has considerable homology to human angiotensin 1 converting enzyme. This secreted protein catalyzes the cleavage of angiotensin I into angiotensin 1-9, and angiotensin II into the vasodilator angiotensin 1-7.

Assay Data

Tris-Bis PAGE




Recombinant human ACE2 on Tris-Bis PAGE under reduced condition. The purity is greater than 95%.

HPLC Data MAU 1006040200

The purity of human ACE2 is greater than 95% as determined by SEC-HPLC.

ELISA Data

Biotinylated nCov-S1, His Tag ELISA



Immobilized biotinylated n-COV S1, His Tag at 1ug/ml (100ul/Well). Dose response curve for ACE2, Fc Tag with the EC50 of 81ng/ml determined by ELISA.