Recombinant 2019 nCOV Spike RBD (His Tag)

bs-46003P




Description	
Protein Sequence	2019-nCoV S protein RBD is expressed with a His tag and Avi at the C-terminal (Arg319-Asn532).
Source	Mammalian Expression System
Accession	QHD43416.1
Mol wt	The protein has a predicted MW of 27 kDa. Due to glycosylation, the protein migrates to 36-40KDa based on Bis-Tris PAGE result.
Endotoxin	Less than 1EU per ug by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE
	> 95% as determined by HPLC
Activity assay	Not tested.
Formulation and Storage	
Formulation	Lyophilized powder (Lyophilized from 0.22um filtered solution in 20mM PB (pH 7.4). Normally 5% trehalose is added as protectant before lyophilization.)
Storage	The product should be stored at -70 $^{\circ}\mathrm{C}$ or -20 $^{\circ}\mathrm{C}$
Background	
	The spike protein (S) of coronavirus (CoV) attaches the virus to its cellular receptor, angiotensin- converting enzyme 2 (ACE2). A defined receptor-binding domain (RBD) on S mediates this interaction. The S protein plays key parts in the induction of neutralizing-antibody and T-cell

Assay Data


Cat. No.





responses, as well as protective immunity.

Recombinant 2019-nCoV S protein RBD on Tris-Bis PAGE under reduced condition. The purity is greater than 95%.



The purity of 2019-nCoV S protein RBD is greater than 95% as determined by SEC-HPLC