Biotinylated Recombinant SARS Spike RBD (His&Avi

Tag)



Cat. No. bs-46008P-Biotin

Description	
Protein Sequence	Biotinylated SARS Spike RBDwith a His tag and Avi at the C-terminal (Arg306-Phe527).
Source	Mammalian Expression System
Accession	P59594
Mol wt	The protein has a predicted MW of 27.9 kDa. Due to glycosylation, the protein migrates to 36-46KDa based on Bis-Tris PAGE result.
Endotoxin	Less than 1EU per ug by the LAL method.
Purity	> 95% as determined by Bis-Tris PAGE
	>95% as determined by HPLC
Activity assay	Not tested.
Formulation an	d Storage
Formulation	Lyophilized powder (Lyophilized from 0.22um filtered solution in 20mM PB (pH 7.4). Normally 5% trehalose is added as protectant before lyophilization)
Storage	The product should be stored at -70 $^{\circ}$ C or -20 $^{\circ}$ C.
Background	
	The spike protein (S) of coronavirus (CoV) attaches the virus to its cellular receptor, angiotensin- converting enzyme 2 (ACE2). A defined receptor-binding domain (RBD) on S mediates this interaction. The S protein plays key parts in the induction of neutralizing-antibody and T-cell responses, as well as protective immunity.
Assess Data	