## **Recombinant SARS Nucleocapsid Protein (His tag)**


Cat. No. bs-49002P



Description	
<b>Protein Sequence</b>	SARS Nucleocapsid Protein full length with a His tag in N terminus (1-422).
Source	Escherichia coil Expression System
Accession	
Mol wt	46kD
Endotoxin	Not tested.
Purity	≥90% as determined by SDS-PAGE
Application	Recommended for sandwich immunoassays in ELISA and CLIA. Each laboratory should determine an optimum working titer for use in its particular application.
Activity assay	Not tested.
Formulation an	d Storage
Format	Liquid
Concentration	≥0.5 mg/ml
Buffer	sterile PBS, pH7.4
Storage	Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles.
Background	
	The nucleocapsid (N) protein of SARS-coronavirus (SARS-CoV) is the key protein for the formation of

the helical nucleocapsid during virion assembly. The nucleocapsid (N) protein of SARS-CoV enters the host cell together with the viral RNA and interferes with several cellular processes. Some of these processes involve interactions between SARS-CoV N protein and host cell proteins. It has also been demonstrated that the SARS-CoV N protein can bind to DNA in vitro. These interactions might have a role in the pathology of SARS. The N protein may be of potential value in vaccine development for specific prophylaxis and treatment against SARS.





SDS-PAGE for recombinant SARS Nucleocapsid Protein