

## **Product Datasheet**

Name: Mouse Anti-SARS N (nucleocapsid) Monoclonal Antibody

**Description:** Hybridoma clones have been derived from hybridization of myeloma cells with

spleen cells of BALB/c mouse immunized with recombinant SARS N (nucleocapsid).

Catalog No.	Isotype	Clone No.	Usage	Buffer
bsm-49131M	lgG2b	8G8A	Capture /Detection	10mM PBS (pH7.4)

**Specificity:** Mab react with recombinant antigen SARS N (nucleocapsid)

**Host:** Mouse

Clonality: Monoclonal

Format: Liquid

**Concentration:** ≥1 mg/ml

**Purification:** ≥90% (SDS-PAGE)

Preservative: 0.1% Proclin300

**Application:** Recommended for sandwich immunoassays in ELISA and CLIA. Each laboratory

should determine an optimum working titer for use in its particular application.

**Storage:** Store at -20 °C for three years. Avoid repeated freeze/thaw cycles.

**Background:** The nucleocapsid (N) protein of SARS-coronavirus (SARS-CoV) is the key protein for the formation of the helical nucleocapsid during virion assembly. The nucleocapsid (N) protein of SARS-CoV enters the host cell together with the viral RNA and interferes with several cellular processes. Some of these processes involve interactions between SARS-CoV N protein and host-cell proteins. It has also been demonstrated that the SARS-CoV N protein can bind to DNA in vitro. These interactions might have a role in the pathology of SARS. The N protein may be of potential value in vaccine development for specific prophylaxis and treatment against SARS.

**Note:** This product as supplied is intended for research or further manufacturing use only.