

\$400-901-9800
sales@bioss.com.cn
support@bioss.com.cn

Recombinant human OSM (Active, Tag free)

Cat. Number:	ck0012
Quantity size:	10µg / 50ug / 500µg
Protein Sequence:	26-252aa, Tag free, full length mature protein
Swiss-Prot:	P13725
Gene ID:	5008
Source:	Human cells derived
Structure:	Glycosylated monomer
Purity:	>95% by SDS-PAGE
MW:	30kDa
Endotoxin Level:	<0.5EU/ug
Formulation:	Lyophilized from a 0.2µm filtered solution in PBS without carrier protein
	Animal and Xeno free
Activity Assay:	The activity was measured by its ability to stimulate the proliferation of human TF-1
	cells (human erythroleukemic indicator cell line).
Reconstitution:	Briefly centrifuge the vial before opening. It is recommended to reconstitute the
	protein in sterile PBS containing 0.1% endotoxin-free recombinant human serum
	albumin.
Stability & Storage:	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. In general: 12
	months from date of receipt, -20 to -80°C as supplied. 1 month, 2 to 8°C under
	sterile conditions after reconstitution. 3 months, -20 to -80 $^\circ\!C$ under sterile
	conditions after reconstitution.
Description:	Oncostatin M (OSM) is a 30 kDa secreted cytokine belonging to the Interleukin 6
	family. OSM inhibits the proliferation of a number of tumor cell lines. It also
	stimulates proliferation of AIDS-KS cells. It regulates cytokine production, including
	IL-6, G-CSF and GM-CSF from endothelial cells. OSM plays crucial roles in
	inflammation, neuroprotection, hematopoiesis, metabolism and development.

Important Note: This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.