

**线虫靶标蛋白抗体**

产品编号	产品名称	应用
bs-8770R	C-actin	WB
bs-8771R	C-gapdh	WB
bs-8772R	C-tubulin-alpha	WB
bs-8773R	C-tubulin-beta	WB
bs-8774R	Histone H3	WB
bsm-8775M	Histone H3	WB

**拟南芥靶标蛋白抗体**

产品编号	产品名称	应用
bs-6988R	RuBisCO large subunit	WB,ELISA,IHC-F,IF
bs-2219R	Pectinesterase inhibitor 18	WB,ELISA,IHC-P,IHC-F,IF
bs-2222R	cellulase	WB,ELISA,IHC-P,IHC-F,IF
bs-2221R	ATEXPA10	WB,ELISA,IHC-P,IHC-F,IF
bs-2247R	AHA3	ELISA,IHC-P,IHC-F,IF
bs-4516R	CPN10	ELISA,IHC-P,IHC-F
bs-0245R	ELP1	IHC-P,IHC-F,IF
bs-13098R	eRF1	ELISA,IHC-P,IHC-F,IF,ICC
bs-0401R	Metal ion transporter	WB,ELISA,IHC-P,IHC-F,IF
bs-2249R	ACA11	ELISA,IHC-P,IHC-F,IF
bs-4517R	LOC103437570	WB,ELISA,IHC-P,IHC-F
bs-4518R	LOC103437570	ELISA,IHC-P,IHC-F,IF