

2024 Jan-Bioss 产品文献引用精选

1. Complement in breast milk modifies offspring gut microbiota to promote infant health

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-2673R	Rabbit Anti-Human C5b-9 antibody	38242132	64.5	FC	Staphylococcus lentus

2. Targeting Erbin-mitochondria axis in platelets/megakaryocytes promotes B cell-mediated antitumor immunity

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-5838R	Rabbit anti-ESAM antibody	38232736	29	IF,CoIP	Mouse,Human

3. EGFR-driven lung adenocarcinomas coopt alveolar macrophage metabolism and function to support EGFR signaling and growth

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1231R-FITC	Rabbit Anti- α ABCG1/FITC Conjugated antibody	38241033	28.2	FC	Mouse

4. Pro-Regenerative Glycopeptide Hydrogel Activates Type 2 Immune Response for Wound Healing via Macrophage-T Cell Crosstalk

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1452R	rabbit polyclonal anti-GATA3 antibody	10.1002/adfm.202307711	19	IF	Mouse

5. The therapeutic effect of gold clusters in an EAE model of multiple sclerosis by regulating T cell differentiation

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1317R	anti-stat1 antibody	10.1016/j.nantod.2023.102128	17.4	WB	Mouse
bs-0296G-FITC	goat anti-mouse IgG H&L/FITC	10.1016/j.nantod.2023.102128	17.4	IHC	Mouse

6. Subcellular Nanobionic Liposome with High Zeta Potential Enhances Intravesical Adhesion and Drug Delivery

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
-----	-------------	-----------	--------------	--------------	-----------------

bs-1035R	rabbit polyclonal anti-CD86 antibody	38252681	17.1	IF	Mouse
bs-2527R	rabbit polyclonal anti-CD163 antibody	38252681	17.1	IF	Mouse

7.A Dual Role of Mesenchymal Stem Cell Derived Small Extracellular Vesicles on TRPC6 Protein and Mitochondria to Promote Diabetic Wound Healing

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-10423R	anti-Collagen I antibody	38290527	17.1	WB	Rat

8.Reprogramming mechanism dissection and trophoblast replacement application in monkey somatic cell nuclear transfer

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0335Rs	rabbit anti-monkey IgG whole serum	38228612	16.6	Other	Monkey

9.A drug-free cardiovascular stent functionalized with tailored collagen supports in-situ healing of vascular tissues

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1035R	rabbit polyclonal anti-CD86 antibody	38272886	16.6	IF	Mouse

10.LNP-RNA-engineered adipose stem cells for accelerated diabetic wound healing

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-3337R	anti-Phospho-PKR (Thr446+Thr451)antibody	38272900	16.6	FC	Mouse

11.Rescuing SERCA2 pump deficiency improves bone mechano-responsiveness in type 2 diabetes by shaping osteocyte calcium dynamics

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-25532R	anti-Osx antibody	38291059	16.6	WB	Mouse

12.Serine synthesis sustains macrophage IL-1 β production via NAD⁺-dependent protein acetylation

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-6741R	Rabbit polyclonal anti-ASC antibody	38266638	16	IF	Mouse

13. Intermittent cytomegalovirus infection alters neurobiological metabolism and induces cognitive deficits in mice

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-2527R-AF680	AF680 anti-CD163 antibody	38182037	15.1	FC	Mouse

14. Lubricating MXenzyme-based hybrid hydrogel reverses oxidative damage to alleviate osteoarthritis

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-10725R	anti-Sox9 antibody	10.1016/j.cej.2024.148815	15.1	WB	Rat

15. Inorganic Phosphate as "Bioenergetic Messenger" Triggers M2-Type Macrophage Polarization

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-19804R	anti-SLC20a1 antibody	38247159	15.1	IF	Rat

16. Anti-inflammasome bio-heterojunction (AI-bioHJ): Revolutionizing diabetic wound healing with in situ self-transformation and programmed gas therapy

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-10021R	anti-NLRP3 antibody	10.1016/j.cej.2024.149014	15.1	WB,IHC	Rat

17. Expanding the insight of ecological risk on the novel chiral pesticide mefentrifluconazole: Mechanism of enantioselective toxicity to earthworms (*Eisenia fetida*)

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0523R	anti-Sox2 antibody	38271877	13.6	WB	Earthworm
bs-3578R	anti-TCTP antibody	38271877	13.6	WB	Earthworm

18. The role of mTORC1/TFEB axis mediated lysosomal biogenesis and autophagy impairment in fluoride neurotoxicity and the intervention effects of resveratrol

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-3501R	anti-synapsin-1 (Syn1) antibody	38335616	13.6	WB	Rat

19. Relieving Macrophage Dysfunction by Inhibiting SREBP2 Activity: A Hypoxic Mesenchymal Stem Cells-Derived Exosomes Loaded Multifunctional Hydrogel for Accelerated Diabetic Wound Healing

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-2072R	anti-mouse/human/rat/chicken iNOS antibody	38247194	13.3	FC	Mouse

20.CYB5R3 functions as a tumor suppressor by inducing ER stress-mediated apoptosis in lung cancer cells via the PERK-ATF4 and IRE1 α -JNK pathways

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-12162R	anti-CYB5R3 antibody	38253797	12.8	WB,IHC,IF	Human,Mouse

21.Pituitary neuroendocrine tumors with PIT1/SF1 co-expression show distinct clinicopathological and molecular features

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1912R	anti-alpha-subunit antibody	38228887	12.7	IHC	Human

22.Association between NTCP hepatic expression and inflammation/fibrosis as well as gender-specific differences in chronic HBV-infected patients

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1958R	anti-NTCP antibody	38258306	12.7	WB	Human

23.Keratinocyte-to-macrophage communication exacerbate psoriasiform dermatitis via LRG1-enriched extracellular vesicles

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0061R	anti- β -actin antibody	38250043	12.4	WB	Human,Mouse
bs-0638R	anti-TGF β R1 antibody	38250043	12.4	IF	Mouse
bs-0295G-HRP	HRP goat anti-rabbit antibody IgG	38250043	12.4	WB	Human,Mouse
bs-0296G-HRP	HRP goat anti-mouse antibody IgG	38250043	12.4	WB	Human,Mouse

24.Integrative analysis of spatial and single-cell transcriptome data from human pancreatic cancer reveals an intermediate cancer cell population associated with poor prognosis

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-11163R	anti-SPRR3 antibody	38297291	12.3	IHC	Human

25. Gut microbiome-derived butyrate inhibits the immunosuppressive factors PD-L1 and IL-10 in tumor-associated macrophages in gastric cancer

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-20373R	anti-IL-10 antibody	38197259	12.2	IHC	Mouse

26. The zinc finger protein DHHC09 S-acylates the kinase STRK1 to regulate H2O2 homeostasis and promote salt tolerance in rice

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0639R	anti-GFP antibody	38180963	11.6	WB	Rice
bs-2247R	anti-AHA3 antibody	38180963	11.6	WB	Rice

27. Targeting ZDHHC21/FASN axis for the treatment of diffuse large B-cell lymphoma

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
AK230	Free Aliphatic Acid Assay Kit	38195819	11.4	Other	Other

28. Insights into the Acceleration Mechanism of Intracellular N and Fe Co-doped Carbon Dots on Anaerobic Denitrification Using Proteomics and Metabolomics Techniques

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
AK300	NAD ⁺ /NADH assay kit	38268063	11.4	Other	Other

29. Nano-modified viruses prime the tumor microenvironment and promote the photodynamic virotherapy in liver cancer

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0648R	rabbit anti-CD8 antibody	38163894	11	IHC	Mouse
bsm-52535R	rabbit anti-MHC II antibody	38163894	11	IHC	Mouse
bsm-60761R	rabbit anti-CD206 antibody	38163894	11	IHC	Mouse

30. Chemo-photothermal nanoplatform with diselenide as the key for ferroptosis in colorectal cancer

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bsm-60738R	anti-Ki67 antibody	38224739	10.8	IHC	Mouse

31. Risedronate-functionalized manganese-hydroxyapatite amorphous particles: A potent adjuvant for subunit vaccines and cancer immunotherapy

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0283P-FITC	OVA-FITC	38244843	10.8	Other	Other

32. Chronic arsenic exposure-provoked biotoxicity involved in liver-microbiota-gut axis disruption in chickens based on multi-omics technologies

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0189R	anti- α -SMA antibody	38237767	10.7	WB	Chicken
bs-0718R	anti-Smad2 antibody	38237767	10.7	WB	Chicken
bsm-33207M	anti-TNF- α antibody	38237767	10.7	WB	Chicken
bs-0465R	anti-NF- κ B antibody	38237767	10.7	WB	Chicken
bs-23709R	anti-Lamin B antibody	38237767	10.7	WB	Chicken

33. Double-layered N-S1 protein nanoparticle immunization elicits robust cellular immune and broad antibody responses against SARS-CoV-2

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-4819R	Rabbit Anti-CD45 antibody	38291444	10.2	IF	Mouse

34. Novel-Ink-Based Direct Ink Writing of Ti6Al4v Scaffolds with Sub-300 μ m Structural Pores for Superior Cell Proliferation And Differentiation

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
BA00208	CCK-8	38180708	10	Other	Other
C01-06002	4% paraformaldehyde	38180708	10	Other	Other

35. A Versatile Intestine-on-Chip System for Deciphering the Immunopathogenesis of Inflammatory Bowel Disease

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1086R	rabbit polyclonal anti-human laminin α 5 antibody	38253407	10	IF	Human

36. GDF15 activates AMPK and inhibits gluconeogenesis and fibrosis in the liver by attenuating the TGF- β 1/SMAD3 pathway

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-13386R	anti-G6Pase antibody	38176644	9.8	WB	Mouse

37.Zinc finger protein ZNF638 regulates triglyceride metabolism via ANGPTL8 in an estrogen dependent manner

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bsm-49019M	anti-Angptl8 antibody	38211696	9.8	ELISA	Mouse

38.Manganese disrupts the maturation and degradation of axonal autophagosome leading to hippocampal synaptic toxicity in mice via the activation of LRRK2 on phosphorylation of Rab10

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bsm-51545M	anti-LRRK2 antibody	38224893	9.8	WB	Mouse

39.Atrazine exposure promotes cardiomyocyte pyroptosis to exacerbate cardiotoxicity by activating NF- κ B pathway

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-6319R	Rabbit anti-IL-1 β antibody	38224882	9.8	IHC	Mouse
bs-6741R	Rabbit anti-ASC antibody	38224882	9.8	IHC	Mouse

40.Salmonella manipulates the host to drive pathogenicity via induction of interleukin 1 β

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-6982R	Anti-Neutrophil Elastase antibody	38236896	9.8	IHC	Mouse

41.Photo-induced universal modification of small-diameter decellularized blood vessels with a hemocompatible peptide improves in vivo patency

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0646R-FITC	FITC-conjugated Anti-CD34 antibody	10.1016/j.actbio.2024.01.012	9.7	IF	Rat
bs-0189R-PE	PE anti- α -SMA antibody	10.1016/j.actbio.2024.01.012	9.7	IF	Rat
bs-4754R-PE	PE anti-vWF antibody	10.1016/j.actbio.2024.01.012	9.7	IF	Rat

42.ROS-Suppression Nanoplatform Combined Activation of STAT3/Bcl-2 Pathway for Preventing Myocardial Infarction in Mice

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bsk12015	Mouse IL-1 β ELISA kit	38288981	9.5	ELISA	Mouse
bsk13005	Mouse IL-6 ELISA kit	38288981	9.5	ELISA	Mouse

bsk12002	Mouse TNF- α ELISA kit	38288981	9.5	ELISA	Mouse
----------	-------------------------------	----------	-----	-------	-------

43.MMP14high macrophages orchestrate progressive pulmonary fibrosis in SR-Ag-induced hypersensitivity pneumonitis

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bsm-33187M	mouse anti-ACTA2 antibody	38218353	9.3	IF,WB	Mouse
bs-0126R	rabbit anti-HSP70 antibody	38218353	9.3	WB	Mouse

44.Transient immune activation without loss of intraepidermal innervation and associated Schwann cells in patients with complex regional pain syndrome

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-15197R	rabbit anti-C5a(Complement component 5a) antibody	38233858	9.3	IF	Human

45.Dioscin Alleviates Periodontitis by Inhibiting NLRP3 Inflammasome Activation via Regulation of K+ Homeostasis and Mitochondrial Function

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1278R	and anti-8-hydroxy-2'-deoxyguanosine (8-OHdG) antibody	38385066	9.2	IF,IHC	Mouse

46.TRIM22 induces cellular senescence by targeting PHLPP2 in hepatocellular carcinoma

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-3233R	Anti-P-IKK β Y188 antibody	38199981	9	WB	Human

47.Mechanical stiffness promotes skin fibrosis via Piezo1-Wnt2/Wnt11-CCL24 positive feedback loop

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-8568R	Anti-Wnt11 antibody	38267432	9	WB	Human

48.Oil mistparticulate matter exposure induces hyperlipidemia-related inflammation via microbiota/SCFAs/GPR43 axis inhibition and TLR4/NF- κ B activation

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
-----	-------------	-----------	--------------	--------------	-----------------

bs-0816R	Anti-Lysozyme 1 antibody	38199482	8.9	IHC	Rat
bsm-51653M	Anti-regenerating islet-derived protein 3-gamma (Reg3γ) antibody	38199482	8.9	WB	Rat

49.Targeted delivery of rosuvastatin enhances treatment of HHcy-induced atherosclerosis using macrophage membrane-coated nanoparticles

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-24563R	Anti-NLRP3 antibody	10.1016/j.jpha.2024.01.005	8.8	WB	Mouse
bs-23418R	Anti-ABCA1 antibody	10.1016/j.jpha.2024.01.005	8.8	IF	Mouse

50.A multifunctional bimetallic nanoplatform for enhancing radiotherapy via oxidative damage and ferroptosis of cervical cancer

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-2130R	Anti-Ki67 antibody	10.1016/j.seppur.2024.126567	8.6	IHC	Mouse

51.Prickle1-driven basement membrane deposition of the iPSC-derived embryoid bodies is separable from the establishment of apicobasal polarity

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-11259R	anti-Rab6a antibody	38185785	8.5	IHC	Mouse

52.Single-cell RNA-seq of in vitro expanded cells from cranial neural crest reveals a rare odontogenic sub-population

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-8557R	anti-DSPP antibody	38196265	8.5	IHC	Mouse

53.Modulation of apoptosis and Inflammasome activation in chondrocytes: co-regulatory role of Chlorogenic acid

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0032R	anti-Bcl-2 pAb antibody	38169388	8.4	WB	Chicken

54.Faecalibacterium prausnitzii as a potential Antiatherosclerotic microbe

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
-----	-------------	-----------	--------------	--------------	-----------------

bs-1432R	anti-CD68 antibody	38243314	8.4	IF	Mouse
----------	--------------------	----------	-----	----	-------

55.MYH1G-AS is a chromatin-associated lncRNA that regulates skeletal muscle development in chicken

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-9762R	anti-FGF18 antibody	38177995	8.3	IF,WB	Chicken
bs-12653R	anti-SMARCA5 antibody	38177995	8.3	IF,WB	Chicken
bs-6936R	anti-MYF5 antibody	38177995	8.3	WB	Chicken
bs-23966R	anti-SMAD4 antibody	38177995	8.3	WB	Chicken
bs-3419R	anti-phosphorylated (p-) SMAD2 antibody	38177995	8.3	WB	Chicken
bs-3425R	anti-p-SMAD3 antibody	38177995	8.3	WB	Chicken

56.Dual-functional 3D-printed porous bioactive scaffold enhanced bone repair by promoting osteogenesis and angiogenesis

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0195R	anti-CD31 antibody	10.1016/j.mtbio.2024.100943	8.2	IHC	Rabbit
bs-20398R	anti-HIF-1 antibody	10.1016/j.mtbio.2024.100943	8.2	IHC	Rabbit
bs-1134R	anti-Runx-2 antibody	10.1016/j.mtbio.2024.100943	8.2	IHC	Rabbit
bs-4917R	anti-osteocalcin (OCN) antibody	10.1016/j.mtbio.2024.100943	8.2	IHC	Rabbit

57.Drug-free and multifunctional sodium bicarbonate/hyaluronic acid hybrid dressing for synergistic healing of infected wounds

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0295G-HRP	goat anti-rabbit IgG/HRP	38191113	8.2	IHC	Mouse

58.Preparation and antitumor study of intelligent injectable hydrogel: Carboxymethyl chitosan–aldehyde gum Arabic composite graphene oxide hydrogel

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
C01-01002	phosphate-buffered saline (PBS)	10.1016/j.ijbiomac.2024.129429	8.2	Other	Other
C7077	RPMI 1640 medium	10.1016/j.ijbiomac.2024.129429	8.2	Other	Other

C7074	foetal bovine serum	10.1016/j.ijbiom ac.2024.129429	8.2	Other	Other
-------	---------------------	------------------------------------	-----	-------	-------

59.Regulation of the colon-targeted release rate of lactoferrin by constructing hydrophobic ethyl cellulose/pectin composite nanofibrous carrier and its effect on anti-colon cancer activity

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0757R	anti-cdk2 antibody	38242414	8.2	WB	Human
bsm-60698R	anti-p21 antibody	38242414	8.2	WB	Human
bsm-33411M	anti-bcl-2 antibody	38242414	8.2	WB	Human
bs-0127M	anti-bax antibody	38242414	8.2	WB	Human
bsm-52462R	anti-p-jnk antibody	38242414	8.2	WB	Human
bs-2592R	anti-jnk antibody	38242414	8.2	WB	Human
bsm-33284M	anti-caspase 3 antibody	38242414	8.2	WB	Human

60.Gelatin-based dynamic response antioxidant, anti-inflammatory multifunctional hydrogel for enhanced diabetic wound repair

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1035R-FITC	Fluorescein isothiocyanate (FITC)-conjugated anti-CD86 antibody	38253143	8.2	FC	Mouse
bs-4727R-APC	allophycocyanin (APC)-conjugated anti-CD206 antibody	38253143	8.2	FC	Mouse

61.Nanotexture and crystal phase regulation for synergistic enhancement in re-endothelialization on medical pure titanium surface

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
BH0190	mouse anti-CD31 antibody	10.1016/j.mtbio .2024.100968	8.2	IF	Human
bs-10196R	rabbit anti- α -SMA antibody	10.1016/j.mtbio .2024.100968	8.2	IF	Human
bs-0195R	anti-CD31 antibody	10.1016/j.mtbio .2024.100968	8.2	IF	Rabbit

62.Programmed nanocarrier loaded with paclitaxel and dual-siRNA to reverse chemoresistance by synergistic therapy

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
C05-02001	bicinchoninic acid (BCA) protein assay kit	38290632	8.2	Other	Other

63.The Hippo-YAP signaling pathway drives CD24-mediated immune evasion in esophageal squamous cell carcinoma via macrophage phagocytosis

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-20824R	anti-CK5/6 antibody	38168654	8	IHC	Mouse

64.RNA 5-methylcytosine writer NSUN5 promotes hepatocellular carcinoma cell proliferation via a ZBED3-dependent mechanism

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-13553R	rabbit anti-ZBED3 antibody	38182896	8	WB	Human

65.A pH-Sensitive Nanoparticle as Reactive Oxygen Species Amplifier to Regulate Tumor Microenvironment and Potentiate Tumor Radiotherapy

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-3440R	Rabbit anti-phospho-TBK1 (Ser172) antibody	10.2147/IJN.S436160	8	WB	Mouse
bs-3195R	Rabbit anti-phospho-IRF3 (Ser396) antibody	10.2147/IJN.S436160	8	WB	Mouse
bs-7497R	Rabbit anti-TBK1 antibody	10.2147/IJN.S436160	8	WB	Mouse
bs-2993R	Rabbit anti-IRF3 antibody	10.2147/IJN.S436160	8	WB	Mouse

66.Combined Placental Mesenchymal Stem Cells with Guided Nanoparticles Effective Against Diabetic Nephropathy in Mouse Model

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-10498R	rabbit anti-CD3 antibody	38293609	8	IF	Mouse

67.Astragaloside IV negatively regulates Gpr97-TPL2 signaling to protect against hyperhomocysteine-exacerbated sepsis associated acute kidney injury

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-10924R	anti-Gpr97 antibody	10.1016/j.phymed.2024.155346	7.9	WB	Mouse,Rat

68.Gubi Zhitong formula alleviates osteoarthritis in vitro and in vivo via regulating BNIP3L-mediated mitophagy

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
-----	-------------	-----------	--------------	--------------	-----------------

bs-41346R	anti-ACAN antibody	10.1016/j.phym ed.2023.155279	7.9	FC,WB	Rat
BA00208	CCK-8	10.1016/j.phym ed.2023.155279	7.9	Other	Other

69.Catalpolaglycone disrupts mitochondrial thermogenesis by specifically binding to a conserved lysine residue of UCP2 on the proton leak tunnel

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1533R	anti-COX4 antibody	10.1016/j.phym ed.2024.155356	7.9	IF	Mouse

70.Catalpol attenuates ischemic stroke by promoting neurogenesis and angiogenesis via the SDF-1 α /CXCR4 pathway

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-20322R	Rabbit anti-CD31 antibody	10.1016/j.phym ed.2024.155362	7.9	IF	Rat
bsm-33192M	anti-GAP-43 antibody	10.1016/j.phym ed.2024.155362	7.9	WB	Rat

71.Gegen Qinlian tablets delay Alzheimer's disease progression via inhibiting glial neuroinflammation and remodeling gut microbiota homeostasis

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0982R	anti-p-NF- κ B p65 antibody	10.1016/j.phym ed.2024.155394	7.9	IHC	Rat

72.Gualou-Xiebai-Banxia-Tang regulates liver-gut axis to ameliorate Metabolic Syndrome in HFD-fed mice

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-55193R	anti-SCD1 rabbit polyclonal antibody	10.1016/j.phym ed.2023.155320	7.9	WB	Mouse

73.Artesunate improves glucose and lipid metabolism in db/db mice by regulating the metabolic profile and the MAPK/PI3K/Akt signalling pathway

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-6417R	anti-p-PI3K antibody	10.1016/j.phym ed.2024.155382	7.9	WB	Mouse

74.Ebola virus sequesters IRF3 in viral inclusion bodies to evade host antiviral immunity

Skus	ProductNames	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-8335R	anti-STING antibody (Rabbit polyclonal)	38285487	7.7	IF	Human

75.Corylifol A suppresses osteoclastogenesis and alleviates ovariectomy-induced bone loss via attenuating ROS production and impairing mitochondrial function

Skus	ProductNames	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-3666R	anti-Calmodulin antibody	38244329	7.5	WB	Mouse
bs-8519R	anti-p-Calmodulin antibody	38244329	7.5	WB	Mouse

76.MIIP downregulation drives colorectal cancer progression through inducing peri-cancerous adipose tissue browning

Skus	ProductNames	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0061R	anti-β-actin antibody	38245780	7.5	WB	Mouse,Human

77.JAK1/JAK2 degraders based on PROTAC for topical treatment of atopic dermatitis

Skus	ProductNames	Pubmed id	Journal_tier	Applications	CrossReactivity
bsm-54138R	Anti-JAK1 antibody	38262152	7.5	WB,IHC	Mouse
bs-23003R	anti-JAK2 antibody	38262152	7.5	WB,IHC	Mouse

78.Dysregulated lncRNAs regulate human umbilical cord mesenchymal stem cell differentiation into insulin-producing cells by forming a regulatory network with mRNAs

Skus	ProductNames	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-12205R	Mouse anti-Sox17 antibody	38273351	7.5	IF	Human
bs-0924R	rabbit anti-MAFA antibody	38273351	7.5	IF	Human
bs-10379R	rabbit anti-GLUT2 antibody	38273351	7.5	IF	Human

79.Isoflavones isolated from chickpea sprouts alleviate ovariectomy-induced osteoporosis in rats by dual regulation of bone remodeling

Skus	ProductNames	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-20624R	anti-OPG antibody	38290254	7.5	IHC	Rat

BA00208	Cell Counting Kit-8 (CCK-8)	38290254	7.5	Other	Other
bs-1134R	anti-RUNX2 antibody	38290254	7.5	WB	Rat
bs-0578R	anti-Collagen I antibody	38290254	7.5	WB	Rat
bs-6292R	anti-Alkaline Phosphatase antibody	38290254	7.5	WB	Rat

80.Culture and expansion of murine proximal airway basal stem cells

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-7166R	anti-Mucin 5AC antibody	38287366	7.5	IF	Mouse

81.Surface Oxygen Deficiency Enabled Spontaneous Antiprotein Fouling in WO3 Nanosheets for Biosensing in Biological Fluids

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0297P	human Immunoglobulin G (IgG)	38174654	7.4	Other	Other

82.SERS–Temperature Dual-Mode T-type Lateral Flow Strip for Accurate Detection of Free and Total Prostate-Specific Antigens in Blood

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-43141P	Carcinoembryonic antigen (CEA)	38176009	7.4	Other	Other
bs-43140P	CA125	38176009	7.4	Other	Other
bs-43121P	alpha-fetoprotein (AFP)	38176009	7.4	Other	Other

83.ANGPTL4 accelerates ovarian serous cystadenocarcinoma carcinogenesis and angiogenesis in the tumor microenvironment by activating the JAK2/STAT3 pathway and interacting with ESM1

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
C05-01002	phenylmethylsulfonyl fluoride (PMSF)	38212795	7.4	Other	Other
C05-03001	SDS–PAGE loading buffer	38212795	7.4	Other	Other

84.Metabolic reprogramming based on RNA sequencing of gemcitabine-resistant cells reveals the FASN gene as a therapeutic for bladder cancer

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-10090R	Rabbit anti-BMP6 IgG antibody	38218866	7.4	IHC,IF	Human

bs-7081R	Rabbit anti-CARD10 IgG antibody	38218866	7.4	IHC	Human
----------	---------------------------------	----------	-----	-----	-------

85.Silver Nanocluster-Tagged Electrochemiluminescence Immunoassay with a Sole and Narrow Triggering Potential Window

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-43141P	human carcinoembryonic antigen (CEA)	38235596	7.4	Other	Other
bs-43140P	cancer antigen 125 (CA125)	38235596	7.4	Other	Other
bs-43121P	human alpha fetoprotein (AFP)	38235596	7.4	Other	Other
bs-41609P	human prostate specific antigen (PSA)	38235596	7.4	Other	Other

86.Engineered in vivo and in vitro tumor model recapitulates vasculogenic mimicry signatures in melanoma

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bsm-60738R	rabbit anti-Ki-67 antibody	10.1002/btm2.10648	7.4	IF	Mouse

87.Targeting ERS-mitophagy in hippocampal neurons to explore the improvement of memory by tea polyphenols in aged type 2 diabetic rats

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-22173R	anti-PINK1 antibody	38286317	7.4	WB	Rat
bs-2469R	anti-PERK antibody	38286317	7.4	WB	Rat
bs-3330R	anti-p-PERK antibody	38286317	7.4	WB	Rat
bs-4842R	anti-p-eIF2 α antibody	38286317	7.4	WB	Rat

88.PITPNC1 Suppress CD8+ T cell immune function and promote radioresistance in rectal cancer by modulating FASN/CD155

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-2525R	anti-CD155 antibody	38291470	7.4	IHC	Mouse,Human
bs-0388R	anti-IFN γ antibody	38291470	7.4	IF	Human

89.Age-related exacerbation of lung damage after trauma is associated with increased expression of inflammasome components

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
-----	-------------	-----------	--------------	--------------	-----------------

bs-6982R	rabbit anti-mouse neutrophil elastase (NE) antibody	10.3389/fimmu.2023.1253637	7.3	IHC	Mouse
----------	---	----------------------------	-----	-----	-------

90.Human PD-L1 overexpression decreases xenogeneic human T-cell immune responses towards porcine kidneys

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-4941R	anti-CD274 Polyclonal Antibody	38352884	7.3	IHC	Pig

91.A new look at angiogenesis inhibition of geniposide in experimental arthritis by blocking angiopoietin-2 exocytosis

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0468R	Anti-CD31 antibody	38185885	7.2	IHC	Rat

92.Calcium-enriched Carbon Nanoparticles Loaded with Indocyanine Green for Near-Infrared Fluorescence Imaging-Guided Synergistic Calcium Overload, Photothermal Therapy, and Glutathione-Depletion-Enhanced Photodynamic Therapy

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
AK102	GSH assay kit	10.1039/D3TB02690F	7	Other	Other

93.Maintaining the Mitochondrial Quality Control System Was a Key Event of Tanshinone IIA against Deoxynivalenol-Induced Intestinal Toxicity

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-10008R	Anti-Claudin-1 antibody	38247545	7	WB	Pig
bs-2788R	Anti-Claudin-3 antibody	38247545	7	WB	Pig
bs-1329R	Anti-ZO-1 antibody	38247545	7	WB	Pig
bs-1172R	Anti-N-cad antibody	38247545	7	WB	Pig
bs-10011R	Anti-Occludin antibody	38247545	7	WB	Pig

94.Dietary supplementation of benzoic acid and essential oils combination enhances intestinal resilience against LPS stimulation in weaned piglets

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0061R	Anti-β-actin antibody	38238856	7	WB	Pig

95.Sodium Butyrate Alleviates Free Fatty Acid-Induced Steatosis in Primary Chicken Hepatocytes via Regulating the ROS/GPX4/Ferroptosis Pathway

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
-----	-------------	-----------	--------------	--------------	-----------------

bs-41331R	Anti-Transferrin receptor (TFRC) antibody	10.3390/antiox13020140	7	WB	Chicken
bs-5907R	Anti-FTH1 antibody	10.3390/antiox13020140	7	WB	Chicken
bs-19051R	Anti-NCOA4 antibody	10.3390/antiox13020140	7	WB	Chicken
bs-10411R	Anti-COX-2 antibody	10.3390/antiox13020140	7	WB	Chicken
bs-3884R	Anti-GPX4 antibody	10.3390/antiox13020140	7	WB	Chicken
bs-41373R	Anti-GAPDH antibody	10.3390/antiox13020140	7	WB	Chicken

96. Neuroprotection Afforded by an Enriched Mediterranean-like Diet Is Modified by Exercise in a Rat Male Model of Cerebral Ischemia

SKU	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-10412R	anti-KDR antibody	10.3390/antiox13020138	7	FC	Rat

97. Oxidative Stress Plays an Important Role in Glutamatergic Excitotoxicity-Induced Cochlear Synaptopathy: Implication for Therapeutic Molecules Screening

SKU	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-6313R	anti-4HNE antibody	10.3390/antiox13020149	7	WB	Mouse

98. Generation of a new therapeutic d-peptide that induces the differentiation of acute myeloid leukemia cells through A TLR-2 signaling pathway

SKU	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1019R	anti-TLR-2 antibody	38272890	7	CoIP	Human

99. Revitalizing Photoaging Skin through Eugenol in UVB-Exposed Hairless Mice: Mechanistic Insights from Integrated Multi-Omics

SKU	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bsm-33058M	anti-p53 antibody	10.3390/antiox13020168	7	IHC	Mouse

100. Bone Marrow Mesenchymal Stem Cell-Derived Exosomes Ameliorate Aging-Induced BTB Impairment in Porcine Testes by Activating Autophagy and Inhibiting ROS/NLRP3 Inflammasomes via the AMPK/mTOR Signaling Pathway

SKU	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
-----	-------------	-----------	--------------	--------------	-----------------

bs-1329R	anti-ZO-1 Antibody	10.3390/antiox1 3020183	7	WB,IHC	Pig
bs-10011R	anti-Occludin Antibody	10.3390/antiox1 3020183	7	WB	Pig
bs-21509R	anti-Claudin-11 Antibody	10.3390/antiox1 3020183	7	WB	Pig
bs-0061R	anti- β -actin Antibody	10.3390/antiox1 3020183	7	WB	Pig
bs-0295G-HRP	HRP-conjugated goat anti-rabbit secondary antibody	10.3390/antiox1 3020183	7	WB	Pig

101.Selenium deficiency exacerbates ROS/ER stress mediated pyroptosis and ferroptosis induced by bisphenol A in chickens thymus

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-14287R	anti-GSDMD antibody	10.1016/j.jes.20 24.01.002	6.9	IF	Chicken

102.Vitamin D3 reverses immune tolerance and enhances the cytotoxicity of effector T cells in coal pneumoconiosis

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1183R	anti-IL-17 antibody	38218105	6.8	IHC	Mouse
bs-23217R	anti-NF- κ B antibody	38218105	6.8	IF	Mouse

103.Harnessing Cytokine-Induced killer cells to accelerate diabetic wound healing: an approach to regulating post-traumatic inflammation

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bsm-33284M	anti-C-caspase 3 antibody	10.1093/rb/rbad 116	6.7	WB	Mouse

104.Alpha-synuclein dynamics bridge Type-I Interferon response and SARS-CoV-2 replication in peripheral cells

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bsm-54277R	rabbit anti- α -syn mAb antibody	38191441	6.7	IF	Human

105.A strategy for mechanically integrating robust Hydrogel-Tissue hybrid to promote the anti-calcification and endothelialization of bioprosthetic heart valve

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
------	-------------	-----------	--------------	--------------	-----------------

bs-20403R	Rabbit anti-CD68 Antibody	10.1093/rb/rbae003	6.7	IHC	Pig
bs-10498R	Rabbit anti-CD3 Antibody	10.1093/rb/rbae003	6.7	IHC	Pig

106. Bile acids induce liver fibrosis through the NLRP3 inflammasome pathway and the mechanism of FXR inhibition of NLRP3 activation

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-10196R	anti- α -SMA antibody	38172440	6.6	IHC	Human, Mouse
bs-0812R	anti-IL-1 β antibody	38172440	6.6	IHC	Human, Mouse

107. Periodontitis aggravates COPD through the activation of $\gamma\delta$ T cell and M2 macrophage

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1035R	anti-CD86 antibody	38214520	6.4	IF	Mouse

108. Gut microbiota dysbiosis and decreased levels of acetic and propionic acid participate in glucocorticoid-induced glycolipid metabolism disorder

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
C2064	Bouin's fixating solution and ethanol	38226811	6.4	Other	Other

109. The relation between BTK expression and iron accumulation of myeloid cells in multiple sclerosis

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-3577R	anti-Dmt-1 antibody	38254312	6.4	WB	Human
bs-10348R	anti-hepcidin antibody	38254312	6.4	WB	Human

110. IP3RPEP6, a novel peptide inhibitor of IP3 receptor channels that does not affect connexin-43 hemichannels

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1050R	rabbit polyclonal anti-Cx40 antibody	38240350	6.3	WB	Human

111. Amyotrophic lateral sclerosis associated disturbance of iron metabolism is blunted by swim training-role of AKT signaling pathway

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
-----	-------------	-----------	--------------	--------------	-----------------

bs-3579R	anti-SLC40A1 antibody	38171451	6.2	WB	Human
----------	-----------------------	----------	-----	----	-------

112.P75NTR regulates autophagy through the YAP-mTOR pathway to increase the proliferation of interfollicular epidermal cells and promote wound healing in diabetic mice

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1158P	advanced glycosylation end products (AGEs)	38176461	6.2	Other	Other

113.A portable colorimetric sensing platform for rapid and sensitive quantification of dichlorvos pesticide based on Fe-Mn bimetallic oxide nanozyme-participated highly efficient chromogenic catalysis

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
D13034	acetylcholinesterase (AChE)	10.1016/j.aca.2024.342243	6.2	Other	Other

114.Oral administration of LfcinB alleviates DSS-induced colitis by improving the intestinal barrier and microbiota

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bsm-33117M	anti-NF-κB p65 antibody	10.1039/D3FO05236B	6.1	WB	Mouse
bs-3485R	anti-p-p65 antibody	10.1039/D3FO05236B	6.1	WB	Mouse
bs-2592R	anti-JNK antibody	10.1039/D3FO05236B	6.1	WB	Mouse
bsm-52462R	anti-p-JNK antibody	10.1039/D3FO05236B	6.1	WB	Mouse
bsm-33337M	anti-ERK antibody	10.1039/D3FO05236B	6.1	WB	Mouse
bs-1522R	anti-P-ERK antibody	10.1039/D3FO05236B	6.1	WB	Mouse
bs-0637R	anti-p38 antibody	10.1039/D3FO05236B	6.1	WB	Mouse
bs-5476R	anti-p-p38 antibody	10.1039/D3FO05236B	6.1	WB	Mouse
bs-34023R	anti-ZO-1 antibody	10.1039/D3FO05236B	6.1	WB	Mouse
bsm-60433R	anti-Claudin-1 antibody	10.1039/D3FO05236B	6.1	WB	Mouse
bs-10011R	anti-occludin antibody	10.1039/D3FO05236B	6.1	WB	Mouse

115.Gut microbiota dysbiosis contributes to α-synuclein-related pathology associated with C/EBPβ/AEP signaling activation in a mouse model of Parkinson's disease

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0295G-AF594	goat-anti-rabbit IgG H&L/AF594 antibody	38227539	6.1	IF	Mouse

116.Peanut oil odor enhances the immunomodulatory effect on immunosuppressed mice by regulating cAMP signaling pathway via brain-spleen axis

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0036R	Anti-phospho-CREB antibody	10.1039/D3FO03629D	6.1	WB	Mouse
bs-18129R	anti-phospho-IKBA antibody	10.1039/D3FO03629D	6.1	WB	Mouse

117.Purple Sweet Potato Polysaccharide Exerting an Anti-inflammatory Effect via a TLR-Mediated Pathway by Regulating Polarization and Inhibiting the Inflammasome Activation

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-10208R	anti-p-IRAK-4 antibody	38233194	6.1	WB	Mouse
bs-3192R	anti-p-IRAK-1 antibody	38233194	6.1	WB	Mouse
bs-3438R	anti-p-TAK1 antibody	38233194	6.1	WB	Mouse
bs-3229R	anti-p-IKK-α antibody	38233194	6.1	WB	Mouse
bs-5399R	anti-p-IKK-β antibody	38233194	6.1	WB	Mouse
bs-3722R	anti-p-MEK1 + MEK2 antibody	38233194	6.1	WB	Mouse
bs-3274R	anti-p-MEK3 + MKK6 antibody	38233194	6.1	WB	Mouse
bs-3440R	anti-p-TBK1 antibody	38233194	6.1	WB	Mouse
bs-4114R	anti-IKK-ε antibody	38233194	6.1	WB	Mouse
bs-9278R	anti-p-IRF-3 antibody	38233194	6.1	WB	Mouse
bs-3153R	anti-p-c-Fos antibody	38233194	6.1	WB	Mouse
bs-8583R	anti-p-IKK-ε antibody	38233194	6.1	WB	Mouse
bs-0295G-HRP	HRP-labeled goat anti-rabbit IgG	38233194	6.1	WB	Mouse

118.Methionine Promotes Milk Synthesis through the BRCC36-BRG1-mTOR Signaling Axis in Mammary Epithelial Cells

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0466R	anti-β-casein antibody	38240727	6.1	WB	Mouse
bs-6972R	Anti-RNA polymerase II antibody	38240727	6.1	ChIP	Mouse

119.Insight into the Mechanism Underlying the Reduction of Digestibility and IgG/IgE Binding Ability in Ovalbumin during Different High-Temperature Conduction Modes-Induced Glycation

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0295G-HRP	goat anti-rabbit IgG-HRP conjugate	38275225	6.1	ELISA	Rabbit

120. Pancreatic β cell derived extracellular vesicles containing surface preproinsulin are involved in glucose stimulated insulin secretion

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0862R	Anti-insulin antibody	38286207	6.1	WB	Human

121. Activation of autophagy by Citri Reticulatae Semen extract ameliorates amyloid-beta-induced cell death and cognition deficits in Alzheimer's disease

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-23379R	anti-A β 1-42 (CT) rabbit antibody	10.4103/NRR.NRR-D-23-00954	6.1	IF	Mouse

122. Preclinical Investigations on Anti-fibrotic Potential of Long-Term Oral Therapy of Sodium Astragaloside in Animal Models of Cardiac and Renal Fibrosis

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bsm-33287M	anti-TGF- β 1 antibody	38357273	6	WB	Rat

123. Effect of Trichinella spiralis-Derived Antigens on Nonalcoholic Fatty Liver Disease Induced by High-Fat Diet in Mice

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
C05-06002	5%BSA	10.1021/acsptsci.3c00276	6	Other	Other
bs-1074R	Anti-Nrf2 antibody	10.1021/acsptsci.3c00276	6	WB	Mouse
bs-2013R	Anti-phospho-Nrf2 antibody	10.1021/acsptsci.3c00276	6	WB	Mouse
bsm-52830R	Anti-NQO-1 antibody	10.1021/acsptsci.3c00276	6	WB	Mouse

124. A Cataract-Causing Mutation in the TRPM3 Cation Channel Disrupts Calcium Dynamics in the Lens

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-9046R	polyclonal Anti-TRPM3 Antibody	38334649	6	IF	Mouse

125. Sex-Dependent Variations in Hypothalamic Fatty Acid Profile and Neuropeptides in Offspring Exposed to Maternal Obesity and High-Fat Diet

SKU	Product Name	Pubmed ID	Journal Tier	Applications	Cross Reactivity
bs-0922r	Anti-NGN-3 antibody	10.3390/nu16030340	5.9	WB	Mouse

126. PDGF-BB accelerates TSCC via fibroblast lactates limiting miR-26a-5p and boosting mitophagy

SKU	Product Name	Pubmed ID	Journal Tier	Applications	Cross Reactivity
bs-2698R	Anti-MCT4 antibody	38169376	5.8	IHC	Human
bs-18205R	Anti-LDHA antibody	38169376	5.8	IHC	Human
bs-9455R	Anti-Hexokinase II antibody	38169376	5.8	IHC	Human
bs-10249R	Anti-MCT1/SLC16A3 antibody	38169376	5.8	IHC	Human

127. HDAC6 inhibitor ACY-1215 enhances STAT1 acetylation to block PD-L1 for colorectal cancer immunotherapy

SKU	Product Name	Pubmed ID	Journal Tier	Applications	Cross Reactivity
bs-0480R	anti-IFN- γ antibody	38231305	5.8	WB,IF	Mouse

128. Unraveling T cell exhaustion in the immune microenvironment of osteosarcoma via single-cell RNA transcriptome

SKU	Product Name	Pubmed ID	Journal Tier	Applications	Cross Reactivity
bs-4963R	anti-MYC antibody	38280005	5.8	IHC	Human

129. DEAD-box helicase 1 inhibited CD8+ T cell antitumor activity by inducing PD-L1 expression in hepatocellular carcinoma

SKU	Product Name	Pubmed ID	Journal Tier	Applications	Cross Reactivity
bs-0648R	anti-CD8 antibody	38243657	5.7	IHC	Mouse

130. The predictive value of BTG1 for the response of newly diagnosed acute myeloid leukemia to decitabine

SKU	Product Name	Pubmed ID	Journal Tier	Applications	Cross Reactivity
BA00208	Cell Counting Kit-8 (CCK-8) assay	38254153	5.7	Other	Other
C03-05002	red blood cell lysate	38254153	5.7	Other	Other
bs-0061R	anti- β -actin antibody	38254153	5.7	WB	Human

131.Engineered three-dimensional bioactive scaffold for enhanced bone regeneration through modulating transplanted adipose derived mesenchymal stem cell and stimulating angiogenesis

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0295G-FITC	FITC goat anti-rabbit secondary antibody	38344289	5.7	IF	Rabbit
bs-0195R	anti-CD31 antibody	38344289	5.7	IF	Rabbit

132.Ulva lactuca polysaccharide inhibits hepatocellular carcinoma growth by induces the expression of CD5L and activates complement cascade

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-2673R	anti-C5b-9 antibody	10.1016/j.jff.2023.105994	5.6	IHC	Mouse

133.miR-129 Regulates Yak Intramuscular Preadipocyte Proliferation and Differentiation through the PI3K/AKT Pathway

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-10278R	Anti-phospho-PIK3CG antibody	38203803	5.6	WB	Bovine
bs-0686R	Anti-PTEN antibody	38203803	5.6	WB	Bovine
bs-6951R	Anti-AKT1+2+3 antibody	38203803	5.6	WB	Bovine
bs-5193R	Anti-phospho-AKT1+2+3 antibody	38203803	5.6	WB	Bovine
bs-5209R	Anti-phospho-AKT3 antibody	38203803	5.6	WB	Bovine

134.Dietary selenomethionine reduced oxidative stress by resisting METTL3-mediated m6A methylation level of Nrf2 to ameliorate LPS-induced liver necroptosis in laying hens

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0812R	anti-IL-1 β antibody	38176622	5.6	WB	Chicken
bs-0061R	anti- β -actin antibody	38176622	5.6	WB	Chicken

135.Melatonin alleviates necrotizing enterocolitis by reducing bile acid levels through the SIRT1/FXR signalling axis

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-17529R	anti-OST- α antibody	38176339	5.6	WB	Human

136.TSP-1 increases autophagy level in cartilage by upregulating HSP27 which delays progression of osteoarthritis

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
------	-------------	-----------	--------------	--------------	-----------------

bs-20665R	anti-HSP27 antibody	38183909	5.6	WB,IHC	Rat
bs-1105R	anti-Caspase12 antibody	38183909	5.6	WB,IHC	Rat

137. Vitamin D3 promotes gastric cancer cell autophagy by mediating p53/AMPK/mTOR signaling

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-14318R	anti-phospho-AMPK alpha 1 (Ser356) polyclonal antibody	38259281	5.6	WB	Human

138. Involvement of Yes-Associated Protein 1 Activation in the Matrix Degradation of Human-Induced-Pluripotent-Stem-Cell-Derived Chondrocytes Induced by T-2 Toxin and Deoxynivalenol Alone and in Combination

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bsm-52214R	anti-Phospho-YAP (Ser127) antibody	10.3390/ijms25020878	5.6	WB	Human

139. A thermoprotective probiotic function by thermostable lactic acid bacteria and its causal structure

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-20399R	rabbit polyclonal anti-HIF-1 alpha antibody	10.1016/j.jff.2024.106001	5.6	IHC	Chicken

140. By regulating the IP3R/GRP75/VDAC1 complex to restore mitochondrial dynamic balance, selenomethionine reduces lipopolysaccharide-induced neuronal apoptosis

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-20189R	anti-MCU antibody	38219075	5.6	IF,WB	Chicken
C7066	JC-1 $\Delta\psi$ m detection kit	38219075	5.6	Other	Other
bs-0295G-HRP	HRP-labeled goat anti-rabbit IgG	38219075	5.6	WB	Chicken
bs-0557R	anti-MFN1 antibody	38219075	5.6	WB	Chicken
bs-2988R	anti-MFN2 antibody	38219075	5.6	WB	Chicken
bs-11764R	anti-Opa1 antibody	38219075	5.6	WB	Chicken

141. Tilapia (*Oreochromis niloticus*) oligopeptide TBP-1 inhibits hepatocellular carcinoma metastasis by suppressing inflammation and epithelial-mesenchymal transition

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-41200P	human TGF- β 1	10.1016/j.jff.2024.106020	5.6	Other	Other

142.6-formylindolo[3, 2-b]carbazole alters gut microbiota and prevents the progression of ankylosing spondylitis in mice

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bsm-33036M	anti-β-actin antibody	38244515	5.6	WB	Mouse

143.Panaxydol extracted from Panax ginseng inhibits NLRP3 inflammasome activation to ameliorate NASH-induced liver injury

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-41373R	anti-GAPDH antibody	38262161	5.6	WB	Mouse

144.The HIF-1α/EGF/EGFR Signaling Pathway Facilitates the Proliferation of Yak Alveolar Type II Epithelial Cells in Hypoxic Conditions

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-10007R	anti-EGFR antibody	10.3390/ijms25031442	5.6	WB	Bovine
bs-2010R	anti-EGF antibody	10.3390/ijms25031442	5.6	WB	Bovine
bs-2006R	anti-PCNA antibody	10.3390/ijms25031442	5.6	WB	Bovine
bs-20386R	anti-Bax antibody	10.3390/ijms25031442	5.6	WB	Bovine
bsm-33036M	anti-β-actin antibody	10.3390/ijms25031442	5.6	WB	Bovine
bs-0295G-HRP	Goat Anti-Rabbit IgG H&L / HRP	10.3390/ijms25031442	5.6	WB	Bovine

145.Blocking brown adipocyte β3-adrenoceptor attenuates blood-spinal cord barrier impairment and chronic postsurgical pain in a rat model of preoperative stress

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-10966R	anti-β-actin antibody	38278068	5.6	WB	Rat

146.Effects of unburned tobacco smoke on inflammatory and oxidative mediators in the rat prefrontal cortex

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0982R	rabbit polyclonal anti-NF-κB(phospho-Ser 536) antibody	38333013	5.6	WB	Rat
bs-20355R	rabbit polyclonal anti-Nf-κB antibody	38333013	5.6	WB	Rat

147. Marine fungus *Aspergillus* c1. sp metabolite activates the HSF1/PGC-1 α axis, inducing a thermogenic program for treating obesity

SKU	Product Name	Pubmed ID	Journal Tier	Applications	Cross Reactivity
bs-3757R	anti-HSF1 antibody	38333010	5.6	WB	Mouse

148. Traumatic Brain Injury Induces Nociceptin/Orphanin FQ and Nociceptin Opioid Peptide Receptor Expression within 24 Hours

SKU	Product Name	Pubmed ID	Journal Tier	Applications	Cross Reactivity
bs-0181R	anti-NOP receptor antibody	10.3390/ijms25031658	5.6	WB	Rat

149. The Postbiotic Butyrate Mitigates Gut Mucosal Disruption Caused by Acute Ethanol Exposure

SKU	Product Name	Pubmed ID	Journal Tier	Applications	Cross Reactivity
bs-0195R	anti-CD31 antibody	10.3390/ijms25031665	5.6	IHC	Mouse

150. Astragaloside IV alleviates macrophage senescence and d-galactose-induced bone loss in mice through STING/NF- κ B pathway

SKU	Product Name	Pubmed ID	Journal Tier	Applications	Cross Reactivity
bs-23797R	anti-p16 antibody	38290207	5.6	IF, IHC	Mouse

151. GPR137 inactivates Hippo signaling to promote gastric cancer cell malignancy

SKU	Product Name	Pubmed ID	Journal Tier	Applications	Cross Reactivity
bs-4635R	anti-Phospho-MST1(pMST1-Thr183) antibody	38163861	5.5	WB, CoIP	Human
bs-3504R	anti-MST1 antibody	38163861	5.5	WB, CoIP	Human
bs-4663R	anti-MST2 antibody	38163861	5.5	WB	Human
bs-3245R	anti-phospho-LATS1 (pLATS1-Thr1079) antibody	38163861	5.5	WB, CoIP	Human
bs-16270R	anti-GPR137 antibody	38163861	5.5	WB	Human
bs-2007R	anti-PCNA antibody	38163861	5.5	WB	Human

152. Salmonella typhimurium Vaccine Candidate Delivering Infectious Bronchitis Virus S1 Protein to Induce Protection

SKU	Product Name	Pubmed ID	Journal Tier	Applications	Cross Reactivity
-----	--------------	-----------	--------------	--------------	------------------

bs-0296G-FITC	fluorescein isothiocyanate (FITC)-labeled goat anti-mouse IgG antibody	38275762	5.5	IF	Chicken
---------------	--	----------	-----	----	---------

153. Upregulated Selenoprotein I during LPS-induced B cell activation promotes lipidomic changes and is required for effective differentiation into IgM secreting plasma B cells

SKU	Product Name	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-3850R	Anti-BCMA/CD269 polyclonal antibody	38289835	5.5	FC	Mouse

154. Protective effects of *Tinospora cordifolia* miers extract against hepatic and neurobehavioral deficits in thioacetamide-induced hepatic encephalopathy in rats via modulating hyperammonemia and glial cell activation

SKU	Product Name	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0649R	anti-CD68 antibody	38176666	5.4	WB	Rat
bs-6319R	anti-pro IL-1 β antibody	38176666	5.4	WB	Rat

155. Qingda granules alleviate cerebral ischemia/reperfusion injury by inhibiting TLR4/NF- κ B/NLRP3 signaling in microglia

SKU	Product Name	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1363R	anti-IBA-1 antibody	38184025	5.4	IHC	Rat

156. *Withania somnifera* root extract inhibits MGO-induced skin fibroblast cells dysfunction via ECM-integrin interaction

SKU	Product Name	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-10589R	rabbit anti-COL1 antibody	38185262	5.4	WB	Human

157. Enhancing Liver Delivery of Gold Nanoclusters via Human Serum Albumin Encapsulation for Autoimmune Hepatitis Alleviation

SKU	Product Name	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-22386R	anti-p-STAT3 antibody	10.3390/pharmaceutics16010110	5.4	WB	Mouse
bs-0081R	anti-Caspase3 antibody	10.3390/pharmaceutics16010110	5.4	WB	Mouse
bs-3206R	anti-p-JAK2 antibody	10.3390/pharmaceutics16010110	5.4	WB	Mouse

bs-23003R	anti-JAK2 antibody	10.3390/pharmaceutics16010110	5.4	WB	Mouse
-----------	--------------------	-------------------------------	-----	----	-------

158.The accumulation of modular serine protease mediated by a novel circRNA sponging miRNA increases Aedes aegypti immunity to fungus

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0295G-FITC	FITC-conjugated goat anti-rabbit	38233907	5.4	Other	Aedes aegypti
bs-0296G-AF594	Alexa Fluor 594-conjugated goat anti-mouse	38233907	5.4	Other	Aedes aegypti

159.Confirmation of pain-related neuromodulation mechanism of Bushen Zhuangjin Decoction on knee osteoarthritis

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0295G-HRP	HRP goat anti-rabbit IgG	38266947	5.4	WB	Rat

160.Porcine deltacoronavirus nsp5 antagonizes type I interferon signaling by cleaving IFIT3

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bsm-33346M	mouse anti-Flag tag antibody	38289117	5.4	IF,WB	Human,Pig

161.EEBR induces Caspase-1-dependent pyroptosis through the NF-κB/NLRP3 signalling cascade in non-small cell lung cancer

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0061R	anti-β-actin antibody	38214430	5.3	WB	Human,Mouse
bs-10442R	anti-Caspase-1 antibody	38214430	5.3	WB,IF	Human,Mouse
bs-4986R	anti-IL-18 antibody	38214430	5.3	IF	Mouse
bs-25615R	anti-IL-1β antibody	38214430	5.3	IF	Mouse

162.Transcriptomic dynamics governing serotonergic dysregulation in the dorsal raphe nucleus following mild traumatic brain injury

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-2590R	rabbit polyclonal anti-REST antibody	38246304	5.3	IHC	Mouse

163.CDK12 regulates angiogenesis of advanced prostate cancer by IGFBP3

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0295G-HRP	Goat Anti-Rabbit IgG Antibody (H+L) HRP	38186306	5.2	CoIP	Human

164.Tumor-associated macrophages and PD-L1 in prostate cancer: a possible key to unlocking immunotherapy efficacy

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-2527R	rabbit anti-CD163 antibody	38189834	5.2	IHC	Human

165.MICA+ Tumor Cell Upregulated Macrophage-Secreted MMP9 via PROS1-AXL Axis to Induce Tumor Immune Escape in Advanced Hepatocellular Carcinoma (HCC)

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0832R	anti-MICA antibody	10.3390/cancers16020269	5.2	IHC	Human

166.Puerarin Attenuates Iron Overload-Induced Ferroptosis in Retina through a Nrf2-Mediated Mechanism

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-2075R	anti-HO-1 antibody	38196088	5.2	WB	Mouse

167.Tongbian decoction restores intestinal microbiota and activates 5-hydroxytryptamine signaling: implication in slow transit constipation

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
C02-04002	DAPI	38287961	5.2	Other	Other

168.S100A9 Induces Macrophage M2 Polarization and Immunomodulatory Role in the Lesion Site After Spinal Cord Injury in Rats

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-8481R	rabbit polyclonal anti-MHC II antibody	38206470	5.1	WB	Rat

169.Discovery of 2,4-diarylaminopyrimidine derivatives bearing sulfonamide moiety as novel FAK inhibitors

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
-----	-------------	-----------	--------------	--------------	-----------------

BA00204	Cell cycle kit	10.1016/j.bioorg .2024.107134	5.1	Other	Other
---------	----------------	----------------------------------	-----	-------	-------

170.Mechanism of multi-organ compensation under different iodine intake in pregnant rats: results from a repeated-measures study of iodine metabolism

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bsm-52262R	rabbit anti-GAPDH (Loading Control) antibody	38170273	5	WB	Rat

171.Nicotinamide phosphoribose transferase facilitates macrophage-mediated pulmonary fibrosis through the Sirt1-Smad7 pathway in mice

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-3393R	anti-Phospho-Sirt1 antibody	38280463	5	WB	Mouse

172.Assessing trans-endothelial transport of nanoparticles for delivery to abdominal aortic aneurysms

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1310R	anti-Integrin α V β 3 antibody	38192169	4.9	IF	Rat

173.The artificial amino acid change in the sialic acid-binding domain of the hemagglutinin neuraminidase of newcastle disease virus increases its specificity to HCT 116 colorectal cancer cells and tumor suppression effect

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-4529R	anti-NDV HN antibody	38178138	4.8	IHC	Mouse

174.Lactate accumulation induced by Akt2-PDK1 signaling promotes pulmonary fibrosis

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
BA00208	CCK-8 assay	38226859	4.8	Other	Other

175.SMPDL3B is palmitoylated and stabilized by ZDHHC5, and its silencing aggravates diabetic retinopathy of db/db mice: Activation of NLRP3/NF- κ B pathway

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-12642R	rabbit-sourced SMPDL3B antibody	38266744	4.8	IHC,IF	Human,Mouse

176.Establishment and analysis of immortalized chicken skeletal muscle satellite cell lines1

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0296G-Cy3	Cy3-goat anti-mouse IgG (H+L)	10.1016/j.jia.2024.01.034	4.8	IF	Chicken
C-0005	goat serum	10.1016/j.jia.2024.01.034	4.8	Other	Other

177.Roflupram alleviates autophagy defects and reduces mutant hSOD1-induced motor neuron damage in cell and mouse models of amyotrophic lateral sclerosis

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-3602R	anti-ULK1 antibody	38218579	4.7	WB	Mouse

178.Genetic and Biological Characteristics of Duck-Origin H4N6 Avian Influenza Virus Isolated in China in 2022

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0310G-HRP	HRP-conjugated goat anti-chicken antibody	10.3390/v16020207	4.7	Other	Duck-Origin H4N6 Avian Influenza Virus

179.Astragalus Polysaccharide Alleviates Ulcerative Colitis by Regulating the Balance of mTh17/mTreg Cells through TIGIT/CD155 Signaling

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-2525R	anti-rabbit CD155 antibody	10.3390/molecules29010241	4.6	WB,IF	Mouse
bs-6181R	anti-rabbit CD113 antibody	10.3390/molecules29010241	4.6	WB,IF	Mouse

180.Notch intracellular domains form transcriptionally active heterodimeric complexes on sequence-paired sites

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bsm-0966M-Biot	biotinylated anti-HA antibody	38168761	4.6	ChIP	Human

181.Colorectal cancer with low SLC35A3 is associated with immune infiltrates and poor prognosis

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
-----	-------------	-----------	--------------	--------------	-----------------

bs-0766R	rabbit anti-CD4 antibody	38172565	4.6	IHC	Human
----------	--------------------------	----------	-----	-----	-------

182. Proteomics and phosphoproteomics to study Tuina reverses capsule fibrosis in frozen shoulder: a research report based on rats

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-2715R	anti-Thbs1 antibody	38182704	4.6	WB	Rat
bs-0061R	anti-β-Actin antibody	38182704	4.6	WB	Rat

183. Hsa_circ_0081065 exacerbates IH-induced EndMT via regulating miR-665/HIF-1α signal axis and HIF-1α nuclear translocation

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0195R	anti-CD31 antibody	38195914	4.6	WB	Human

184. Transcription factor FoxO1 regulates myoepithelial cell diversity and growth

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-7111R	anti-Eda2r antibody	38212454	4.6	IF, WB	Mouse

185. Ritonavir Has Reproductive Toxicity Depending on Disrupting PI3K/PDK1/AKT Signaling Pathway

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-3351R	anti-phospho-PTEN antibody	10.3390/toxics12010073	4.6	WB	Dog

186. Phomopsterone B Alleviates Liver Fibrosis through mTOR-Mediated Autophagy and Apoptosis Pathway

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-4563R	anti-bcl-2 antibody	10.3390/molecules29020417	4.6	WB	Human

187. Acute Kidney Injury Caused by Rhabdomyolysis Is Ameliorated by Serum Albumin-Based Supersulfide Donors through Antioxidative Pathways

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1278R	anti-8-OHdG Polyclonal Antibody	38256961	4.6	IHC	Mouse

188.All-trans retinoic acid alleviates collagen-induced arthritis and promotes intestinal homeostasis

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-10008R	anti-claudin-1 Antibody	38245637	4.6	IHC	Rat
bs-34023R	anti-ZO-1 Antibody	38245637	4.6	IHC	Rat

189.Xiao Cheng Qi Decoction, an Ancient Chinese Herbal Mixture, Relieves Loperamide-Induced Slow-Transit Constipation in Mice: An Action Mediated by Gut Microbiota

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-2127R	anti-5-HT4 receptor Antibody	10.3390/ph17020153	4.6	WB	Mouse

190.SN-38, an active metabolite of irinotecan, inhibits transcription of nuclear factor erythroid 2-related factor 2 and enhances drug sensitivity of colorectal cancer cells

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0295PC	normal rabbit IgG	38270247	4.6	Other	Other

191.Unique E2-binding specificity of artificial RING fingers in cancer cells

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-8348R	rabbit anti-UbcH5b antibody	38291082	4.6	WB	Human

192.IRF4 induces M1 macrophage polarization and aggravates ulcerative colitis progression by the Bcl6-dependent STAT3 pathway

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0908R	anti-JAK2 Antibody	38164749	4.5	WB	Mouse

193.Combined exposure to lead and microplastics increased risk of glucose metabolism in mice via the Nrf2/NF-κB pathway

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1074R	Rabbit Anti-Nrf2 antibody	38180308	4.5	WB	Mouse
bs-2075R	Anti-HO-1 antibody	38180308	4.5	WB	Mouse
bs-3648R	Anti-Kelch-like ECH-associated protein 1 (Keap-1) antibody	38180308	4.5	WB	Mouse

194. Electroacupuncture improves articular microcirculation and attenuates cartilage hypoxia in a male rabbit model of knee osteoarthritis

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0101R	Anti-PKM2 antibody	10.1016/j.jtcme.2024.01.002	4.5	IHC,WB	Rabbit

195. Dendrobium officinale polysaccharide alleviates thiacloprid-induced kidney injury in quails via activating the Nrf2/HO-1 pathway

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1074R	Anti-Nrf2 antibody	38224485	4.5	WB	Chicken
bs-0032R	Anti-B-cell lymphoma gene 2 (Bcl-2) antibody	38224485	4.5	WB	Chicken
bs-28034R	Anti-Bcl-2-associated X protein (Bax) antibody	38224485	4.5	WB	Chicken
bs-2075R	Anti-HO-1 antibody	38224485	4.5	WB	Chicken
bsm-33284M	Anti-caspase-3 antibody	38224485	4.5	WB	Chicken

196. A novel lipid metabolism-based risk model associated with immunosuppressive mechanisms in diffuse large B-cell lymphoma

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-9102R	anti-ARSK antibody	38254162	4.5	IHC,IF	Human

197. Identification of genes regulated by trait sensitivity to negative feedback and prolonged alcohol consumption in rats

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-2126R	rabbit anti-5-HT3A antibody	38172401	4.4	WB	Rat

198. FOXD1 promotes proliferation, migration, and epithelial-mesenchymal transformation in esophageal squamous cell carcinoma by targeting SNAI1

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1165R	anti-β-catenin antibody	10.1016/j.procbio.2024.01.006	4.4	WB	Human
bs-1172R	anti-N-cadherin antibody	10.1016/j.procbio.2024.01.006	4.4	WB	Human
bs-1371R	anti-SNAI1 antibody	10.1016/j.procbio.2024.01.006	4.4	WB	Human
bs-1519R	anti-E-cadherin antibody	10.1016/j.procbio.2024.01.006	4.4	WB	Human

bs-0756R	anti-vimentin antibody	10.1016/j.procbio.2024.01.006	4.4	WB	Human
----------	------------------------	-------------------------------	-----	----	-------

199.Dermal FOXO3 activity in response to Wnt/ β -catenin signaling is required for feather follicle development of goose embryos (*Anser cygnoides*)

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0127R	anti-Bax protein antibody	38330682	4.4	WB	Goose
bs-0032R	anti-Bcl2 protein antibody	38330682	4.4	WB	Goose
bs-2006R	anti-PCNA protein antibody	38330682	4.4	WB	Goose
bs-23663R	anti- β -catenin protein antibody	38330682	4.4	WB	Goose
bs-1843R	anti-LEF1 protein antibody	38330682	4.4	WB	Goose

200.Role of Nrf2/HO-1 pathway on inhibiting activation of ChTLR15/ChNLRP3 inflammatory pathway stimulated by *E. tenella* sporozoites

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1074R	anti-Nrf2 antibody	10.1016/j.psj.2024.103445	4.4	WB	Chicken

201.Inflammatory responses and barrier disruption in the trachea of chicks following *Mycoplasma gallisepticum* infection: a focus on the TNF- α -NF- κ B/MLCK pathway

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-2513R	anti-p-I κ B α antibody	38225621	4.4	WB	Chicken
bs-0465R	anti-p65 antibody	38225621	4.4	WB	Chicken
bs-0982R	anti-p-p65 antibody	38225621	4.4	WB	Chicken
bsm-33207M	anti-TNF- α antibody	38225621	4.4	WB	Chicken

202.Sodium butyrate alleviates free fatty acid-induced steatosis in primary chicken hepatocytes via the AMPK/PPAR α pathway

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1115R	anti-AMP-activated protein kinase α 1 (AMPK α 1) antibody	10.1016/j.psj.2024.103482	4.4	WB	Chicken
bs-5551R	anti-p-AMPK α 1 antibody	10.1016/j.psj.2024.103482	4.4	WB	Chicken

bs-23779R	anti-carnitine palmitoyl-transferase 1A (CPT1A) antibody	10.1016/j.psj.2024.103482	4.4	WB	Chicken
-----------	--	---------------------------	-----	----	---------

203.Pathogenic Characteristics of an Infection with Canine Influenza Virus and Streptococcus equi subsp. zooepidemicus Alone or in Combination in Mice

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-8853R	anti-Smad 2/3 (Thr8) antibody	10.1155/2024/2237621	4.3	IF,WB	Mouse

204.Exploring the prognostic potential of m6A methylation regulators in low-grade glioma: implications for tumor microenvironment modulation

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-18553R	anti-ZCCHC4 antibody	38173044	4.2	IHC	Human

205.Myogenic Differentiation and Immunomodulatory Properties of Rat Adipose-Derived Mesenchymal Stem/Stromal Cells

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0521R-FITC	FITC anti-CD44 antibody	10.3390/biology13020072	4.2	ICC	Rat
bs-0646R-FITC	FITC anti-CD34 antibody	10.3390/biology13020072	4.2	ICC	Rat
bs-0295P-FITC	Rabbit IgG isotype control	10.3390/biology13020072	4.2	Other	Other
bs-2442R-FITC	FITC anti-MyoD antibody	10.3390/biology13020072	4.2	ICC	Rat

206.Identification and Verification of Hub Mitochondrial Dysfunction Genes in Osteoarthritis Based on Bioinformatics Analysis

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-10738R	anti-SLC19A2 antibody	38292759	4.1	WB	Human

207.Parapyruvate Induces Neurodegeneration in C57BL/6JNarl Mice via Inhibition of the α -Ketoglutarate Dehydrogenase Complex

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-20444R	anti-Tau antibody	38343944	4.1	IHC	Mouse

208.Adipocyte subpopulations mediate lipolysis and obesity induced insulin resistance

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-4888R	anti-p-PPAR γ S273 antibody	38174979	4	WB	Mouse

209.Dysfunctional Bladder Morphology and Functional Impairments Are Identified in the Alzheimer's Disease APP NL-G-F/NL-G-F Murine Model

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-21455R	anti- β 2AR antibody	38160353	4	WB	Mouse

210.Fermented Natto powder alleviates obesity by regulating LXR pathway and gut microbiota in obesity rats

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-2536R	anti-SREBP-2 antibody	38192042	4	WB	Rat
bs-21429R	anti-CYP7A1 antibody	38192042	4	WB	Rat
bsm-52822R	anti-HMGCR antibody	38192042	4	WB	Rat
bs-41373R	anti-GAPDH antibody	38192042	4	WB	Rat
bs-18451R	anti-liver X receptors (LXR) antibody	38192042	4	WB	Rat

211.Sodium butyrate restricts neutrophils migration and NETs formation through reducing macrophage-derived CXCL16 in calculous cholecystitis

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1441R	anti-CXCL16 antibody	10.1016/j.heliyon.2024.e25189	4	IHC	Human
bs-4943R	anti-MPO antibody	10.1016/j.heliyon.2024.e25189	4	IF	Mouse

212.Tanshinone IIA attenuates osteoarthritis via inhibiting aberrant angiogenesis in subchondral bone

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-10589R	anti-Col II antibody	38253247	3.9	IF	Mouse
bs-1223R	anti-Aggregan antibody	38253247	3.9	IF	Mouse
bs-10250R	anti-MMP13 antibody	38253247	3.9	IF	Mouse
bs-3573R	anti-ADAMTS5 antibody	38253247	3.9	IF	Mouse

213.Berberine alleviates ischemia reperfusion injury induced AKI by regulation of intestinal microbiota and reducing intestinal inflammation

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0921R	anti-Sirt1 antibody	38291383	3.9	IF	Rat

214.microRNA-622 upregulates cell cycle process by targeting FOLR2 to promote CRC proliferation

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-9595R	anti-FOLR2 antibody	38166756	3.8	IHC	Mouse

215.ALDH3A1 upregulation inhibits neutrophils N2 polarization and halts oral cancer growth

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0982R	anti-p-NF- κ B antibody	38225738	3.8	IHC,WB	Human
bs-10123R	anti-IKK antibody	38225738	3.8	WB	Human
bs-3237R	anti-p-IKK antibody	38225738	3.8	WB	Human
bs-2513R	anti-p-I κ B antibody	38225738	3.8	WB	Human

216.LOX-1 regulation of H-type vascular endothelial cell regeneration in hyperglycemia

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-2044R	anti-LOX-1 antibody	38244081	3.8	WB,IF	Rat
bs-20399R	anti-HIF-1 α antibody	38244081	3.8	WB	Rat
bs-34032R	anti-VEGF antibody	38244081	3.8	WB	Rat

217.Jiaohong pills attenuate neuroinflammation and amyloid- β protein-induced cognitive deficits by modulating the mitogen-activated protein kinase/nuclear factor kappa-B pathway

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0982R	anti-phospho-p65-NF- κ B(p-p65) antibody	38177948	3.7	WB	Mouse
bs-5476R	anti-phospho-p38 MAPK (p-p38) antibody	38177948	3.7	WB	Mouse

218.Hypoxic stress caused apoptosis of MDBK cells by p53/BCL6-mitochondrial apoptosis pathway

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bsm-60223R	anti-BCL6 antibody	38178593	3.7	WB	Bovine
bsm-33036M	anti- β -actin antibody	38178593	3.7	WB	Bovine
bs-0013R	anti-Cytc antibody	38178593	3.7	WB	Bovine
bsm-33199M	anti-Caspase-3 antibody	38178593	3.7	WB	Bovine
bs-0049R	anti-Caspase-9 antibody	38178593	3.7	WB	Bovine

219. Carbon disulfide induces accumulation of TDP-43 in the cytoplasm and mitochondrial dysfunction in rat spinal cords

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-19882R	anti-Phospho-Parkin(Ser 65) antibody	38185999	3.7	CoIP	Rat

220. MiR-92a-2-5p suppresses esophageal squamous cell carcinoma cell proliferation and invasion by targeting PRDX2

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
BA00208	CCK-8 solution	38211680	3.7	Other	Other

221. Integration of RRBS and RNA-seq unravels the regulatory role of DNMT3A in porcine Sertoli cell proliferation

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
C5027	RIPA lysis buffer	38264208	3.7	Other	Other

222. A novel assay to measure low-density lipoproteins binding to proteoglycans

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-5178R	Rabbit anti-chondroitin sulfate antibody	38295021	3.7	ICC	Mouse

223. Bombesin receptor-activated protein homolog deficiency altered the pattern of pathological changes of psoriasis - like skin lesion in mice

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-2508R	anti-CD11C antibody	38169666	3.6	IF	Mouse

224. Initial characterization of a transgenic mouse with overexpression of the human D1-dopamine receptor in the heart

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1007R	rabbit polyclonal anti-DRD1 antibody	38177456	3.6	WB	Mouse

225. DNAJA1 positively regulates amino acid-stimulated milk protein and fat synthesis in bovine mammary epithelial cells

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
-----	-------------	-----------	--------------	--------------	-----------------

bs-0466R	anti- β -casein antibody	10.1002/cbf.3918	3.6	WB	Bovine
----------	--------------------------------	------------------	-----	----	--------

226.Neuroprotective effects of maternal melatonin administration in early-onset placental insufficiency and fetal growth restriction

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-21511R	anti-MTR1B antibody	38225450	3.6	IHC	Sheep

227.Study of iodine transport and thyroid hormone levels in the human placenta under different iodine nutritional status

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-10900R	rabbit anti-GAPDH (Loading Control) antibody	38221821	3.6	WB	Human

228.Q-switched 1064 nm Nd: YAG laser restores skin photoageing by activating autophagy by TGF β 1 and ITGB1

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0415R	anti-TIMP1 antibody	10.1111/exd.15006	3.6	IF	Mouse

229.Propranolol inhibits EMT and metastasis in breast cancer through miR-499-5p-mediated Sox6

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bsm-33170M	anti-vimentin antibody	38294713	3.6	WB	Mouse
D10393	Hematoxylin	38294713	3.6	Other	Other

230.Transcription Factor SATB2 Regulates Skeletal Muscle Cell Proliferation and Migration via HDAC4 in Pigs

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-23102R	anti-KI67 antibody	10.3390/genes15010065	3.5	WB	Pig

231.Reuterin isolated from the probiotic Lactobacillus reuteri promotes periodontal tissue regeneration by inhibiting Cx43-mediated the intercellular transmission of endoplasmic reticulum stress

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-8680R	anti-IRE1 antibody	38193526	3.5	WB	Other

232. Interferon-tau protects bovine endometrial epithelial cells against inflammatory injury by regulating the PI3K/AKT/ β -catenin/FoxO1 signaling axis

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-28034R	anti-bax antibody	38220437	3.5	WB	Bovine
bsm-33411M	anti-bcl-2 antibody	38220437	3.5	WB	Bovine
bs-3485R	anti-p-p65 antibody	38220437	3.5	WB	Bovine
bs-20160R	anti-p65 antibody	38220437	3.5	WB	Bovine
bs-20595R	anti-TLR4 antibody	38220437	3.5	WB	Bovine

233. Thioredoxin-2 (TXN2) suppresses hydrogen peroxide-activated nuclear factor kappa B signaling via alleviating oxidative stress in bovine adipocytes

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-4256R	anti-TXN2 antibody	38246558	3.5	WB	Bovine

234. Demystifying the Landscape of Endometrial Immune Microenvironment in Luteal-Phase from Cuprotosis: Implications for the Mechanism and Treatment of RPL

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-13008R	anti-DLST antibody	38253297	3.5	IF	Mouse
bs-18298R	anti-LIPT1 antibody	38253297	3.5	IF	Mouse

235. Renal regenerative capacity related to stem cell reserve in nephrectomized rats

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0778R	anti-CD90 antibody	38206410	3.4	IF	Rat

236. Effects of hypoxia in the diabetic corneal stroma microenvironment

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-11054R	anti-Keratocan antibody	38224848	3.4	WB	Human

237. Polygalasaponin F ameliorates middle cerebral artery occlusion-induced focal ischemia / reperfusion injury in rats through inhibiting TXNIP/NLRP3 signaling pathway

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-10445R	Rabbit Anti-NSE antibody	10.1016/j.jneur oim.2023.57828 1	3.3	IF	Rat

238.Effect of 4-nonylphenol (4-NP) on sperm function: Insights into the PI3K/PDK1/AKT signaling pathway during capacitation

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-3351R	anti-phospho-PTEN (Ser380 + Thr382 + Thr383) polyclonal antibody	38246476	3.3	WB	Pig

239.Proteomic Characteristics of the Prefrontal Cortex and Hippocampus in Mice with Chronic Ketamine-Induced Anxiety and Cognitive Impairment

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0222R	anti-NMDAR2B antibody	38266908	3.3	WB	Mouse

240.Fructose aggravates copper-deficiency-induced cardiac remodeling by inhibiting SERCA2a

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0296G-AF594	Alexa Fluor 594-conjugated Goat Anti-Mouse IgG Antibody (H+L)	38271051	3.3	IF	Rat

241.Special Staining and Protein Expression of VEGF/EGFR and P53/NF-κB in Cryptorchid Tissue of Erhualian Pigs

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-40296G-HRP	Goat Anti-Mouse IgG (H+L)/HRP	10.3390/life14010100	3.2	IHC	Pig
bs-40295G-HRP	Goat Anti-Rabbit IgG (H+L)/HRP	10.3390/life14010100	3.2	IHC	Pig
bs-40295G-IRDye	IRDye® 800CW Goat Anti-Rabbit IgG H&L	10.3390/life14010100	3.2	IF	Pig

242.Downregulation of transcription 1 hinders the replication of Dabie bandavirus by promoting the expression of TLR7, TLR8, and TLR9 signaling pathway

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0296G-HRP	HRP goat anti-mouse antibody	38194758	3.2	WB	Human

243. Metformin Increases the Response of Cholangiocarcinoma Cells to Gemcitabine by Suppressing Pyruvate Kinase M2 to Activate Mitochondrial Apoptosis

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-1519R	anti-E-cadherin antibody	38170336	3.1	IF	Human

244. Role of TFEB-autophagy lysosomal pathway in palmitic acid induced renal tubular epithelial cell injury

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-2188R	anti-GAPDH antibody	10.1016/j.bbrc.2024.149472	3.1	WB	Human

245. Huaier Regulates Oxaliplatin Resistance in Colorectal Cancer by Regulating Autophagy and Inhibiting the Wnt/ β -catenin Signaling Pathway

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-1700R	anti-Wnt3a antibody	31345004	3.1	WB	Human

246. The serine protease 2 gene regulates lipid metabolism through the LEP/ampk α 1/SREBP1 pathway in bovine mammary epithelial cells

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-10501R	anti-PRSS2 antibody	38271832	3.1	WB	Bovine
bs-1402R	anti-SREBP1 antibody	38271832	3.1	WB	Bovine
bs-10344R	anti-AMPK α 1 antibody	38271832	3.1	WB	Bovine
bsm-52132R	anti-P-AMPK α 1 (Ser496) antibody	38271832	3.1	WB	Bovine
bsm-33036M	anti- β -actin antibody	38271832	3.1	WB	Bovine

247. Caveolin 1 Regulates the Tight Junctions between Sertoli Cells and Promotes the Integrity of Blood-Testis Barrier in Yak via the FAK/ERK Signaling Pathway

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
IHC0111	Histone H3 Ready-To-Use IHC Kit	10.3390/ani14020183	3	IHC	Bovine
C-0003	30-diaminobenzidine (DAB) kit	10.3390/ani14020183	3	Other	Other
bs-0295G-HRP	horseradish peroxidase (HRP)-conjugated goat anti-rabbit	10.3390/ani14020183	3	WB	Bovine
bs-0368G-HRP	Goat Anti-Mouse IgM / HRP	10.3390/ani14020183	3	WB	Bovine

248. Immunofluorescence study of cytoskeleton in endothelial cells induced with malaria sera

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0715R	anti- β -tubulin antibody	38183117	3	IF	Human

249. The antimicrobial peptide Abaecin alleviates colitis in mice by regulating inflammatory signaling pathways and intestinal microbial composition

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bsm-33117M	anti-p65 antibody	38242174	3	WB	Mouse
bs-3485R	anti-p-p65 antibody	38242174	3	WB	Mouse
bs-0637R	anti-p38 antibody	38242174	3	WB	Mouse
bs-5476R	anti-p-p38 antibody	38242174	3	WB	Mouse
bs-2592R	anti-JNK antibody	38242174	3	WB	Mouse
bsm-52462R	anti-p-JNK antibody	38242174	3	WB	Mouse
bsm-33337M	anti-ERK antibody	38242174	3	WB	Mouse
bs-1522R	anti-p-ERK antibody	38242174	3	WB	Mouse
bs-10011R	anti-occludin antibody	38242174	3	WB	Mouse
bs-0790R	anti-claudin-1 antibody	38242174	3	WB	Mouse
bs-34023R	anti-ZO-1 antibody	38242174	3	WB	Mouse

250. Mitochondrial Acid 5 Increases Ram Sperm Quality by Improving Mitochondrial Function during Storage at 4 °C

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-3955R	anti-MT-ND6 antibody	10.3390/ani14030368	3	WB	Sheep

251. C6 Ceramide Inhibits Canine Mammary Cancer Growth and Metastasis by Targeting EGR3 through JAK1/STAT3 Signaling

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-6448R	anti-EGR3 antibody	10.3390/ani14030422	3	WB	Dog

252. SP1 and KROX20 Regulate the Proliferation of Dermal Papilla Cells and Target the CUX1 Gene

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-17135R	rabbit polyclonal anti-SOX18 antibody	10.3390/ani14030429	3	IF	Sheep

253. Inhibition of FoxO1 alleviates polycystic ovarian syndrome by reducing inflammation and the immune response

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-2537R	rabbit anti-FoxO1 antibody	38189995	2.9	IHC, WB	Rat

254. A deep-red xanthene-based highly sensitive fluorescent probe for detection of hypochlorite

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
D10171	Glutathione (GSH)	38286607	2.9	Other	Other
D10209	cysteine (Cys)	38286607	2.9	Other	Other

255. Cuproptosis mediates copper-induced testicular spermatogenic cell death

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-0295G-HRP	goat anti-rabbit IgG heavy and light chains/horseradish peroxidase (H&L/HRP)	38284772	2.9	WB	Mouse
bs-0296G-HRP	goat anti-mouse IgG H&L/HRP	38284772	2.9	WB	Mouse

256. Evaluation of cell proliferation marker CDC-7 in salivary adenoid cystic carcinoma

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-5148R	anti-CDC-7 antibody	38183819	2.8	IHC	Human

257. Cell-permeable PI3 kinase competitive peptide inhibits KIT mutant mediated tumorigenesis of gastrointestinal stromal tumor (GIST)

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-0294D-HRP	HRP conjugated donkey anti-goat IgG	38206538	2.8	IP, WB	Mouse, Monkey

258. Establishment and characterization of cisplatin-resistant cell lines from canine mammary gland tumors

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-11208R	anti-Sox4 antibody	10.1016/j.theriogenology.2024.01.017	2.8	WB	Dog

BA00101	Annexin V-FITC/PI kit	10.1016/j.theriogenology.2024.01.017	2.8	Other	Other
---------	-----------------------	--------------------------------------	-----	-------	-------

259. Carvedilol induces pyroptosis through NLRP3-caspase1-ASC inflammasome by nuclear migration of NF-κB in prostate cancer models

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-0465R	anti-NF-κB antibody	38270665	2.8	WB	Human
bsm-52110R	anti-Histone-H3.3 antibody	38270665	2.8	WB	Human

260. Follicular fluid exosomes inhibit expression of BTG2 and promote glucose uptake in granulosa cells by delivering miR-21-5p

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-12710R	anti-p-IRS1 antibody	38301506	2.8	WB	Pig

261. Schizandrin A attenuates early brain injury following subarachnoid hemorrhage through suppressing neuroinflammation

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bsk12002	Mouse TNF-α ELISA kit	38285214	2.8	ELISA	Mouse
bsk12015	Mouse IL-1β ELISA kit	38285214	2.8	ELISA	Mouse
bsk12004	Mouse IL-6 ELISA kit	38285214	2.8	ELISA	Mouse

262. MHC-I in the hippocampus promotes comorbid depressive symptoms in bone cancer pain via the upregulation of microglial TREM2/DAP12 signaling

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-12630R	rabbit anti-DAP12 antibody	38176616	2.7	IHC	Mouse

263. Precise construction of hyaluronic acid (HA)-coated drug-loaded iron oxide metal-organic frameworks: Investigation of cell death mechanism in neuroblastoma carcinoma cells

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
S0088	Hyaluronic acid kit	10.1557/s43578-023-01244-x	2.7	Other	Other
D13126	hyaluronidase	10.1557/s43578-023-01244-x	2.7	Other	Other

264. GSDMD suppresses keratinocyte differentiation by inhibiting FLG expression and attenuating KCTD6-mediated HDAC1 degradation in atopic dermatitis

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-16100R	anti-FLG2 antibody	38250727	2.7	WB	Human

265.Short-Time Alternating Current Electrical Stimulation and Cell Membrane-Related Components

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-10927R	anti-aquaporins 5 antibody	10.3390/app14020812	2.7	IF	Human
bs-2506R	anti-aquaporins 7 antibody	10.3390/app14020812	2.7	IF	Human
bs-2060R	anti-aquaporins 9 antibody	10.3390/app14020812	2.7	IF	Human

266.Immortalization of American miniature horse-derived fibroblast by cell cycle regulator with normal karyotype

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0917R	rabbit polyclonal anti-Brdu antibody	38288466	2.7	Other	Horse
bs-8533R	anti-vimentin antibody	38288466	2.7	IF	Horse
bs-20702R	anti-desmin antibody	38288466	2.7	IF	Horse

267.Memantine versus Ginkgo biloba Extract: A Comparative Study on Cognitive Dysfunction Treatment in a Novel Rat Model

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
V5301	anti-A β 1-42 antibody	38286405	2.7	IHC	Rat

268.Scutellarin alleviates tensile stress-induced proliferation and migration of venous smooth muscle cells via mediating the p38 MAPK pathway

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-5476R	Rabbit Polyclonal anti-p-P38 MAPK (Thr180) Antibody	38211409	2.6	WB	Rat

269.Spirulina platensis as a growth booster for broiler; Insights into their nutritional, molecular, immunohistopathological, and microbiota modulating effects

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-2072R	polyclonal rabbit anti-iNOS Antibody	38183085	2.6	IHC	Chicken

270.Role of cholinergic innervation in biliary remnants of patients with biliary atresia

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
PV-0024	SP-0024 Histostain-Plus Kit	38259596	2.6	IHC	Human

271.Tissue and Cell Effects of probiotic supplementation on chronic inflammatory process modulation in colorectal carcinogenesis

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-6309R	polyclonal anti-IL-6 Antibody	10.1016/j.tice.2023.102293	2.6	IHC	Mouse

272.Study on the mechanism of action of Wu Mei Pill in inhibiting rheumatoid arthritis through TLR4-NF- κ B pathway

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1287R	anti-I κ B- α antibody	38218891	2.6	WB	Rat
bs-2513R	anti-p-I κ B- α antibody	38218891	2.6	WB	Rat
bs-20705R	anti-MMP-2 antibody	38218891	2.6	WB	Rat
bs-7059R	anti-MMP-9 antibody	38218891	2.6	WB	Rat
bs-43009R	anti-TIMP-1 antibody	38218891	2.6	WB	Rat
bs-0416R	anti-TIMP-2 antibody	38218891	2.6	WB	Rat
bs-1840R	anti-Lamin B antibody	38218891	2.6	WB	Rat

273.Aseptic loosening around total joint replacement in humans is regulated by miR-1246 and miR-6089 via the Wnt signalling pathway

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-12435R	Rabbit Polyclonal anti-Human FRAT2 IgG antibody	38287447	2.6	WB	Human

274.Complex Topology of Ubiquitin Chains Mediates Lysosomal Degradation of MrgC Proteins

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0061R	Rabbit anti- β -Actin antibody	38291169	2.6	WB	Mouse

275.ECM1-associated miR-1260b promotes osteogenic differentiation by targeting GDI1

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
------	-------------	-----------	--------------	--------------	-----------------

bs-7158R	anti-COL1A1 antibody	38266317	2.5	WB	Human
----------	----------------------	----------	-----	----	-------

276.Ser252Trp mutation in fibroblast growth factor receptor 2 promotes branching morphogenesis in mouse salivary glands

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1374R	anti-bone morphogenetic protein (BMP) 4 antibody	38246420	2.4	IHC	Mouse

277.Molecular Network Mechanism Analysis of Urine Stem Cells Against Retinal Aging

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
BA00208	Cell counting Kit-8	38273154	2.4	Other	Other

278.Hesperidin Protects Against High-Fat Diet-Induced Lipotoxicity in Rats by Inhibiting Pyroptosis

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0755R	anti-GAPDH antibody	38294790	2.4	WB	Human

279.Norepinephrine protects against cochlear outer hair cell damage and noise-induced hearing loss via α 2A-adrenergic receptor

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0498R	anti- β 1 antibody	38291397	2.4	WB	Rat
bs-1063R	anti- β 3 antibody	38291397	2.4	WB	Rat

280.The effects of low ambient temperature on steroidogenesis and mitochondrial functions in the testes of wild ground squirrels (*Spermophilus dauricus*)

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bsm-52289R	rabbit monoclonal anti-Cleaved Caspase-3 antibody	38228267	2.3	IHC	Squirrels
bs-0081R	rabbit polyclonal anti-Caspase-3 antibody	38228267	2.3	IHC	Squirrels
bs-23683R	rabbit polyclonal anti-voltage-dependent anion channel (VDAC) antibody	38228267	2.3	IHC	Squirrels
bs-6431R	rabbit polyclonal anti-LHR antibody	38228267	2.3	IHC	Squirrels
bs-0895R	anti-FSHR antibody	38228267	2.3	IHC	Squirrels

bs-23982R	anti-GATA4 antibody	38228267	2.3	IHC	Squirrels
bs-20388R	anti-StAR antibody	38228267	2.3	IHC	Squirrels
bs-3608R	anti-P450scc antibody	38228267	2.3	IHC	Squirrels
bs-3853R	anti-P450c17 antibody	38228267	2.3	IHC	Squirrels
bs-0114M	mouse polyclonal anti-P450arom antibody	38228267	2.3	IHC	Squirrels
bs-23685R	anti-MFN2 antibody	38228267	2.3	IHC	Squirrels

281.Cytosolic mtDNA–cGAS–STING axis contributes to sepsis-induced acute kidney injury via activating the NLRP3 inflammasome

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
js-0368G-AF488	Goat Anti-Mouse IgM / AF488	38238499	2.3	IF	Human

282.Rosacea treatment with mussel adhesive protein delivered via microneedling: In vivo and clinical studies

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-23926R	Anti-TRPV1 antibody	38284129	2.3	IHC	Mouse

283.Modified Buyang Huanwu Decoction alleviates diabetic liver injury via inhibiting oxidative stress in db/db mice

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-5661R	mouse anti-NFκB-P65 antibody	38322549	2.2	WB	Mouse
bs-5514R	rabbit anti-p-IκBα antibody	38322549	2.2	WB	Mouse

284.Immunofluorescence analysis of human eosinophils

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
js-1782R-AF488	Alexa Fluor 488 anti-transferrin receptor (TfnRc) polyclonal antibody	38272178	2.2	IF	Human
js-0295P-AF488	isotype rabbit IgG conjugated to Alexa Fluor 488	38272178	2.2	IF	Human

285.Identification of uromodulin deposition in the stroma of perinephric fibromyxoid nephrogenic adenoma by mass spectrometry

SKU	Product Name	Pubmed ID	Journal Tier	Applications	Cross Reactivity
bs-41397R	rabbit polyclonal anti-uromodulin antibody	38289139	2.2	IHC	Human

286. Colchicine ameliorates myocardial injury induced by coronary microembolization through suppressing pyroptosis via the AMPK/SIRT1/NLRP3 signaling pathway

SKU	Product Name	Pubmed ID	Journal Tier	Applications	Cross Reactivity
bs-10442R	anti-caspase-1 p20 antibody	38172692	2.1	IF, WB	Rat

287. Coenzyme Q10 mitigates macrophage mediated inflammation in heart following myocardial infarction via the NLRP3/IL1 β pathway

SKU	Product Name	Pubmed ID	Journal Tier	Applications	Cross Reactivity
bs-6319R	rabbit anti-IL-1 β polyclonal antibody	38281937	2.1	IF, WB	Mouse

288. Loureirin A Promotes Cell Differentiation and Suppresses Migration and Invasion of Melanoma Cells via WNT and AKT/mTOR signaling Pathways

SKU	Product Name	Pubmed ID	Journal Tier	Applications	Cross Reactivity
bs-15510R	Anti-TRP1 antibody	38199251	2	WB	Mouse, Human

289. Expression of the SARS-CoV-2 Receptor ACE2 and Protease TMPRSS2 in Ocular Hypertension Eyes of Nonhuman Primate and Human

SKU	Product Name	Pubmed ID	Journal Tier	Applications	Cross Reactivity
bs-23443R	rabbit anti-ACE2 antibody	38212998	2	IF	Monkey

290. Icaritin regulates EMT and stem cell-like character in breast cancer through modulating lncRNA NEAT1/TGF β /SMAD2 signaling pathway

SKU	Product Name	Pubmed ID	Journal Tier	Applications	Cross Reactivity
bs-11208R	Anti-SOX4 antibody	38220208	2	WB	Human

291. The Mechanism of 5-Fluorouracil-Induced Hyperpigmentation in HRM-2 Hairless Mice: Focus on the Increase of Blood Vessels

SKU	Product Name	Pubmed ID	Journal Tier	Applications	Cross Reactivity
bs-0819R	Anti-tyrosinase antibody	38281766	2	IHC	Mouse

bs-4754R	Anti-vWF antibody	38281766	2	IHC	Mouse
----------	-------------------	----------	---	-----	-------

292.Evaluation of the Immune Effect of a Triple Vaccine Composed of Fowl Adenovirus Serotype 4 Fibre-2 Recombinant Subunit, Inactivated Avian Influenza (H9N2) Vaccine and Newcastle Disease Vaccine Against Respective Pathogenic Virus Challenge in Chickens

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0432R-HRP	HRP-conjugated rabbit anti-chicken IgY antibody	10.1016/j.japr.2024.100410	1.9	ELISA	Chicken

293.B2M is a Biomarker Associated With Immune Infiltration In High Altitude Pulmonary Edema

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-2072R	Anti-Inducible nitric oxide synthase (NOS2) antibody	37165489	1.8	IF	Rat

294.Rreb1 is a key transcription factor in Sertoli cell maturation and function and spermatogenesis in mouse

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-18856R	Anti-RREB1 antibody	38248872	1.7	IF	Mouse

295.Ethanol extract of Abelmoschus manihot suppresses endoplasmic reticulum stress in contrast-induced nephropathy

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
C-0003	DAB solution	10.4103/2221-1691.393577	1.7	Other	Other

296.Safety of Three-Dimensional Conformal Radiotherapy on the Sheep with Pulmonary Cystic Echinococcosis: A Preliminary Study

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0086R	anti-TGF- β antibody	38233676	1.5	IHC	Sheep
bs-4605R	anti-MMP2 antibody	38233676	1.5	IHC	Sheep
bs-4593R	anti-MMP9 antibody	38233676	1.5	IHC	Sheep

297.Ivabradine Alleviates Experimental Autoimmune Myocarditis-Mediated Myocardial Injury

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-20619R	anti-MMP-9 antibody	38296563	1.5	WB	Rat

298. Draining Lymph Node Metastasis Model for Assessing the Dynamics of Antigen-Specific CD8

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-0292P	BSA-V (bovine serum albumin)	10.3791/65646	1.2	Other	Other

299. Chikungunya virus infection IVa Protects against H9N2 Avian Influenza Virus Infection by Inhibiting Inflammation and TLR4/NF- κ B Signaling Pathways in vitro and in vivo

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-0755R	Anti-GAPDH antibody	10.1177/09731296231216165	0.7	WB	Human, Mouse
bs-20594R	Anti-TLR4 antibody	10.1177/09731296231216165	0.7	WB, IHC	Human, Mouse
bs-1047R	Anti-MyD88 antibody	10.1177/09731296231216165	0.7	WB	Human, Mouse
bs-20160R	Anti-NF κ B p65 antibody	10.1177/09731296231216165	0.7	WB, IHC, IF	Human, Mouse
bs-5661R	Anti-phospho-NF κ B p65 (Thr435) antibody	10.1177/09731296231216165	0.7	WB	Human, Mouse
bs-0812R	Anti-IL-1 Beta antibody	10.1177/09731296231216165	0.7	IHC	Mouse
bs-0782R	Anti-IL-6 antibody	10.1177/09731296231216165	0.7	IHC	Mouse
bs-0780R	Anti-IL-8 antibody	10.1177/09731296231216165	0.7	IHC	Mouse