

2024 Feb-Bioss 产品文献引用精选

1.A novel peptide PDHK1-241aa encoded by circPDHK1 promotes ccRCC progression via interacting with PPP1CA to inhibit AKT dephosphorylation and activate the AKT-mTOR signaling pathway

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bsm-33070M	anti-Ki-67 antibody	38360682	37.3	IHC	Mouse

2.Multi-omic profiling of clear cell renal cell carcinoma identifies metabolic reprogramming associated with disease progression

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1050R	anti-human CX40 antibody	38361033	30.8	IHC	Human

3.Molecular Mechanisms of Intracellular Delivery of Nanoparticles Monitored by an Enzyme-Induced Proximity Labeling

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-13999R	anti-CORO7 antibody	38300384	26.6	IF	Mouse
bs-13490R	anti-GOLIM4 antibody	38300384	26.6	IF	Mouse

4.CD4+ T cells produce IFN-I to license cDC1s for induction of cytotoxic T-cell activity in human tumors

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-2789R	rabbit anti-human TAP1 polyclonal antibody	38383773	24.1	FC	Human
bs-2374R	rabbit anti-human TAP2 polyclonal antibody	38383773	24.1	FC	Human

5.Manipulation of Surface Electrical Charge on Nanocomposite Membranes Confers Wide Spectrum Bactericidal Effects and Promotes Tissue Regeneration

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-4917R	anti-OCN antibody	10.1002/adfm.202314024	19	IF	Mouse

6.Controlled-release hydrogel loaded with magnesium-based nanoflowers synergize immunomodulation and cartilage regeneration in tendon-bone healing

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bsm-33409M	anti-Col II antibody	38440323	18.9	IF	Rat
bs-1223R	anti-Aggrecan antibody	38440323	18.9	IF	Rat

7.Polydopamine-modified metal-organic frameworks nanoparticles enhance the corrosion resistance and bioactivity of polycaprolactone coating on high-purity magnesium

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
-----	-------------	-----------	--------------	--------------	-----------------

bs-0195R	anti-CD31 antibody	10.1016/j.jma.2023.12.014	17.6	IHC	Rabbit
bs-20624R	anti-OPG antibody	10.1016/j.jma.2023.12.014	17.6	IHC	Rabbit
bs-4917R	anti-OCN antibody	10.1016/j.jma.2023.12.014	17.6	IHC	Rabbit

8.Reinforced Immunogenic Endoplasmic Reticulum Stress and Oxidative Stress via an Orchestrated Nanophotoinducer to Boost Cancer Photoimmunotherapy

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bsR-0292P	Rhodamine B isothiocyanate (RBITC)-labeled BSA	38382065	17.1	Other	Other
D-9106	4,6-diamidino-2-phenylindole (DAPI)	38382065	17.1	Other	Other
bsm-33070M	anti-Ki67 antibody	38382065	17.1	IHC	Mouse

9.Neutrophil Nanodecoys Inhibit Tumor Metastasis by Blocking the Interaction between Tumor Cells and Neutrophils

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1629R	Rabbit anti-CXCR2 antibody	38422392	17.1	WB	Mouse

10.UPP1 promotes lung adenocarcinoma progression through the induction of an immunosuppressive microenvironment

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bsm-33187m	anti- α -SMA antibody	38331898	16.6	IF	Human

11.Molecular co-assembled strategy tuning protein conformation for cartilage regeneration

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-10725R	anti-Sox9 antibody	38374253	16.6	IF	Rabbit
bs-0296R-Cy7	Rabbit Anti-Mouse IgG H&L / Cy7 antibody	38374253	16.6	IF	Mouse

12.Selective Suppression of Integrin-Ligand Binding by Single Molecular Tension Probes Mediates Directional Cell Migration

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-3159R-AF594	anti-phosphorylated FAK in Tyr397 antibody	38311584	15.1	IF	Mouse

13.Hierarchically Assembled Nanofiber Scaffolds with Dual Growth Factor Gradients Promote Skin Wound Healing Through Rapid Cell Recruitment

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
-----	-------------	-----------	--------------	--------------	-----------------

bs-0292P-FITC	FITC-BSA	38326085	15.1	Other	Other
bsm-60235R	Anti-K6 antibody	38326085	15.1	IHC	Mouse

14. Dry powder microneedle-enabled transdermal anti-inflammatory therapy for obesity, diabetes, hyperlipidemia, and fatty liver

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1278R	anti-8-OHdG antibody	10.1016/j.cej.2024.149395	15.1	IF	Mouse

15. 3D printing of personalized magnesium composite bone tissue engineering scaffold for bone and angiogenesis regeneration

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-20735R	anti-FAK antibody	10.1016/j.cej.2024.149444	15.1	WB	Human
bs-3159R	anti-p-FAK antibody	10.1016/j.cej.2024.149444	15.1	WB	Human

16. Myofibrillogenesis Regulator-1 Regulates the Ubiquitin Lysosomal Pathway of Notch3 Intracellular Domain Through E3 Ubiquitin-Protein Ligase Itchy Homolog in the Metastasis of Non-Small Cell Lung Cancer

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-8713R	Anti-Itchy homolog (ITCH) antibody	38342606	15.1	CoIP, WB	Human

17. A collective hyperthermia-facilitated glutathione inhibition strategy for boosted oxidative stress treatment

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-2184R	Rabbit anti-NQO1 antibody	10.1016/j.cej.2024.150056	15.1	WB	Mouse
bsm-33284M	anti-caspase-3 antibody	10.1016/j.cej.2024.150056	15.1	WB	Mouse

18. Multivalent Aptamer-Based Lysosome-Targeting Chimeras (LYTACs) Platform for Mono- or Dual-Targeted Proteins Degradation on Cell Surface

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0437P-FITC	FITC-streptavidin (SA)	38425146	15.1	Other	Other

19. The osteoclastic activity in apical distal region of molar mesial roots affects orthodontic tooth movement and root resorption in rats

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-10802R	anti-TNF- α antibody	38418457	14.9	IHC	Rat

20. Tracking endogenous proteins based on RNA editing-mediated genetic code expansion

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
-----	-------------	-----------	--------------	--------------	-----------------

bs-3930R	rabbit anti-Cav2.1 (CACNA1A) antibody	38302606	14.8	ICC	Mouse
----------	---------------------------------------	----------	------	-----	-------

21. Adipose triglyceride lipase suppresses noncanonical inflammasome by hydrolyzing LPS

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0311R	Anti-biotin antibody	38413746	14.8	WB,IP	Human,Mouse
bs-0294D-AF488	Alexa Fluor Plus 488-conjugated donkey anti-goat IgG antibody	38413746	14.8	IF	Human,Mouse
bs-0295D-AF594	Alexa Fluor 594-conjugated donkey anti-rabbit IgG antibody	38413746	14.8	IF	Human,Mouse

22. Deep simulated annealing for the discovery of novel dental anesthetics with local anesthesia and anti-inflammatory properties

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1278R	anti-8-ohdG antibody	10.1016/j.apsb.2024.01.019	14.5	IF	Rat
bs-0295G-AF647	goat anti-rabbit IgG H&L/AF647 antibody	10.1016/j.apsb.2024.01.019	14.5	IF	Rat
bs-0296G-AF488	goat anti-mouse IgG H&L/AF488 antibody	10.1016/j.apsb.2024.01.019	14.5	IF	Rat

23. Integrated Single-Cell Transcriptomic Atlas of Human Kidney Endothelial Cells

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-23861R	anti-RGCC antibody	38351505	13.6	IF	Human

24. Alginate di-aldehyde modified metal-organic framework nanocarriers as delivery and adjuvant in inactivated pseudorabies vaccination

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
C5085	alum adjuvant	38376908	13.3	Other	Other

25. Sodium butyrate ameliorates high glucose-suppressed neuronal mitophagy by restoring PRKN expression via inhibiting the RELA-HDAC8 complex

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1278R	anti-8-hydroxy-2-deoxyguanosine antibody	38409852	13.3	IF	Mouse
bs-6313R	anti-4-hydroxy-2-noneal antibody	38409852	13.3	IF	Mouse

26. Interorgan communication in neurogenic heterotopic ossification: the role of brain-derived extracellular vesicles

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1996R	anti-L1CAM antibody	38383487	12.7	IHC,WB	Rat
bs-1996R-Bio	anti-L1CAM biotinylated antibody	38383487	12.7	FC	Rat

27. A DNA tetrahedron-based ferroptosis-suppressing nanoparticle: superior delivery of curcumin and alleviation of diabetic osteoporosis

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1158P	advanced glycation end products (AGEs)	38424439	12.7	Other	Other

28. Hnrnpk Protects Against Osteoarthritis through Targeting WWC1 mRNA and Inhibiting Hippo Signaling Pathway

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-3605R	anti-YAP antibody	38414246	12.4	IHC	Mouse

29. Gut microbiota modulate CD8+ T cell immunity in gastric cancer through Butyrate/GPR109A/HOPX

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0480R	anti-IFN- γ antibody	38319728	12.2	WB,IHC,IF	Mouse
bs-1035R	anti-CD86 antibody	38319728	12.2	IF	Mouse

30. Astrocyte-induced mGluR1 activates human lung cancer brain metastasis via glutamate-dependent stabilization of EGFR

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-2590R	Rabbit anti-NRSF/REST polyclonal antibody	38309264	11.8	ICC	Human

31. A tumor suppressor protein encoded by circKEAP1 inhibits osteosarcoma cell stemness and metastasis by promoting vimentin proteasome degradation and activating anti-tumor immunity

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
BA00103	Annexin V-AF647 and PI Apoptosis Detection kit	38383479	11.3	Other	Other

32. An antibacterial biologic patch based on bacterial cellulose for repair of infected hernias

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-6433R	Rabbit Anti-Cathelicidin/Camp antibody	38494213	11.2	IF	Mouse

33.mGluR7 allosteric modulator AMN082 corrects protein synthesis and pathological phenotypes in FXS

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-11085R	rabbit anti-PCDH7 antibody	38374465	11.1	WB	Mouse

34.Enhanced photodynamic therapy efficacy of Ni-doped/oxygen vacancy double-defect Ni-ZnO@C photosensitizer in bacteria-infected wounds based on ROS damage and ATP synthesis inhibition

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
AK300	NADH detection kit	10.1016/j.jmst.2024.01.018	10.9	Other	Other

35.Porous Ti3SiC2 ceramics with improved osteogenic functions via biomineralization as load-bearing bone implants

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-2072R	rabbit anti-iNOS polyclonal antibody	10.1016/j.jmst.2024.01.025	10.9	IF	Mouse

36.Adenosine-modulating synthetic high-density lipoprotein for chemoimmunotherapy of triple-negative breast cancer

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-5900R-PE	PE-anti-ADORA2B antibody	38295994	10.8	FC	Mouse
bs-0189R-AF488	AF488-anti- α -SMA antibody	38295994	10.8	FC	Mouse

37.SiATG5-loaded cancer cell membrane-fused liposomes induced increased uptake of albumin-bound chemotherapeutics by pancreatic cancer cells

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0945P-FITC	Human Serum Albumin-FITC Conjugate (FITC-HSA)	38311244	10.8	Other	Other

38.BPA and low-Se exacerbate apoptosis and mitophagy in chicken pancreatic cells by regulating the PTEN/PI3K/AKT/mTOR pathway

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-2951R	anti-p62 antibody	38311007	10.7	IF	Chicken

39.Ginseng-derived nanoparticles alleviate inflammatory bowel disease via the TLR4/MAPK and p62/Nrf2/Keap1 pathways

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bsm-33036M	anti- β -actin antibody	38302938	10.2	WB	Mouse

40.Spatiotemporal transformable nano-assembly for on-demand drug delivery to enhance anti-tumor immunotherapy

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bsm-33400M-FITC	FITC conjugated anti-collagen I antibody	38434719	10.2	IF	Mouse
bsm-33187M-AF594	AF594 conjugated anti- α -SMA antibody	38434719	10.2	IF	Mouse
bs-0396R-AF488	Alexa Fluor [®] 488 conjugated anti-mouse VCAM-1 antibody	38434719	10.2	IF	Mouse
bs-5758R-FITC	FITC conjugated anti-mouse FAP- α antibody	38434719	10.2	IF	Mouse
bs-0086R	anti-TGF- β antibody	38434719	10.2	WB	Mouse
bs-10196R	anti- α -SMA antibody	38434719	10.2	WB	Mouse

41.Turning Threat to Therapy: A Nanozyme-Patch in Surgical Bed for Convenient Tumor Vaccination by Sustained In Situ Catalysis

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-20633R	Rabbit anti-HMGB-1 polyclonal antibody	38301259	10	IF	Mouse

42.Mn-Based Artificial Mitochondrial Complex "VI" Acts as an Electron and Free Radical Conversion Factory to Suppress Macrophage Inflammatory Response

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-40295G-HRP	goat anti-rabbit IgG H&L antibody	38301194	10	IHC	Mouse

43.Platinum-loaded Cerium Oxide Capable of Repairing Neuronal Homeostasis for Cerebral Ischemia-Reperfusion Injury Therapy

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0782R	anti-IL-6 antibody	38323853	10	WB	Mouse
bsm-33287M	anti-TGF- β antibody	38323853	10	WB	Mouse

44.3D Printing of Cobalt-Incorporated Chlorapatite Bioceramic Composite Scaffolds with Antioxidative Activity for Enhanced Osteochondral Regeneration

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0698R	anti-IL-10 antibody	38363057	10	IF	Mouse
bs-2150R	anti-THF- α antibody	38363057	10	IF	Mouse

45.Environmental dose of 16 priority-controlled PAHs induce endothelial dysfunction: An in vivo and in vitro study

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
-----	-------------	-----------	--------------	--------------	-----------------

bs-3447R	anti-phosphorylated eNOS (Ser1177) antibody	38340817	9.8	WB	Rat,Human
----------	---------------------------------------------	----------	-----	----	-----------

46. Acidosis activates breast cancer ferroptosis through ZFAND5/SLC3A2 signaling axis and elicits M1 macrophage polarization

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-10142R	anti-ZFAND5 antibody	38360142	9.7	WB	Human

47. Silk Fibroin-Coated Nano-MOFs Enhance the Thermal Stability and Immunogenicity of HBsAg

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
BA00121	Cy5.5 antibody conjugation kit	38323561	9.5	Other	Other
C5084	Rehydragel@LV alum adjuvant	38323561	9.5	Other	Other

48. Nanoplasmonic Detection of EGFR Mutations Based on Extracellular Vesicle-Derived EGFR-Drug Interaction

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1523R	anti-CD63 antibody	38335730	9.5	WB	Human
bsm-51699M	anti-TSG101 antibody	38335730	9.5	WB	Human
bs-1693R	anti-Calnexin antibody	38335730	9.5	WB	Human

49. A Multifunctional, Tough, Stretchable, and Transparent Curcumin Hydrogel with Potent Antimicrobial, Antioxidative, Anti-inflammatory, and Angiogenesis Capabilities for Diabetic Wound Healing

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-6319R	anti-Interleukin-1 β (IL-1 β) antibody	38359334	9.5	IHC	Rat
bs-0581R	anti-Interleukin-4 (IL-4) antibody	38359334	9.5	IF	Mouse
bs-1035R	anti-cluster of differentiation 86 (CD86) antibody	38359334	9.5	IF	Mouse
bsm-52791R	anti-cluster of differentiation 206 (CD206) antibody	38359334	9.5	IF	Mouse
bs-0195R	anti-cluster of differentiation 31 (CD31) antibody	38359334	9.5	IHC	Rat

50. A novel mitochondria-targeting DHODH inhibitor induces robust ferroptosis and alleviates immune suppression

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-16882R	anti-LACTB antibody	38423231	9.3	WB	Human,Mouse

51. Deubiquitinase Mym1 regulates neural stem cell proliferation and differentiation by controlling Id4 expression

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-0199R	anti-Gfap antibody	38342917	9	WB	Mouse
bsm-33199M	anti-Active Caspase-3 antibody	38342917	9	IF	Mouse

52. CHOP-mediated Gasdermin E expression promotes pyroptosis, inflammation, and mitochondrial damage in renal ischemia-reperfusion injury

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-0295G-RBITC	goat anti-rabbit IgG-rhodamine antibody	38388468	9	IF	Mouse
bs-0296G-RBITC	goat anti-mouse IgG-rhodamine antibody	38388468	9	IF	Mouse

53. Antibacterial, ROS scavenging and angiogenesis promoting ε-polylysine/gelatin based hydrogel containing CTLP to regulate macrophages for pressure ulcer healing

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-1035R	Rabbit Anti-CD86 antibody	38408382	9	WB	Mouse
bs-0468R	Rabbit Anti-CD31 antibody	38408382	9	WB	Mouse
bs-5193R	Anti-phospho-AKT1 + AKT2 + AKT3 (Thr308+Thr309+Thr305) antibody	38408382	9	WB	Mouse

54. Phthalate drives splenic inflammatory response via activating HSP60/TLR4/NLRP3 signaling axis-dependent pyroptosis

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-2513R	anti-p-IκB Polyclonal antibody	38382728	8.9	WB	Mouse

55. Paeoniflorin ameliorates chronic colitis via the DR3 signaling pathway in group 3 innate lymphoid cells

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bsm-33036M	anti-β-actin antibody	10.1016/j.jpha.2024.01.008	8.8	WB	Mouse, Human
bs-41373R	anti-GAPDH antibody	10.1016/j.jpha.2024.01.008	8.8	WB	Mouse

56. Long-term activation of anti-tumor immunity in pancreatic cancer by a p53-expressing telomerase-specific oncolytic adenovirus

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-2499R	anti-CD69 antibody	38316993	8.8	IHC	Human

57. Carbon dots and covalent organic frameworks based FRET immunosensor for sensitive detection of Escherichia coli O157:H7

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1563R	anti-E. coli O157: H7 antibody	38489878	8.8	Other	Other
C5047	20 × Tris Buffered Saline (TBS buffer)	38489878	8.8	Other	Other

58. miR615-3p inhibited FBLN1 and osteogenic differentiation of umbilical cord mesenchymal stem cells by associated with YTHDF2 in a m6A-miRNA interaction manner

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-10316R	anti-DSPP antibody	38353178	8.5	WB,IF	Human
bs-12359R	anti-DMP1 antibody	38353178	8.5	WB,IF	Human
bs-0470R	anti-OCN antibody	38353178	8.5	IF	Human
bs-0295P	Rabbit IgG	38353178	8.5	Other	Other

59. T-type voltage-gated channels, Na⁺/Ca²⁺-exchanger, and calpain-2 promote photoreceptor cell death in inherited retinal degeneration

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1997R	rabbit anti-NCX2 antibody	38303059	8.4	IHC	Mouse

60. Musashi-2 potentiates colorectal cancer immune infiltration by regulating the post-translational modifications of HMGB1 to promote DCs maturation and migration

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-8000P	LPS	38347600	8.4	Other	Other

61. Oral delivery of miR-146a-5p overexpression plasmid-loaded Pickering double emulsion modulates intestinal inflammation and the gut microbe

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1329R	anti-ZO-1 antibody	38307433	8.2	WB	Mouse
bs-0465R	anti-p65 antibody	38307433	8.2	WB	Mouse

62. Downregulation of VEGFA accelerates AGEs-mediated nucleus pulposus degeneration through inhibiting protective mitophagy in high glucose environments

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1158P	advanced glycation end products (AGEs)	38320636	8.2	Other	Other

63. Whole-transcriptome analysis of longissimus dorsi muscle in cattle-yaks reveals the regulatory functions of ADAMTS6 gene in myoblasts

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-8009R	rabbit Anti-ADAMTS6 antibody	38342263	8.2	WB	Mouse
bs-9610R	rabbit Anti-MyHC antibody	38342263	8.2	WB	Mouse
bsm-52028M	mouse Anti-Cdk4 antibody	38342263	8.2	WB	Mouse
bs-2006R	rabbit Anti-PCNA antibody	38342263	8.2	WB	Mouse
bs-0061R	rabbit Anti-beta-Actin antibody	38342263	8.2	WB	Mouse

64.METTL3 enhances E. coli F18 resistance by targeting IKBKG/NF- κ B signaling via an m6A-YTHDF1-dependent manner in IPEC-J2 cells

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-17608R	rabbit anti-METTL14 antibody	38346619	8.2	WB	Pig
bs-22575R	rabbit anti-IKBKG antibody	38346619	8.2	WB	Porcine

65.Characterization of a novel γ -type carbonic anhydrase, Sjy-CA2, in Saccharina japonica: Insights into carbon concentration mechanism in macroalgae

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bsm-33226M	anti-RbcL antibody	38423426	8.2	WB	Saccharina japonica
bs-0715R	anti- β -tubulin antibody	38423426	8.2	WB	Saccharina japonica

66.Antler thymosin β 10 reduces liver fibrosis via inhibiting TGF- β 1/SMAD pathway

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bsm-33345M	anti-TGF- β 1 antibody	38428779	8.2	IHC,WB	Rat
bs-0638R	anti-TGF- β R1 antibody	38428779	8.2	WB	Rat
bs-0117R	anti-TGF- β R2 antibody	38428779	8.2	WB	Rat

67.Artesunate-Nanoliposome-TPP, a Novel Drug Delivery System That Targets the Mitochondria, Attenuates Cisplatin-Induced Acute Kidney Injury by Suppressing Oxidative Stress and Inflammatory Effects

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-4999R	anti-RAGE antibody	38371457	8	IF,WB	Mouse

68.Hydrogels Containing Chitosan-Modified Gold Nanoparticles Show Significant Efficacy in Healing Diabetic Wounds Infected with Antibiotic-Resistant Bacteria

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-6761R	anti-interleukin-10 (IL-10) antibody	38406603	8	IF	Rat

bsm-33207M	anti-tumor necrosis factor- α (TNF- α) antibody	38406603	8	IF	Rat
bs-0468R	anti-cluster of differentiation 31 (CD31) antibody	38406603	8	IF	Rat

69.BOP1 contributes to the activation of autophagy in polycystic ovary syndrome via nucleolar stress response

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-12877R	Rabbit anti-BOP1 antibody	38409361	8	WB	Mouse

70.Effects of Folate-Chicory Acid Liposome on Macrophage Polarization and TLR4/NF- κ B signaling pathway in ulcerative colitis mouse

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-7058R-AF488	anti-AF488-F4/80 antibody	38503151	7.9	FC	Mouse
bs-4727R-APC	anti-APC-CD206 antibody	38503151	7.9	FC	Mouse
bs-1014R-PE-Cy5.5	anti-PerCP-Cy5.5-CD11b antibody	38503151	7.9	FC	Mouse
bs-1035R-PE	anti-PE-CD86 antibody	38503151	7.9	FC	Mouse

71.Dioscin decreases M2 polarization via inhibiting a positive feedback loop between RBM47 and NF- κ B in glioma

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1162R	anti-CXCL2 antibody	10.1016/j.phymed.2024.155417	7.9	WB	Human

72.Rutin ameliorated lipid metabolism dysfunction of diabetic NAFLD via AMPK/SREBP1 pathway

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0351R	anti-GLUT2 antibody	38394735	7.9	IF	Human

73.Rosmarinic acid mitigates acrylamide induced neurotoxicity via suppressing endoplasmic reticulum stress and inflammation in mouse hippocampus

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-6309R	anti-IL-6 antibody	38394736	7.9	WB	Mouse
bs-1219R	anti-Bip antibody	38394736	7.9	WB	Mouse

74.Andrographolide sensitizes KRAS-mutant colorectal cancer cells to cetuximab by inhibiting the EGFR/AKT and PDGFR β /AKT signaling pathways

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
------	-------------	-----------	--------------	--------------	-----------------

C01-03001	10% normal goat serum	38394734	7.9	Other	Other
-----------	-----------------------	----------	-----	-------	-------

75.Paris polyphylla ethanol extract and polyphyllin I ameliorate adenomyosis by inhibiting epithelial–mesenchymal transition

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-10048R	rabbit anti-Von Willebrand factor polyclonal antibody	38452697	7.9	IF	Human

76.The mechanism of curcumin to protect mouse ovaries from oxidative damage by regulating AMPK/mTOR mediated autophagy

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
D10393	hematoxylin dyeing solution	38471315	7.9	Other	Other
S0159	eosin	38471315	7.9	Other	Other
bs-1353R	anti-BECN1 antibody	38471315	7.9	WB	Mouse,Human
bs-10216R	anti-SOD1 antibody	38471315	7.9	WB	Mouse
bs-0061R	anti-β-actin antibody	38471315	7.9	WB	Mouse,Human

77.The Sanbi Decoction alleviates intervertebral disc degeneration in rats through intestinal flora and serum metabolic homeostasis modulation

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-10581R	anti-MMP13 antibody	38484462	7.9	IHC,WB	Rat
bs-4191R	anti-ADAMTS4 antibody	38484462	7.9	WB	Rat

78.Ginseng-derived nanoparticles alleviate alcohol-induced liver injury by activating the Nrf2/HO-1 signalling pathway and inhibiting the NF-κB signalling pathway in vitro and in vivo

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0061R	anti-β-actin antibody	38458086	7.9	WB	Mouse,Human

79.HER2 phosphorylation induced by TGF-β promotes mammary morphogenesis and breast cancer progression

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bsm-33346M	anti-Flag antibody	38407425	7.8	WB	Human

80.PSTPIP2 ameliorates aristolochic acid nephropathy by suppressing interleukin-19-mediated neutrophil extracellular trap formation

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-2713R	anti-kidney injury molecule-1 (KIM-1) antibody	38314821	7.7	IHC,IF,WB	Mouse
bs-19580R	anti-PSTPIP2 antibody	38314821	7.7	IF,WB,CoIP	Mouse
bs-0061R	anti-β-actin antibody	38314821	7.7	WB	Mouse

81.ADSC-Exos enhance functional recovery after spinal cord injury by inhibiting ferroptosis and promoting the survival and function of endothelial cells through the NRF2/SLC7A11/GPX4 pathway

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-6883R	anti-SLC7A11 antibody	38306845	7.5	WB	Rat,Human
bs-1074R	anti-NRF2 antibody	38306845	7.5	WB	Rat,Human
bs-6313R	anti-4-HNE antibody	38306845	7.5	WB	Rat,Human

82.3D printed PLGA scaffold with nano-hydroxyapatite carrying linezolid for treatment of infected bone defects

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-10423R	anti-Collagen I antibody	38320333	7.5	IHC	Rabbit

83.Serine protease inhibitor from the muscle larval Trichinella spiralis ameliorates non-alcoholic fatty liver disease in mice via anti-inflammatory properties and gut-liver crosstalk

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-10196R	rabbit anti-mouse α -SMA antibody	38325266	7.5	IHC	Mouse
bs-10802R	rabbit anti-mouse TNF- α antibody	38325266	7.5	IHC	Mouse
bs-6319R	rabbit anti-mouse IL-1 β antibody	38325266	7.5	IHC	Mouse

84.Human amniotic mesenchymal stromal cell-derived exosomes promote neuronal function by inhibiting excessive apoptosis in a hypoxia/ischemia-induced cerebral palsy model: A preclinical study

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-34012R	anti-Bcl2 antibody	38394849	7.5	WB	Rat

85.Koumine inhibits IL-1 β -induced chondrocyte inflammation and ameliorates extracellular matrix degradation in osteoarthritic cartilage through activation of PINK1/Parkin-mediated mitochondrial autophagy

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-3573R	anti-Adamts5 antibody	38412715	7.5	WB	Rat

86.The inhibition of YAP Signaling Prevents Chronic Biliary Fibrosis in the Abcb4-/- Model by Modulation of Hepatic Stellate Cell and Bile Duct Epithelium Cell Pathophysiology

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-5791R	anti-Integrin α V β 6 antibody	37307826	7.4	WB	Human

87.ZSWIM4 inhibition improves chemosensitivity in epithelial ovarian cancer cells by suppressing intracellular glycine biosynthesis

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
IHC001	IHC kit	38383406	7.4	Other	Other
D-9106	DAPI	38383406	7.4	Other	Other

88. Microbiota modulates the steroid response to acute immune stress in male mice

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-3898R	anti-CYP11B1 polyclonal antibody	38361932	7.3	WB	Mouse

89. Myeloid-derived suppressor cell-derived osteoclasts with bone resorption capacity in the joints of arthritic SKG mice Authors

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-10037R-PE-Cy7	PE-Cy7 conjugated anti NFκB2 p52 (polyclonal) antibody	10.3389/fimmu.2024.1168323	7.3	FC	Mouse
bs-3562R-PE	PE conjugated anti Relb (polyclonal) antibody	10.3389/fimmu.2024.1168323	7.3	FC	Mouse

90. Urinary miR-185-5p is a biomarker of renal tubulointerstitial fibrosis in IgA nephropathy

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bsm-0978M	anti-GAPDH antibody	38426107	7.3	WB	Mouse, Human

91. Integrative analysis identifies oxidative stress biomarkers in non-alcoholic fatty liver disease via machine learning and weighted gene co-expression network analysis

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-0742R	anti-CDKN1B antibody	38476236	7.3	WB	Mouse

92. D-pinitol alleviates diabetic cardiomyopathy by inhibiting the optineurin-mediated endoplasmic reticulum stress and glycophagy signaling pathway

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-1158P	advanced glycation end-products (AGEs)	38311336	7.2	Other	Other

93. PANoptosis subtypes predict prognosis and immune efficacy in gastric cancer

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
BA00208	CCK8 kit	38347337	7.2	Other	Other

94. Cichoric acid targets RANKL to inhibit osteoclastogenesis and prevent ovariectomy-induced bone loss

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bsm-33117M	anti-p65 antibody	38358727	7.2	IF	Mouse

95. Lactate/GPR81 recruits regulatory T cells by modulating CX3CL1 to promote immune resistance in a highly glycolytic gastric cancer

Skus	ProductNames	Pubmed id	Journal_tier	Applications	CrossReactivity
bsm-52079R	anti-FOXP3 antibody	38419759	7.2	WB,IHC	Human,Mouse
bs-0811R	anti-CX3CL1 antibody	38419759	7.2	WB,IHC,IF	Human,Mouse

96.3D printed fibroblast-loaded hydrogel for scleral remodeling to prevent progression of myopia

Skus	ProductNames	Pubmed id	Journal_tier	Applications	CrossReactivity
bsm-33400M	anti-collagen type I antibody	38362614	7	IF	Mouse

97.Vitamin A regulates mitochondrial biogenesis and function through p38 MAPK-PGC-1 α signaling pathway and alters the muscle fiber composition of sheep

Skus	ProductNames	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1832R	anti-PGC-1 α antibody	38310300	7	WB	Sheep
bsm-33034M	anti- β -tubulin antibody	38310300	7	WB	Sheep

98.Formyl Peptide Receptor 2-Dependent cPLA2 and 5-LOX Activation Requires a Functional NADPH Oxidase

Skus	ProductNames	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0296G-HRP	Goat-anti-mouse antibody	38397818	7	WB	Human
bs-0295G-HRP	goat-anti-rabbit antibody	38397818	7	WB	Human

99.Amphiphilic multi-targeting copolymer micelles efficiently deliver pZNF580 to promote endothelial cell proliferation and migration

Skus	ProductNames	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0061R	anti- β -actin antibody	38412450	7	WB	Human

100.Establishment of a chicken intestinal organoid culture system to assess deoxynivalenol-induced damage of the intestinal barrier function

Skus	ProductNames	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0539R	anti-chromogranin A antibody	38369477	7	IF	Chicken

101.Inhibition of the ILK-AKT pathway by upregulation of PARVB contributes to the cochlear cell death in Fascin2 gene knockout mice

Skus	ProductNames	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-5444R	anti-p-ILK antibody	38374196	7	WB	Mouse

102.A Mixture of Formic Acid, Benzoic Acid, and Essential Oils Enhanced Growth Performance via Modulating Nutrient Uptake, Mitochondrion Metabolism, and Immunomodulation in Weaned Piglets

Skus	ProductNames	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0061R	anti- β -actin antibody	38397844	7	WB	Pig

103. Intrauterine Growth Restriction Affects Colonic Barrier Function via Regulating the Nrf2/Keap1 and TLR4-NF- κ B/ERK Pathways and Altering Colonic Microbiome and Metabolome Homeostasis in Growing-Finishing Pigs

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1074R	anti-Nrf2 antibody	10.3390/antiox13030283	7	WB	Pig
bs-2013R	anti-p-Nrf2 antibody	10.3390/antiox13030283	7	WB	Pig
bs-3648R	anti-Keap1 antibody	10.3390/antiox13030283	7	WB	Pig
bs-0061R	anti- β -actin antibody	10.3390/antiox13030283	7	WB	Pig

104. Oxidative stress-mediated activation of FTO exacerbates impairment of the epithelial barrier by up-regulating IKBKB via N6-methyladenosine-dependent mRNA stability in asthmatic mice exposed to PM_{2.5}

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0465R	Anti-NF- κ B antibody	38325270	6.8	IF, WB	Mouse, Human

105. p53-dependent HIF-1 α /autophagy mediated glycolysis to support Cr(VI)-induced cell growth and cell migration

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0737R	anti-HIF-1 α antibody	38335577	6.8	WB	Mouse, Human

106. The endogenous opioid system in the medial prefrontal cortex mediates ketamine's antidepressant-like actions

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-3724R	anti-phospho-OPRM1 (Ser375) antibody	38346984	6.8	WB	Rat

107. GPR55 is expressed in glutamate neurons and functionally modulates drug taking and seeking in rats and mice

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-7686R	anti-GPR55 antibody	38374108	6.8	Other	Other

108. Mechanism of KAT2A regulation of H3K36ac in manganese-induced oxidative damage to mitochondria in the nervous system and intervention by curcumin

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-6313R	anti-4-HNE antibody	38417317	6.8	WB	Rat, Human

109. Carrier-free self-assembled nanomedicine based on celastrol and galactose for targeting therapy of hepatocellular carcinoma via inducing ferroptosis

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bsm-60428R	anti-ASGPR antibody	38354520	6.7	WB, IF	Human, Mouse

110.Endometrial mesenchymal stromal/stem cells improve regeneration of injured endometrium in mice

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bsm-52059R	rabbit anti-cytokeratin 19 monoclonal antibody	38347646	6.7	IHC	Mouse

111.SS-31 inhibits mtDNA-cGAS-STING signaling to improve POCD by activating mitophagy in aged mice

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1278R	anti-8-OHdG antibody	38411634	6.7	IF	Mouse

112.Optimization of hair follicle spheroids for the hair-on-a-chip

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-3662R	anti-Wnt10b antibody	38372380	6.6	IF	Human

113.Farnesoid X Receptor Protects Murine Lung Against IL-6-promoted Ferroptosis Induced by Poly(I:C)

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-12867R	anti-FXR antibody	38300138	6.4	WB	Mouse,Human

114.Indoleamine 2,3-dioxygenase 1-mediated iron metabolism in macrophages contributes to lipid deposition in nonalcoholic steatohepatitis

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-15493R	anti-Indoleamine 2,3-dioxygenase 1 (IDO1) antibody	38402297	6.3	IHC	Mouse

115.MT-TN mutations lead to progressive mitochondrial encephalopathy and promotes mitophagy

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-4843R	anti-LC3B antibody	38320662	6.2	WB	Human

116.Spermidine ameliorates osteoarthritis via altering macrophage polarization

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bsm-34028M	anti-F4/80 antibody	38367900	6.2	IF	Mouse

117.Three-dimensional DNA biomimetic networks (B-3D Net)-based ratiometric fluorescence platform for cancer-related gene biosensing

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
------	-------------	-----------	--------------	--------------	-----------------

bs-0437P-AP	Streptavidin-labelled phosphatase enzymes (ALP-SA)	38499419	6.2	Other	Other
-------------	----------------------------------------------------	----------	-----	-------	-------

118. Dietary supplementation with α -ionone alleviates chronic UVB exposure-induced skin photoaging in mice

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-2090R	anti-p53 antibody	38328833	6.1	IHC	Mouse

119. Dysbiosis of gut microbiota and metabolites is associated with radiation-induced colorectal fibrosis and is restored by adipose-derived mesenchymal stem cell therapy

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
IHC001	IHC kit (for Rabbit Primary Antibody)	38350495	6.1	Other	Other

120. Nupr1-mediated vascular smooth muscle cell phenotype transformation involved in methamphetamine induces pulmonary hypertension

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-7106R	anti-Nupr1 antibody	38347241	6.1	IHC,IF	Mouse

121. Maternal Phlorizin Intake Protects Offspring from Maternal Obesity-Induced Metabolic Disorders in Mice via Targeting Gut Microbiota to Activate the SCFA-GPR43 Pathway

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bsm-60016R	rabbit anti-MUC2 antibody	38349207	6.1	IF	Mouse
bs-13536R	rabbit anti-GPR43 antibody	38349207	6.1	WB,IHC	Mouse
bs-0061R	anti- β -actin antibody	38349207	6.1	WB	Mouse

122. Modulatory effect of the fermented Polygonatum cyrtonema Hua on the immune homeostasis and gut integrity in a dextran-sulfate-sodium (DSS)-induced colitis model

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1540R	anti-IL-7 receptor (CD127) antibody	38440931	6.1	IF	Mouse
bs-23110R	anti-RAR-related orphan receptor gamma (ROR γ t) antibody	38440931	6.1	IF	Mouse

123. Melatonin restores DNFB-induced dysbiosis of skin microbiota in a mouse model of atopic dermatitis

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0581R	anti-IL-4 antibody	38387700	6.1	IHC	Mouse
bs-0560R	anti-IL-13 antibody	38387700	6.1	IHC	Mouse

124.β-Nicotinamide Mononucleotide Supplementation Prolongs The Lifespan Of Premature Ageing Mice And Protects Colon Function In Ageing Mice

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-20747R	anti-Lgr5 antibody	38445897	6.1	IF,WB	Mouse

125.Liquiritin reduces chondrocyte apoptosis through P53/PUMA signaling pathway to alleviate osteoarthritis

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
BA00208	CCK8 solution	38423170	6.1	Other	Other

126.DPPA2/4 Promote the Pluripotency and Proliferation of Bovine Extended Pluripotent Stem Cells by Upregulating the PI3K/AKT/GSK3β/β-Catenin Signaling Pathway

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0301R	anti-DNMT3B antibody	38474345	6	WB	Bovine

127.Therapeutic Effects of a Novel Aptamer on Coronaviral Infection-Induced Lung Injury and Systemic Inflammatory Responses

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0828R	Rabbit anti-PAR1 antibody	38474386	6	WB	Mouse

128.Advanced Glycation End Products Upregulate CD40 in Human Retinal Endothelial and Müller Cells: Relevance to Diabetic Retinopathy

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0666R	anti-fibronectin antibody	38474393	6	IF	Human
bs-2929R	anti-human CD40 antibody	38474393	6	IF	Human

129.Zinc Oxide–Selenium Nanoparticles for Inhibiting the Proliferation of Porcine Reproductive and Respiratory Syndrome Virus

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-23941R	anti-PRRSV N antibody	10.1021/acsanm.3c05307	5.9	IF	Monkey

130.NIPBL-mediated RAD21 facilitates tumorigenicity by the PI3K pathway in non-small-cell lung cancer

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-2067R	anti-PI3K antibody	38378967	5.9	WB	Human
bs-5570R	anti-phospho-PI3K (Tyr317) antibody	38378967	5.9	WB	Human

131. Transcript-dependent effects of the CALCA gene on the progression of post-traumatic osteoarthritis in mice

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1860R	anti-CRLR antibody	38396204	5.9	IF	Mouse

132. Isopropyl 3-(3,4-dihydroxyphenyl) 2-hydroxypropanoate protects septic myocardial injury via regulating GAS6/Axl-AMPK signaling pathway

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-5180R	anti-Axl antibody	38301968	5.8	WB	Mouse
bs-2257R	anti-SIRT1 antibody	38301968	5.8	WB	Mouse

133. Elucidating the mechanisms underlying astrocyte-microglia crosstalk in hippocampal neuroinflammation induced by acute diquat exposure

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-7690R	anti-P2X4R antibody	38305974	5.8	WB	Mouse

134. Selective inhibition of partial EMT-induced tumour cell growth by cerium valence states of extracellular ceria nanoparticles for anticancer treatment

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1172R	anti-N-cadherin antibody	38382224	5.8	ICC	Human

135. Dexamethasone upregulates macrophage PIEZO1 via SGK1, suppressing inflammation and increasing ROS and apoptosis

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0771R	anti-CaMK2 alpha antibody	38354960	5.8	WB	Mouse
bs-3081R	anti-phospho-CaMK2 alpha (Tyr230) antibody	38354960	5.8	WB	Mouse
bsm-33033M	anti-GAPDH antibody	38354960	5.8	WB	Mouse

136. Apical papilla stem cell-derived exosomes regulate lipid metabolism and alleviate inflammation in the MCD-induced mouse NASH model

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0061R	Anti- β -Actin antibody	38395263	5.8	WB	Mouse

137. Cav3.1 T-type calcium channel blocker NNC 55-0396 reduces atherosclerosis by increasing cholesterol efflux

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bsm-52283R	Anti-SR-B1 antibody	38423188	5.8	WB	Human

138. A novel cartilage-targeting MOF-HMME-RGD sonosensitizer combined with sonodynamic therapy to enhance chondrogenesis and cartilage regeneration

Skus	ProductNames	Pubmed id	Journal_tier	Applications	CrossReactivity
bsm-33409M	anti-Col2a1 antibody	38361795	5.7	IHC	Rat

139. Microbial interventions in yak colibacillosis: Lactobacillus-mediated regulation of intestinal barrier

Skus	ProductNames	Pubmed id	Journal_tier	Applications	CrossReactivity
bsm-33102M	Anti-Cytokeratin 18 antibody	38390621	5.7	ICC	Yak

140. USF1 transcriptionally activates USP14 to drive atherosclerosis by promoting EndMT through NLRC5/Smad2/3 axis

Skus	ProductNames	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-22363R	anti-VE-Cadherin antibody	38424494	5.7	WB,IF	Mouse,Human

141. Role of Spexin in White Adipose Tissue Thermogenesis under Basal and Cold-Stimulated Conditions

Skus	ProductNames	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0061R	anti-βactin antibody	38339044	5.6	WB	Mouse
bs-11740R-AF647	anti-Ebf2-Alexa Fluor 647 antibody	38339044	5.6	FC	Mouse

142. ATP5me alleviates high glucose-induced myocardial cell injury

Skus	ProductNames	Pubmed id	Journal_tier	Applications	CrossReactivity
bsm-33199M	anti-cleaved-caspase3 antibody	38320353	5.6	WB	Rat

143. Del-1 Plays a Protective Role against COPD Development by Inhibiting Inflammation and Apoptosis

Skus	ProductNames	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1278R	Anti-8-hydroxy-2'-deoxyguanosine (8-OHdG) antibody	38396634	5.6	IHC	Mouse

144. Fraxin (7-hydroxy-6-methoxycoumarin 8-glucoside) confers protection against ionizing radiation-induced intestinal epithelial injury in vitro and in vivo

Skus	ProductNames	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-2623R	anti-IL-22 antibody	38335653	5.6	WB	Human
bs-2614R	anti-IL-17F antibody	38335653	5.6	WB	Human

145. Baicalin Protects Broilers against Avian Coronavirus Infection via Regulating Respiratory Tract Microbiota and Amino Acid Metabolism

Skus	ProductNames	Pubmed id	Journal_tier	Applications	CrossReactivity
------	--------------	-----------	--------------	--------------	-----------------

bs-2513R	anti-p-IkB- α antibody	38396786	5.6	WB	Chicken
----------	-------------------------------	----------	-----	----	---------

146.Chitosan Versus Dapagliflozin in a Diabetic Cardiomyopathy Mouse Model

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-6313R	Rabbit anti-4-hydroxynonenal (4-HNE) polyclonal antibody	38396795	5.6	IHC	Mouse

147.IL-35 ameliorates lipopolysaccharide-induced endothelial dysfunction by inhibiting endothelial-to-mesenchymal transition

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-2513R	anti-p-IkB α antibody	38335651	5.6	WB	Mouse,Human

148.Atractylenolide-III attenuates osteoarthritis by repolarizing macrophages through inactivating TLR4/NF- κ B signaling

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-18129R	anti-p-IkB α antibody	38346377	5.6	WB	Mouse
bs-0812R	anti-IL-1 β antibody	38346377	5.6	IHC	Rat

149.Pathogenic Th17 cell-mediated liver fibrosis contributes to resistance to PD-L1 antibody immunotherapy in hepatocellular carcinoma

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bsm-43072M	anti-PD-1 antibody	38350354	5.6	IHC	Mouse
bsm-33187M	anti-alpha smooth muscle actin (α -SMA) antibody	38350354	5.6	IHC	Mouse
bs-1103R	anti-PD-L1 antibody	38350354	5.6	IHC,WB	Mouse
bs-1179R	anti-CTLA-4 antibody	38350354	5.6	IHC	Mouse
bs-7363R	anti-TIM-3 antibody	38350354	5.6	IHC	Mouse
bs-0061R	rabbit anti- β -actin antibody	38350354	5.6	WB	Mouse

150.Endothelium-targeted Ddx24 conditional knockout exacerbates ConA-induced hepatitis in mice due to vascular hyper-permeability

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
C2076	O.C.T (optimal cutting temperature) mounting medium	38354508	5.6	Other	Other

151.A visual, rapid, and sensitive detection platform for Vibrio parahaemolyticus based on RPA-CRISPR/Cas12a and an immunochromatographic test strip

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
-----	-------------	-----------	--------------	--------------	-----------------

bs-0295G	goat anti-rabbit IgG antibody	10.1093/fqsafe/fyae008	5.6	Other	Other
----------	-------------------------------	------------------------	-----	-------	-------

152. Mechanical strain regulates osteogenesis via Antxr1/LncRNA H19/Wnt/ β -catenin axis

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
AK352	Alkaline phosphatase (ALP) assay kit	38358001	5.6	Other	Other

153. Protective effects of functional Nano-Selenium supplementation on spleen injury through regulation of p38 MAPK and NF- κ B protein expression

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-0637R	anti-p38 MAPK antibody	38367461	5.6	WB,IHC	Rat
bs-20355R	anti-NF- κ B antibody	38367461	5.6	WB,IHC	Rat
bs-0636R	anti-p-p38 MAPK antibody	38367461	5.6	WB,IHC	Rat
bs-6309R	anti-IL-6 antibody	38367461	5.6	WB,IHC	Rat
bs-41373R	anti-GAPDH antibody	38367461	5.6	WB,IHC	Rat

154. The Complement Component 4 Binding Protein α Gene: A Versatile Immune Gene That Influences Lipid Metabolism in Bovine Mammary Epithelial Cell Lines

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-5069R	anti-HMGCS1 antibody	38397050	5.6	WB	Bovine
bs-23779R	anti-CPT1A antibody	38397050	5.6	WB	Bovine
bs-3787R	anti-SCD1 antibody	38397050	5.6	WB	Bovine
bs-5060R	anti-FADS1 antibody	38397050	5.6	WB	Bovine

155. Exosomal miR-223-3p from bone marrow mesenchymal stem cells targets HDAC2 to downregulate STAT3 phosphorylation to alleviate HBx-induced ferroptosis in podocytes

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-0297P	Human IgG	38444939	5.6	Other	Other

156. Synergistic Anti-Cancer Effects of ERB-041 and Genistein through Estrogen Receptor Suppression-Mediated PI3K/AKT Pathway Downregulation in Canine Mammary Gland Tumor Cells

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-0686R	rabbit anti-PTEN polyclonal antibody	38473712	5.6	WB,ICC	Dog

157. Knockdown of DNMT1 Induces SLCO3A1 to Promote Follicular Growth by Enhancing the Proliferation of Granulosa Cells in Mammals

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-20387R	anti-STAR antibody	38473715	5.6	WB	Human
bs-6656R	anti-CCNB2 antibody	38473715	5.6	WB	Human

158.Indoleamine 2, 3-dioxygenase 1 activation in macrophage exacerbates hepatic ischemia-reperfusion injury by triggering hepatocyte ferroptosis

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-15493R	anti-Indoleamine 2, 3-dioxygenase 1 (IDO-1) polyclonal antibody	38382261	5.6	IHC	Mouse

159.Modulating Th1/Th2 drift in asthma-related immune inflammation by enhancing bone mesenchymal stem cell homing through targeted inhibition of the Notch1/Jagged1 signaling pathway

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1448R	anti-Jagged1 antibody	38387192	5.6	WB	Rat

160.Pterostilbene Induces Apoptosis from Endoplasmic Reticulum Stress Synergistically with Anticancer Drugs That Deposit Iron in Mitochondria

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-9043R	anti-PDZD8 antibody	38473857	5.6	IP,WB	Human

161.Salvia officinalis L. exerts oncostatic effects in rodent and in vitro models of breast carcinoma

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bsm-33277M	anti-caspase-3 antibody	38464730	5.6	IHC	Rat

162.Combined metabolomic and proteomic analysis of sepsis related acute liver injury and its pathogenesis research

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-9879R	anti-Epb42 antibody	38412671	5.6	WB	Mouse

163.The Upregulation of Leucine-Rich Repeat Containing 1 Expression Activates Hepatic Stellate Cells and Promotes Liver Fibrosis by Stabilizing Phosphorylated Smad2/3

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-8853R	anti-p-Smad2/3 (Thr8) antibody	38473980	5.6	IP,WB	Mouse,Human
bs-23445R	anti-TGFβR I antibody	38473980	5.6	WB	Human

164.Dexamethasone Inhibits White Adipose Tissue Browning

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0061R	anti-β-ACTIN antibody	38473960	5.6	WB	Rat

165.Effect of high NEFA concentration on lipid metabolism disorders in hepatocytes based on lipidomics

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
-----	-------------	-----------	--------------	--------------	-----------------

BA00208	CCK8 analysis kit	38482059	5.6	Other	Other
---------	-------------------	----------	-----	-------	-------

166. Effect of the *Pseudopleuronectes americanus*-derived Pleurocidin on DSS-induced Ulcerative colitis in mice and its preliminary molecular mechanisms

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1329R	Rabbit anti-ZO-1/TJP1 antibody	38422770	5.6	IF,WB	Mouse
bs-10011R	Rabbit anti-Occludin antibody	38422770	5.6	WB	Mouse
bs-1428R	Rabbit anti-CLDN1 antibody	38422770	5.6	WB	Mouse
bs-41373R	Rabbit anti-GAPDH antibody	38422770	5.6	WB	Mouse
bs-5662R	Rabbit anti-phospho-NFKB p65 antibody	38422770	5.6	WB	Mouse
bs-20160R	Rabbit anti-NFKB p65 antibody	38422770	5.6	WB	Mouse
bs-2637R	Rabbit anti-ERK1 + ERK2 antibody	38422770	5.6	WB	Mouse
bs-1522R	Rabbit anti-phospho-ERK1 + 2 antibody	38422770	5.6	WB	Mouse
bs-2592R	Rabbit anti-JNK1 + JNK2 + JNK3 antibody	38422770	5.6	WB	Mouse
bs-4163R	Rabbit anti-phospho-JNK1 + JNK2 + JNK3 antibody	38422770	5.6	WB	Mouse
bs-0636R	Rabbit anti-Phospho-p38 MAPK antibody	38422770	5.6	WB	Mouse
bs-0637R	Rabbit anti-p38 MAPK antibody	38422770	5.6	WB	Mouse
bs-0295G-HRP	Goat Anti-Rabbit IgG H&L / HRP antibody	38422770	5.6	WB	Mouse

167. The therapeutic effect of anti-CD19 antibody on DHEA-induced PCOS mice

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0061R	rabbit anti- β -actin antibody	38428145	5.6	WB	Mouse,Human
bsm-33033M	mouse anti-GAPDH antibody	38428145	5.6	WB	Human

168. Construction of a prognostic prediction model in liver cancer based on genes involved in integrin cell surface interactions pathway by multi-omics screening

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0032R	anti-Bcl-2 antibody	38374893	5.5	WB	Human
bs-2203R	anti-ITGAV antibody	38374893	5.5	WB	Human

bs-0637R	anti-MAPK antibody	38374893	5.5	WB	Human
----------	--------------------	----------	-----	----	-------

169.Mechanism study of cross presentation of exogenous antigen induced by cholera toxin-like chimeric protein

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0295G-HRP	HRP-conjugated goat anti-rabbit IgG antibody	38320931	5.5	ELISA	Mouse

170.TERT mediates the U-shape of glucocorticoids effects in modulation of hippocampal neural stem cells and associated brain function

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0233R	anti-TERT antibody	38421107	5.5	WB	Mouse

171.Inhibiting acid-sensing ion channel exerts neuroprotective effects in experimental epilepsy via suppressing ferroptosis

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-2586R	rabbit anti-ASIC1a polyclonal antibody	38357854	5.5	IF	Human

172.Pathological features and molecular signatures of early olfactory dysfunction in 3xTg-AD model mice

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-2392R	anti-pS404 antibody	38366763	5.5	WB	Mouse

173.Melatonin alleviates oxidative stress damage in mouse testes induced by bisphenol A

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
S0152	hematoxylin	38440074	5.5	Other	Other
C02-02003	ethylene diamine tetraacetic acid antigen retrieval buffer	38440074	5.5	Other	Other

174.Integrated network pharmacology, metabolomics, and transcriptomics of Huanglian-Hongqu herb pair in non-alcoholic fatty liver disease

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-3614R	anti-PPAR α antibody	38325669	5.4	WB	Rat

175.Calycosin enhances treg differentiation for alleviating skin inflammation in atopic dermatitis

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-23110R	anti-ROR γ t antibody	38331120	5.4	IHC	Mouse

176. Mechanisms underlying the synergistic effects of chuanxiong combined with Chishao on treating acute lung injury based on network pharmacology and molecular docking combined with preclinical evaluation

SKU	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1840R	anti-Lamin B1 antibody	38342157	5.4	WB	Mouse
bsm-33039M	anti- α -tubulin antibody	38342157	5.4	WB	Mouse

177. Anji white tea relaxes precontracted arteries, represses voltage-gated Ca²⁺ channels and voltage-gated K⁺ channels in the arterial smooth muscle cells: Comparison with green tea main component (–)-epigallocatechin gallate

SKU	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-16865R	rabbit anti-Kv1.2 antibody	38346524	5.4	IF	Rat

178. Shaoyao decoction improves damp-heat colitis by activating the AHR/IL-22/STAT3 pathway through tryptophan metabolism driven by gut microbiota

SKU	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1329R	anti-ZO-1 antibody	38342152	5.4	WB,IF	Mouse

179. Study on therapeutic mechanism of total salvianolic acids against myocardial ischemia-reperfusion injury based on network pharmacology, molecular docking, and experimental study

SKU	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1183R	Rabbit Anti-IL-17 polyclonal antibody	38360382	5.4	WB	Rat
bs-2606R	Rabbit Anti-IL-17RA polyclonal antibody	38360382	5.4	WB	Rat

180. Inhibition of p38 MAPK/NF- κ B p65 signaling pathway activity by rare ginsenosides ameliorates cyclophosphamide-induced premature ovarian failure and KGN cell injury

SKU	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-2130R	rabbit anti-Ki67 antibody	38382656	5.4	IHC	Rat
bs-0637R	anti-p38MAPK antibody	38382656	5.4	WB	Rat,Human
bs-0636R	anti-phospho-p38MAPK antibody	38382656	5.4	WB	Rat,Human
bs-5662R	anti-phospho-NF- κ Bp65 antibody	38382656	5.4	WB	Rat,Human

181. Investigation of healing strategies in a rat corneal opacity model with polychromatic light and stem cells injection

SKU	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bsm-33400M	anti-collagen type I (COL1) antibody	38422971	5.4	IHC	Rat
bs-0948R	anti-collagen type III (COL3) antibody	38422971	5.4	IHC	Rat

bsm-33345M	anti-transforming growth factor beta 1 (TGF- β 1) antibody	38422971	5.4	IHC	Rat
bs-6319R	anti-interleukin 1 beta (IL-1 β) antibody	38422971	5.4	IHC	Rat
bs-11054R	anti-keratocan (KERA) antibody	38422971	5.4	IHC	Rat

182.The rhizomes of *Atractylodes macrocephala* Koidz improve gastrointestinal health and pregnancy outcomes in pregnant mice via modulating intestinal barrier and water-fluid metabolism

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1253R	anti-AQP3 antibody	38403003	5.4	IHC	Mouse
bs-1329R	anti-ZO-1 antibody	38403003	5.4	IF	Mouse

183.Neuroprotective effects of Shenghui decoction via inhibition of the JNK/p38 MAPK signaling pathway in an AIC13-induced zebrafish (*Danio rerio*) model of Alzheimer's disease

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-20594R	rabbit anti-TLR4 antibody	38423408	5.4	WB	Zebrafish

184.Investigating the molecular mechanism of Qizhu anticancer prescription in inhibiting hepatocellular carcinoma based on high-resolution mass spectrometry and network pharmacology

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0233R	anti-TERT antibody	38417600	5.4	WB	Mouse
bs-0255R	anti-ESR2 antibody	38417600	5.4	WB	Mouse
bs-0217R	anti-FGF2 antibody	38417600	5.4	IHC,WB	Mouse
bs-4029R	anti-CA9 antibody	38417600	5.4	IHC,WB	Mouse
bs-10274R	anti-IL2 antibody	38417600	5.4	WB	Mouse
bs-0563R	anti-ABCB1 antibody	38417600	5.4	WB	Mouse
bsm-51481M	anti-ESR1 antibody	38417600	5.4	WB	Mouse

185.Gegen Qinlian decoction ameliorates TNBS-induced ulcerative colitis by regulating Th2/Th1 and Tregs/Th17 cells balance, inhibiting NLRP3 inflammasome activation and reshaping gut microbiota

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1329R	anti-ZO-1 antibody	38428658	5.4	WB	Mouse

186.MTMR7 suppresses the phenotypic switching of vascular smooth muscle cell and vascular intimal hyperplasia after injury via regulating p62/mTORC1-mediated glucose metabolism

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-17888R	anti-MTMR7 antibody	38342025	5.3	IF	Mouse
bs-2130R	anti-Ki-67 antibody	38342025	5.3	IHC	Mouse

187.Nerve growth factor promote osteogenic differentiation of dental pulp stem cells through MEK/ERK signalling pathways

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bsm-52176R	rabbit anti-phospho-MEK (p-MEK) antibody	38333908	5.3	WB	Human

188.PI3K δ inhibition alleviates the brain injury during cerebral ischemia reperfusion via suppressing pericyte contraction in a TNF- α dependent manner

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-10297R	anti-TRPV2 antibody	38365134	5.3	IF,WB	Mouse

189.Ndfip1 protected dopaminergic neurons via regulating mitochondrial function and ferroptosis in Parkinson's disease

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0061R	rabbit-anti- β -actin antibody	38365133	5.3	WB	Mouse,Human

190.Neuroimmune pathways involvement in neurodegeneration of R6/2 mouse model of Huntington's disease

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-21460R	rabbit anti-CD47 polyclonal antibody	38444595	5.3	IHC	Mouse

191.AAV-mediated VEGFA overexpression promotes angiogenesis and recovery of locomotor function following spinal cord injury via PI3K/Akt signaling

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-6417R	anti-p-PI3K antibody	38401852	5.3	WB	Mouse

192.MicroRNA-124 negatively regulates STAT3 to alleviate hypoxic-ischemic brain damage by inhibiting oxidative stress

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1278R	rabbit anti-8-OhdG antibody	38319722	5.2	IHC	Rat

193.Targeting PERK-ATF4-P21 axis enhances the sensitivity of osteosarcoma HOS cells to Mpp α -PDT

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-3330R	rabbit anti-phospho-PERK polyclonal antibody	38319715	5.2	IF,WB	Human
bs-0295G-FITC	FITC-conjugated goat anti-rabbit secondary antibody	38319715	5.2	IF	Human
bs-0295G-AF594	goat anti-rabbit Fluor 594 secondary antibody	38319715	5.2	IF	Mouse

C51003	phosphatase inhibitor cocktail	38319715	5.2	Other	Other
PIC008	protease inhibitor cocktail	38319715	5.2	Other	Other
bs-0061R	anti-β-actin antibody	38319715	5.2	WB	Mouse
bs-0755R	anti-GAPDH antibody	38319715	5.2	WB	Human
C05-07004	ECL Plus WB Substrate	38319715	5.2	Other	Other

194. Microfluidic synthesis of dual-encoded microbeads for multiplex immunoassay

SKU	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0297P	Human IgG	10.1016/j.colsurf a.2024.133416	5.2	Other	Other
bs-0297R	rabbit anti-human IgG antibody	10.1016/j.colsurf a.2024.133416	5.2	Other	Other
bs-0296R	rabbit anti-mouse IgG antibody	10.1016/j.colsurf a.2024.133416	5.2	Other	Other
bs-0294R	rabbit anti-goat IgG antibody	10.1016/j.colsurf a.2024.133416	5.2	Other	Other
bs-0297R-AF488	Alexa Fluor 488-conjugated rabbit anti-human IgG antibody	10.1016/j.colsurf a.2024.133416	5.2	Other	Other

195. Screening of Active Substances Regulating Alzheimer's Disease in Ginger and Visualization of the Effectiveness on 6-Gingerol Pathway Targets

SKU	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
A5300P	Aβ1-42	38397589	5.2	Other	Other

196. Isolation, identification, recombination analysis and pathogenicity experiment of a PRRSV recombinant strain in Sichuan Province, China

SKU	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-23941R-FITC	fluorescein isothiocyanate-labeled anti-PRRSV N polyclonal antibody	38450173	5.2	IF	Pig

197. Hierarchical zeolite coatings featuring a spatial gradient architecture for sequentially-controlled bisphosphonate release in the modulation of osteogenic-osteoclastic balance

SKU	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
C01-01001	Phosphate-buffered saline (PBS)	10.1016/j.micro meso.2024.1130 60	5.2	Other	Other
C2055	4% paraformaldehyde	10.1016/j.micro meso.2024.1130 60	5.2	Other	Other

C8013	phalloidin-FITC	10.1016/j.micro meso.2024.1130 60	5.2	Other	Other
-------	-----------------	-----------------------------------------	-----	-------	-------

198. Multi-omic analysis of precocious puberty girls: pathway changes and metabolite validation

Skus	ProductNames	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-10369R	anti-GnRH antibody	38487340	5.2	IHC	Rat

199. Casticin ameliorates osteoarthritic cartilage damage in rats through PI3K/AKT/HIF-1 α signaling

Skus	ProductNames	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-5570R	anti-p-PI3K antibody	38309612	5.1	IHC,WB	Rat

200. Low-Shear Stress Promotes Atherosclerosis via Inducing Endothelial Cell Pyroptosis Mediated by IKK ϵ /STAT1/NLRP3 Pathway

Skus	ProductNames	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-8583R	anti-p-IKK ϵ Ser172 antibody	38315275	5.1	IF	Mouse

201. The Synergistic Effect Study of Lipopolysaccharide (LPS) and A53T- α -Synuclein: Intranasal LPS Exposure on the A53T- α -Synuclein Transgenic Mouse Model of Parkinson's Disease

Skus	ProductNames	Pubmed id	Journal_tier	Applications	CrossReactivity
D-9106	DAPI	38367134	5.1	Other	Other

202. Ceftriaxone Modulates Ubiquitination of α -Amino-3-Hydroxy-5-Methyl-4-Isoxazole Propionic Acid Receptors to Improve Long-Term Potentiation Impairment Induced by Exogenous β -Amyloid in a Glutamate Transporter-1 Dependent Manner

Skus	ProductNames	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0107P	A β 1-42 (A β)	38374316	5.1	Other	Other

203. Farrerol Alleviates Cerebral Ischemia-Reperfusion Injury by Promoting Neuronal Survival and Reducing Neuroinflammation

Skus	ProductNames	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0295G-AF555	goat anti-rabbit IgG-AlexaFluor-555 conjugate antibody	38376762	5.1	IF	Mouse

204. High Dose of Estrogen Protects the Lungs from Ischemia-Reperfusion Injury by Downregulating the Angiotensin II Signaling Pathway

Skus	ProductNames	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-19605R	anti-Cavin-2 antibody	38386131	5.1	WB,IF	Mouse

205. The preventive effects of Lactobacillus casei 03 on Escherichia coli-induced mastitis in vitro and in vivo

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-20159R	anti-NF-κB p65 antibody	38395896	5.1	WB	Bovine
bs-0982R	anti-NF-κB phospho-p65 antibody	38395896	5.1	WB	Bovine
bs-18129R	anti-phospho-IκBα antibody	38395896	5.1	WB	Bovine
bsm-33036M	anti-β-actin antibody	38395896	5.1	WB	Bovine

206. One-pot synthesis of triple responsive composite hydrogel based on hyaluronic acid and its targeted controlled drug release

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
D13032	Hyaluronidase (HAase, 400 Unit/mg)	10.1016/j.reactfunctpolym.2024.105867	5.1	Other	Other

207. Genome-wide search links senescence-associated secretory proteins with susceptibility for coronary artery disease in mouse and human

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-4595R	anti-COL4A1 antibody	38416803	5.1	IF, WB	Mouse
bs-2095R	anti-ITGA1 antibody	38416803	5.1	IF, WB	Mouse
bsm-60676R	anti-LRP1 antibody	38416803	5.1	IF, WB	Mouse
bs-0294D-Cy5	Cy5-conjugated donkey anti-goat IgG antibody	38416803	5.1	IF	Mouse
bs-0295D-Cy3	Cy3-conjugated donkey anti-rabbit IgG antibody	38416803	5.1	IF	Mouse

208. Inhibition of AMPK activation in Echinococcus granulosus sensu stricto limits the parasite's glucose metabolism and survival

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-80295G-HRP	Goat Anti-Rabbit IgG H&L / HRP antibody	38349157	4.9	WB	ococcus granu
bs-0296G-HRP	Goat Anti-Mouse IgG H&L / HRP antibody	38349157	4.9	WB	ococcus granu

209. Total flavonoids of Astragalus protects glomerular filtration barrier in diabetic kidney disease

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-10657R	anti-PI3 Kinase p110 beta antibody	38365794	4.9	WB	Mouse

210. Fabrication of vascularized tissue-engineered bone models using triaxial bioprinting

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0470R	anti-Osteocalcin Polyclonal antibody	38411369	4.9	IF	Human, Mouse

211. Machine learning combined with single-cell analysis reveals predictive capacity and immunotherapy response of T cell exhaustion-associated lncRNAs in uterine corpus endometrial carcinoma

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0061R	anti-β-Actin antibody	38311301	4.8	WB	Human

212. MiR-21-5p ameliorates Busulfan-induced testicular dysfunction and maintains spermatogenesis1

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-5971R	anti-ZBTB16 antibody	10.1016/j.jia.2024.02.004	4.8	IF	Goat

213. Rap1GAP exacerbates myocardial infarction by regulating the AMPK/SIRT1/NF-κB signaling pathway

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0489P	bromo-deoxyuridine	38320624	4.8	Other	Other

214. Zishen Yutai Pills Promote Angiogenesis at the Maternal-Fetal Interface in Recurrent Spontaneous Abortion Mice by Regulating miR-187/VEGF Axis

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
BA00208	CCK8 kit	38370565	4.8	Other	Other

215. LncRNA CHROMR/miR-27b-3p/MET axis promotes the proliferation, invasion and contributes to rituximab resistance in diffuse large B-cell lymphoma

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
BA00208	CCK-8	38367665	4.8	Other	Other

216. Vitamin D status alters genes involved in ovarian steroidogenesis in muskrat granulosa cells

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-5057R	Rabbit anti-DHCR7 antibody	38402945	4.8	IHC	Muskrat
bs-3900R	Rabbit anti-CYP2R1 antibody	38402945	4.8	IHC	Muskrat
bs-14146R	Rabbit anti-CYP27B1 antibody	38402945	4.8	IHC	Muskrat
bs-3899R	Rabbit anti-CYP24A1 antibody	38402945	4.8	IHC	Muskrat

bs-10618R	Rabbit anti-VDR antibody	38402945	4.8	IHC	Muskrat
bs-3570R	Rabbit anti-StAR antibody	38402945	4.8	IHC	Muskrat
bs-23683R	Rabbit anti-VDAC antibody	38402945	4.8	IHC	Muskrat
bs-10099R	Rabbit anti-P450scc antibody	38402945	4.8	IHC	Muskrat
bs-3906R	Rabbit anti-3 β HSD antibody	38402945	4.8	IHC	Muskrat
bs-3853R	Rabbit anti-P450c17 antibody	38402945	4.8	IHC	Muskrat

217.Coxsackievirus and adenovirus receptor inhibits tilapia lake virus infection via binding to viral segment 8 and 10 encoded protein

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0296G-HRP	HRP-conjugated Goat anti-Mouse antibody	38341116	4.7	WB	ochromis nilot
bs-0295G-HRP	HRP-conjugated Goat anti-Rabbit antibody	38341116	4.7	WB	ochromis nilot

218.Platinum-based neoadjuvant chemotherapy upregulates STING/IFN pathway expression and promotes TILs infiltration in NSCLC

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1103R	anti-PD-L1 antibody	38425343	4.7	IHC	Human
bs-23731R	anti-IFN- β antibody	38425343	4.7	IHC	Human

219.Long non-coding RNA KB-1460A1.5 promotes ferroptosis by inhibiting mTOR/SREBP-1/SCD1-mediated polyunsaturated fatty acid desaturation in glioma

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-6313R	anti-4-HNE antibody	38422369	4.7	IHC	Mouse

220.Identification and preliminary analysis of hub genes associated with bladder cancer progression by comprehensive bioinformatics analysis

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
IHC001	immunohistochemistry detection system kit (Rabbit)	38307969	4.6	Other	Other

221.Comparison of preparation methods of rat kidney single-cell suspensions

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0527R-PE	PE-conjugated anti-CD10 antibody	38307992	4.6	FC	Rat

222. Anti-psoriasis molecular targets and active components discovery of Optimized Yinxieling Formula via affinity-purified strategy

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1193R	anti-IL-23 antibody	38342565	4.6	IHC	Mouse

223. Acquired curved hair is caused by fusion of multiple hair matrix cells

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1948R	rabbit anti-Wnt5a polyclonal antibody	38431439	4.6	IHC	Human
bs-1700R	rabbit anti-Wnt3a polyclonal antibody	38431439	4.6	IHC	Human

224. The GRHL3-regulated long non-coding RNA Inc-DC modulates keratinocytes differentiation by interacting with IGF2BP2 and up-regulating ZNF750

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-6327R	anti-FLG antibody	38383230	4.6	WB	Human
bs-23062R	anti-IVL antibody	38383230	4.6	WB,IF	Human
bsm-52051R	anti-KRT1 antibody	38383230	4.6	WB	Human

225. Piceatannol Alleviates Deoxynivalenol-Induced Damage in Intestinal Epithelial Cells via Inhibition of the NF- κ B Pathway

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1329R	anti-ZO-1 antibody	38398607	4.6	WB	Pig

226. Calcitriol promotes M2 polarization of tumor-associated macrophages in 4T1 mouse mammary gland cancer via the induction of proinflammatory cytokines

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-2987R	anti-VDR antibody	38355711	4.6	WB	Human

227. Ripa-56 protects retinal ganglion cells in glutamate-induced retinal excitotoxic model of glaucoma

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
C05-02001	BCA protein assay kit	38360971	4.6	Other	Other

228. Forensic significance of intracardiac expressions of Nrf2 in acute myocardial ischemia

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-2673R	rabbit anti-human C5b-9 polyclonal antibody	38374168	4.6	IHC	Human
bs-1074R	rabbit anti-human Nrf2 polyclonal antibody	38374168	4.6	IHC	Human

229. Chromium Picolinate Regulates Bone Metabolism and Prevents Bone Loss in Diabetic Rats

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0431R	Rabbit anti-OPG antibody	38474436	4.6	IHC	Rat
bs-0747R	Rabbit anti-RANKL antibody	38474436	4.6	IHC	Rat
bs-1134R	Rabbit anti-RUNX 2 antibody	38474436	4.6	IHC	Rat

230. An Inulin-Type Fructan CP-A from Codonopsis pilosula Alleviated 5-Fluorouracil-Induced Intestinal Mucositis via the ERK/MLCK/MLC2 Pathway and Regulation of Gut Microbiota

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1428R	anti-Claudin-1 antibody	10.3390/ph17030297	4.6	IHC,WB	Rat
bs-1571R	anti-F-actin antibody	10.3390/ph17030297	4.6	IHC,WB	Rat

231. Regulator of Ribosome Synthesis 1 (RRS1) Stabilizes GRP78 and Promotes Breast Cancer Progression

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0297R	Rabbit Anti-Human IgG H&L antibody	38474562	4.6	CoIP,WB	Human

232. Regulation of the Nur77-P2X7r Signaling Pathway by Nodakenin: A Potential Protective Function against Alcoholic Liver Disease

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-10196R	anti- α -SMA antibody	38474588	4.6	Other	Other
bs-0578R	anti-collagen type I (collagen-I) antibody	38474588	4.6	Other	Other

233. SDF-1 involvement in orthodontic tooth movement after tooth extraction

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-4938R	rabbit anti-SDF-1/CXCL12 polyclonal antibody	38424199	4.6	IF	Rat
bs-0812R	anti-IL-1 β antibody	38424199	4.6	IHC	Rat
bs-0782R	anti-IL-6 antibody	38424199	4.6	IHC	Rat
bs-0747R	anti-RANKL antibody	38424199	4.6	IHC	Rat
bs-1611R	anti-Cathepsin K antibody	38424199	4.6	IHC	Rat

234.The Humanization and Maturation of an Anti-PrPc Antibody

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0297P	native IgG	10.3390/bioengineering11030242	4.6	Other	Other

235.Hsa_circ_0008133 contributes to lung cancer progression by promoting glycolysis metabolism through the miR-760/MEX3A axis

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0101R	anti-pyruvate kinase M1/2 (PKM2) antibody	38317294	4.5	WB	Human
bs-4593R	anti-matrix metalloproteinase 9 (MMP9) antibody	38317294	4.5	WB	Human
bs-0061R	anti- β -actin (loading control) antibody	38317294	4.5	WB	Mouse,Human
bs-0295G-HRP	goat anti-rabbit IgG heavy and light/horseradish peroxidase (H&L/HRP) antibody	38317294	4.5	WB	Mouse,Human
C05-07004	BiossECL Plus WB Substrate	38317294	4.5	Other	Other
bsm-33070M	anti-Ki-67 antibody	38317294	4.5	IHC	Mouse

236.Bioinformatics and Integrative Experimental Method to Identifying and Validating Co-Expressed Ferroptosis-Related Genes in OA Articular Cartilage and Synovium

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0519R	anti-activating transcription factor 3 (ATF3) antibody	38370466	4.5	IHC,WB	Mouse
bs-0741R	anti-cyclin-dependent kinase inhibitor 1A (CDKN1A) antibody	38370466	4.5	IHC,WB	Mouse

237.The in vitro and in vivo antiviral effects of IGF1R inhibitors against respiratory syncytial virus infection

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-4985R	rabbit anti-IGF1R antibody	38299600	4.4	WB	Human
bs-41373R	anti-GAPDH antibody	38299600	4.4	WB	Human
bs-1264R	rabbit anti-RSV G protein antibody	38299600	4.4	IF	Mouse,Human

238.Ferroptosis as a Potential Therapeutic Target for Reducing Inflammation and Corneal Scarring in Bacterial Keratitis

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0480R	rabbit anti-IFN- γ antibody	38381413	4.4	IHC	Mouse

239.M2 macrophages promote subconjunctival fibrosis through YAP/TAZ signalling

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-13455R	anti-fibronectin antibody	38335557	4.4	IF, WB	Human
bs-3605R	anti-YAP1 antibody	38335557	4.4	WB	Human
bs-12367R	anti-TAZ antibody	38335557	4.4	WB	Human

240.Acute stress transiently activates macrophages and chemokines in cervical lymph nodes

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-10234R	rabbit anti-CXCL1 polyclonal antibody	38351242	4.4	IHC	Rat

241.Chlorogenic Acid Ameliorates Intestinal Inflammation by Inhibiting NF- κ B and Endoplasmic Reticulum Stress in Lipopolysaccharide-Challenged Broilers

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-2075R	anti-HO-1 antibody	38442474	4.4	WB	Chicken

242.B-GOS alleviates olanzapine-induced lipid disturbances in mice by enriching Akkermansia and upregulation of PGRMC1-Wnt signaling

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-13536R	anti-GPR43 antibody	38325638	4.3	WB	Mouse

243.Exosomal U2AF2 derived from human bone marrow mesenchymal stem cells attenuates the intervertebral disc degeneration through circ_0036763/miR-583/ACAN axis

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1223R	anti-ACAN antibody	38362337	4.3	WB	Human
bs-7158R	anti-Collagen I antibody	38362337	4.3	WB	Human
bs-10589R	anti-Collagen II antibody	38362337	4.3	WB	Human
bs-1523R	anti-CD63 antibody	38362337	4.3	WB	Human
bs-6934R	anti-CD81 antibody	38362337	4.3	WB	Human

244.The VP1/2 Protein of a New Recombinant PRV Strain Promotes the Infectivity and Pathogenicity of PRV in Northeastern China

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
-----	-------------	-----------	--------------	--------------	-----------------

bs-0295G-AF488	Goat Anti-Rabbit IgG H&L / AF488 antibody	10.1155/2024/1575103	4.3	IF	Pig
bs-0368G-AF594	Goat Anti-Mouse IgM / AF594 antibody	10.1155/2024/1575103	4.3	IF	Pig

245.Effect of CB1 Receptor Deficiency on Mitochondrial Quality Control Pathways in Gastrocnemius Muscle

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0061R	anti-B-ACTIN antibody	38392333	4.2	WB	Mouse

246.Effect of hypoxia-HIF-1 α -periostin axis in thyroid cancer

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0737R	rabbit anti-HIF-1 α polyclonal antibody	38391012	4.2	IHC	Human

247.A one-step process for multi-gradient wettability modification on a polymer surface

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-3514R-FITC	FITC-labelled anti-D-dimer detection monoclonal antibody	38421308	4.2	Other	Other
bs-3514R	anti-D-dimer capture monoclonal antibody	38421308	4.2	Other	Other
bs-41172P	human D-dimer biomarker	38421308	4.2	Other	Other

248.m6A methyltransferase METTL14-mediated RP1-228H13.5 promotes the occurrence of liver cancer by targeting hsa-miR-205/ZIK1

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0295G-HRP	Goat Anti-Rabbit IgG H&L / HRP antibody	38426536	4.2	WB	Human

249.Hyperbaric oxygen preconditioning normalizes scrotal temperature, sperm quality, testicular structure, and erectile function in adult male rats subjected to exertional heat injury

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1329R	Rabbit anti-ZO-1 antibody	38341020	4.1	IHC	Rat

250.Protective effect of bioactive components from Rubi fructus against oxidative damage in human ovarian granulosa cells induced by 2,2-azobis (2-methylpropionamide) dihydrochloride

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-2469R	anti-PERK antibody	38349056	4.1	WB	Human

251. Inhibition of Oxidative Stress-Induced Ferroptosis Can Alleviate Rheumatoid Arthritis in Human

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-3577R	Rabbit anti-DMT1 antibody	10.1155/2024/9943747	4.1	WB	Human
bs-9894R	Rabbit anti-TFR2 antibody	10.1155/2024/9943747	4.1	WB	Human
bs-3884R	Rabbit anti-GPX4 antibody	10.1155/2024/9943747	4.1	WB,IF	Human

252. Yellow leaf green tea modulates AMPK/ACC/SREBP1c signaling pathway and gut microbiota in high-fat diet-induced mice to alleviate obesity

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-2047R	anti-CPT-1 α antibody	38407390	4.1	WB	Mouse
bs-1402R	anti-SREBP-1c antibody	38407390	4.1	WB	Mouse

253. Bushenhuo decoction improves polycystic ovary syndrome by regulating exosomal miR-30a-5p/SOCS3/mTOR/NLRP3 signaling-mediated autophagy and pyroptosis

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-8878R	anti-LC3II/I antibody	38302986	4	WB	Rat
bs-1353R	anti-Beclin1 antibody	38302986	4	WB	Rat
bs-4005R	anti-ATG5 antibody	38302986	4	WB	Rat
bs-3602R	anti-ULK1 antibody	38302986	4	WB	Rat
bs-3495R	anti-p-mTOR antibody	38302986	4	WB	Rat
bs-1992R	anti-mTOR antibody	38302986	4	WB	Rat
bs-24250R	anti-SOCS3 antibody	38302986	4	WB	Other
bs-1365R	anti-TSG101 antibody	38302986	4	WB	Rat
bs-23032R	anti-CD63 antibody	38302986	4	WB	Rat
bsm-52639R	anti-calnexin antibody	38302986	4	WB	Rat
bs-0061R	anti- β -actin antibody	38302986	4	WB	Rat

254. Progesterone and adipoQ receptor 7 (PAQR7) mediate the anti-apoptotic effect of P4 on human granulosa cells and its deficiency reduces ovarian function in female mice

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-17741R	anti-PAQR7 antibody	38317224	4	WB	Mouse, Human

255. Identification of PKM2 as a pyroptosis-related key gene aggravates senile osteoporosis via the NLRP3/Caspase-1/GSDMD signaling pathway

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-1134R	anti-RUNX2 antibody	38342404	4	WB	Human
bs-6741R	anti-ASC antibody	38342404	4	WB	Human

256. Inhibition of PTPN3 Expressed in Activated Lymphocytes Enhances the Antitumor Effects of Anti-PD-1 Therapy in Head and Neck Cancer, Especially in Hypoxic Environments

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-4029R-AF488	Alexa Fluor 488-conjugated CA-IX antibody	38297883	3.9	IF	Mouse

257.Grape Seed Proanthocyanidins Protect Pancreatic β Cells Against Ferroptosis via the Nrf2 Pathway in Type 2 Diabetes

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-6313R	Anti-4 hydroxynonenal (4-HNE) antibody	38367173	3.9	IHC	Rat

258.Identification of a Novel NLRP12 Frameshift Mutation (Val730Glyfs41) by Whole-Exome Sequencing in Patients with Crohn's Disease

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0812R	rabbit anti-IL-1 β antibody	10.1155/2024/5573272	3.9	IF	Human

259.Chronic stress promotes gastric cancer progression via the adrenoceptor beta 2/PlexinA1 pathway

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-21452R	anti-ADRB2 antibody	38331165	3.8	IHC,WB	Human,Mouse
bs-2692R	anti-PlexinA1 antibody	38331165	3.8	IHC,WB	Human,Mouse

260.The effect of c-Fos on the prognosis, proliferation, and invasion of oral squamous cell carcinoma

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-10172R	anti-c-Fos antibody	38326698	3.8	IHC	Human

261.The role of ERp29/FOS/EMT pathway in excessive apoptosis of placental trophoblast cells in intrahepatic cholestasis of pregnancy

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
D-9106	DAPI	38346375	3.8	Other	Other

262.Dual effect of vitamin D3 on breast cancer-associated fibroblasts

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-10618R	rabbit anti-VDR antibody	38360633	3.8	WB	Human

263.Deciphering relationship between depression and microbial molecules based on multi-omics: A case study of Chaigui Granules

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-17244R	anti-SAT1 antibody	10.1016/j.chmed.2023.12.003	3.8	WB	Rat

bs-1294R	anti-ODC antibody	10.1016/j.chmed .2023.12.003	3.8	WB	Rat
----------	-------------------	---------------------------------	-----	----	-----

264. Isolation and characterization of genotype 4 Eurasian avian-like H1N1 influenza virus in pigs suffering from pneumonia

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-0310G-HRP	HRP-conjugated goat anti-chicken antibody	38330852	3.7	ELISA	Chicken

265. Porphyromonas gingivalis GroEL exacerbates orthotopic allograft transplantation vasculopathy via impairment of endothelial cell function

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-0646R-Cy5	Cy5-conjugated anti-CD34 antibody	38311876	3.7	FC	Rat

266. Regulation of Onecut2 by miR-9-5p in Japanese encephalitis virus infected neural stem/progenitor cells

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-19643R	anti-Onecut2 (OC2) antibody	38319106	3.7	WB	Mouse, Human

267. Vaginal pH value can affect the susceptibility to human papillomavirus infection

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-0061R	anti-β-actin antibody	38331746	3.7	WB	Human

268. PI3K/AKT/mTORC1 signalling pathway regulates MMP9 gene activation via transcription factor NF-κB in mammary epithelial cells of dairy cows

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-0982R	anti-p-P65 antibody	38343377	3.7	WB	Cow

269. A genome-wide CRISPR screen identified host genes essential for intracellular Brucella survival

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-2229R	anti-Brucella antibody	38376367	3.7	IF	Human

270. An angiogenesis-associated gene-based signature predicting prognosis and immunotherapy efficacy of head and neck squamous cell carcinoma patients

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-8512R	rabbit anti-MSX1 antibody	38347320	3.6	IHC, WB	Human

271. Lesional Infiltration of Eosinophils, Basophils, and M2 Macrophages Expressing Eotaxin-1 and Eotaxin-3 in Granuloma Faciale, but not in Erythema Elevatum Diutinum

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-2527R	anti-CD163 antibody	38348729	3.6	IF	Human

272. The antipsychotic drug pimozide promotes apoptosis through the RAF/ERK pathway and enhances autophagy in breast cancer cells

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0032R	anti-Bcl-2 antibody	38356266	3.6	WB	Human
bsm-33036M	anti-β-actin antibody	38356266	3.6	WB	Human

273. miR-10b-5p promotes tumor growth by regulating cell metabolism in liver cancer via targeting SLC38A2

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-12125R	anti-SLC38A2 antibody	38390840	3.6	WB	Mouse, Human

274. Pterostilbene protects against H₂O₂-induced oxidative stress by regulating GAS6/Axl signaling in HL-1 cells

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-7549R	anti-GAS6 antibody	38403920	3.6	WB	Mouse
bs-5180R	anti-Axl antibody	38403920	3.6	WB	Mouse
bs-5181R	anti-p-Axl antibody	38403920	3.6	WB	Mouse

275. LncRNA S100BP promotes proliferation and steroid hormone synthesis of granulosa cells by sponging MiR-2285bc-BMP2 in bovine

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0061R	anti-β-actin antibody	38412119	3.6	WB	Bovine

276. Activation of σ-1 receptor mitigates estrogen withdrawal-induced anxiety/depressive-like behavior in mice via restoration of GABA/glutamate signaling and neuroplasticity in the hippocampus

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-2252R	anti-GABA antibody	38485341	3.5	IF	Mouse

277. Regulation of myo-miR-24-3p on the Myogenesis and Fiber Type Transformation of Skeletal Muscle

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-23102R	anti-KI67 antibody	10.3390/genes15030269	3.5	WB	Mouse

278. Role of diacylglycerol O-acyltransferase 1 (DGAT1) in lipolysis and autophagy of adipose tissue from ketotic dairy cows

Skus	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-6417R	anti-p-PI3K antibody	38395404	3.5	WB	Bovine

bs-10657R	anti-PI3K antibody	38395404	3.5	WB	Bovine
-----------	--------------------	----------	-----	----	--------

279.The HSP90AA1 Gene Is Involved in Heat Stress Responses and Its Functional Genetic Polymorphism Are Associated with Heat Tolerance in Holstein Cows

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-41233R	anti-Cleaved caspase-3 antibody	38395401	3.5	WB	Bovine

280.SNAT2-mediated regulation of estrogen and progesterone in the proliferation of goat mammary epithelial cells

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-12125R	anti-SLC38A2 antibody	38393495	3.5	WB	Goat

281.Icariin alleviates diabetic renal interstitial fibrosis aggravation by inhibiting miR-320a-3p targeting BMP6

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-10423R	anti-Collagen I antibody	38485350	3.5	WB	Human

282.Combination of histological and metabolomic assessments to evaluate the potential pharmacological efficacy of saikosaponin D

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
BA00208	cell counting kit-8 (CCK-8)	38354536	3.4	Other	Other

283.Thymoquinone regulates microglial M1/M2 polarization after cerebral ischemia-reperfusion injury via the TLR4 signaling pathway

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1035R	rabbit anti-CD86 antibody	38325603	3.4	IF	Mouse
bsm-52791R	rabbit anti-CD206 antibody	38325603	3.4	IF	Mouse

284.IFNT-induced IRF1 enhances bovine endometrial receptivity by transactivating LIFR

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1458R	anti-LIFR antibody	38432052	3.4	WB,IHC	Bovine

285.Effect of fullereneol C60 on lung and renal tissue in lower extremity ischemia-reperfusion injury in sevoflurane-treated rats

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0127R	anti-BAX antibody	38334145	3.4	IHC	Rat
bs-0608R	anti-ICAM-1 antibody	38334145	3.4	IHC	Rat
bs-4563R	anti-BCL-2 antibody	38334145	3.4	IHC	Rat

286.M2 macrophage-derived exosomes alleviate KCa3.1 channel expression in rapidly paced HL-1 myocytes via the NF-κB (p65)/STAT3 signaling pathway

SKU	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-6675R	anti-KCa3.1 antibody	38334149	3.4	IF	Mouse

287.PI3K/AKT signaling pathway associates with pyroptosis and inflammation in patients with endometriosis

SKU	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-6417R	anti-phosphor-PI3K antibody	38364342	3.4	WB	Human

288.Immunogenic Cell death Induced by Low Temperature Plasma to Activate Antitumor Immunity and Enhance Therapeutic Efficiency

SKU	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-2386R	anti-CD47 antibody	10.1088/1361-6463/ad2832	3.4	FC	Mouse

289.Abietane diterpenoids with anti-neuroinflammation activity from Rosmarinus officinalis

SKU	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-8000P	Lipopolysaccharides (LPS)	38378134	3.4	Other	Other

290.Therapeutic effect and mechanism of action of pterostilbene nano drugs in dry eye models

SKU	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1035R	anti-CD86 antibody	38387712	3.4	IF	Mouse

291.Development of Drug Efficacy Testing Platform for Glomerulonephritis

SKU	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-10233R-FITC	FITC anti-nephrin antibody	10.3390/mi15030317	3.4	IF	Human
bs-6983R-AF647	AF647 anti-Wilms' tumor 1 (WT1) antibody	10.3390/mi15030317	3.4	IF	Human

292.Characterization of lncRNA and mRNA profiles in ciliary body in experimental myopia

SKU	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1165R	anti-CTNNA1 antibody	38430983	3.4	IHC,WB	Pig
bs-0128R	anti-PIK3R1 antibody	38430983	3.4	IHC,WB	Pig
bs-0486R	anti-ITGB1 antibody	38430983	3.4	IHC,WB	Pig

293. Angelica dahurica extract and its effective component bergapten alleviated hepatic fibrosis by activating FXR signaling pathway

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-10196R	anti- α -SMA antibody	38334900	3.3	WB,IF	Mouse, Rat, Human
bs-0578R	anti-collagen I antibody	38334900	3.3	WB,IF	Mouse, Rat, Human

294. LncRNA HOTTIP regulates TLR4 promoter methylation by recruiting H3K4 methyltransferase MLL1 to affect apoptosis and inflammatory response of fibroblast-like synoviocyte in rheumatoid arthritis

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-8533R	anti-vimentin antibody	38363110	3.3	IF	Mouse

295. Proteomic analysis of multiple organ dysfunction induced by rhabdomyolysis

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-0094R	anti-Serpina3 antibody	38403185	3.3	IHC, WB	Rat

296. RNA-seq reveals role of cell-cycle regulating genes in the pathogenicity of a field very virulent infectious bursal disease virus

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-2994R	anti-IRF7 antibody	38362295	3.2	WB	Chicken

297. Effects of myeloperoxidase on inflammatory responses with hypoxia in Citrobacter rodentium-infectious mice

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-1955R	anti-MCP1 antibody	38415976	3.2	WB, IF	Mouse
bs-1162R	anti-MIP2 antibody	38415976	3.2	WB, IF	Mouse

298. PPM1G regulates hepatic ischemia/reperfusion injury through STING-mediated inflammatory pathways in macrophages

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-2994R	anti-IRF7 antibody	38372470	3.2	WB	Mouse
bs-3196R	anti-p-IRF7 antibody	38372470	3.2	WB	Mouse
bs-9278R	anti-p-IRF3 antibody	38372470	3.2	IF	Mouse

299. AgPd NPs-based Lateral flow immunoassay for simultaneous detection of glycocholic acid and alpha-fetoprotein

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-0296G	Goat anti-mouse antibody	38372146	3.1	Other	Human

300. Novel Prognostic Model Construction of Tongue Squamous Cell Carcinoma Based on Apigenin-Associated Genes

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
-----	--------------	-----------	--------------	--------------	------------------

BA00208	Cell Counting Kit-8 (CCK-8)	31345004	3.1	Other	Other
---------	-----------------------------	----------	-----	-------	-------

301. Electroacupuncture Improving Intestinal Barrier Function in Rats with Irritable Bowel Syndrome Through Regulating Aquaporins

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1253R	anti-AQP3 antibody	38421507	3.1	WB	Rat
bs-6786R	anti-AQP8 antibody	38421507	3.1	IHC	Rat

302. TAK-242 improves sepsis-associated acute kidney injury in rats by inhibiting the TLR4/NF- κ B signaling pathway

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
AK288	malondialdehyde (MDA) Assay Kit	38482886	3	Other	Other
AK060	superoxide dismutase (SOD) Assay Kit	38482886	3	Other	Other

303. Alpha-fetoprotein upregulates hepatocellular carcinoma cell-intrinsic PD-1 expression through the LATS2/YAP/TEAD1 pathway

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-4708R	anti-AFP antibody	38395204	3	IF	Human

304. Curcumin attenuates aflatoxin B1-induced ileum injury in ducks by inhibiting NLRP3 inflammasome and regulating TLR4/NF- κ B signaling pathway

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-60331R	anti-mucin 2 (MUC2) antibody	38400893	3	IF	Duck
bsm-33207M	anti-TNF- α antibody	38400893	3	WB	Duck
bs-15581R	anti-p-I κ B antibody	38400893	3	WB	Duck
bs-0982R	anti-p-NF- κ B antibody	38400893	3	WB	Duck
bs-41373R	anti-GAPDH antibody	38400893	3	WB	Duck

305. Inhibiting Liver Autophagy and Promoting Hepatocyte Apoptosis by Schistosoma Japonicum Infection

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0549R	Rabbit anti-collagen III a1 antibody	38393131	2.9	WB	Mouse, Human

306. Audiogenic kindling activates glutamatergic system in the hippocampus of rats with genetic predisposition to audiogenic seizures

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-11365R	Anti-SV2B Polyclonal Antibody	38325559	2.9	WB	Rat

307. BMSCs attenuate radiation-induced brain injury induced hippocampal neuronal apoptosis through a PI3K/Akt/Bax/Bcl-2 signaling pathway

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-0061R	anti- β -actin antibody	38331376	2.9	WB	Rat
bs-6417R	anti-P-PI3K antibody	38331376	2.9	WB	Rat

308. NGF Signaling Exacerbates KOA Peripheral Hyperalgesia via the Increased TRPV1-Labeled Synovial Sensory Innervation in KOA Rats

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bsm-33036M	anti- β -actin antibody	38410126	2.9	WB	Rat
bs-10806R	anti-NGF antibody	38410126	2.9	WB,IF	Rat
bs-0193R	anti-TrKA antibody	38410126	2.9	WB,IHC,IF	Rat
bs-23926R	anti-TRPV1 antibody	38410126	2.9	WB,IHC,IF	Rat
bs-0812R	anti-IL-1 β antibody	38410126	2.9	WB	Rat
bs-3806R	anti-PGP 9.5 antibody	38410126	2.9	WB,IHC,IF	Rat

309. Evaluation of PDL1 positive cancer cell-specific binding activity of recombinant anti-PDL1 scFv

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-41373R	Anti-GAPDH antibody	38377106	2.9	WB	Human

310. The changes of neurogenesis in the hippocampal dentate gyrus of SAMP8 mice and the effects of acupuncture and moxibustion

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-20797R	anti-DCX antibody	38395250	2.9	IF	Mouse

311. SIRT2 regulates apoptosis by inducing mitophagy in sheep cumulus cells

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-0295G-HRP	horseradish peroxidase-linked goat anti-rabbit IgG antibody	38330860	2.8	WB	Sheep

312. Application of sperm motion kinematics and motility-related proteins for prediction of male fertility

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-24834R	anti-RSPH9 antibody	38359560	2.8	ELISA	Pig

313. CTC together with Shh and Nrf2 are prospective diagnostic markers for HNSCC

SKU	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
-----	--------------	-----------	--------------	--------------	------------------

bs-0593R	Anti-Epithelial Cell Adhesion Molecule Antibody	38336617	2.8	Other	Human
bs-4819R	Anti-CD45 Antibody	38336617	2.8	IF	Human

314. Erythropoietin-producing hepatocellular receptor B6 is highly expressed in non-functioning pituitary neuroendocrine tumors and its expression correlates with tumor size

Skus	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-25280R	anti-EphB6 antibody	38341842	2.8	IHC	Human

315. Reduced Proliferative Capacity and Defense against Staphylococcus aureus in Human Nasal Mucosal Epithelium Lacking ZNF365

Skus	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-10009R	anti-E-cadherin antibody	38354721	2.8	WB	Human
bs-0756R	anti-vimentin antibody	38354721	2.8	WB	Human
bs-1371R	anti-Snai1 antibody	38354721	2.8	WB	Human

316. A mutation in LPAR2 activates the miR-939-5p-LPAR2-PI3K/AKT axis to regulate the proliferation and apoptosis of granulosa cells in sheep

Skus	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-10368R	anti-LPAR2 antibody	38368704	2.8	WB	Sheep
bs-20773R	anti-Caspase 9 antibody	38368704	2.8	WB	Sheep
bsm-33411M	anti-Bcl-2 antibody	38368704	2.8	WB	Sheep

317. Effect of estradiol after bacterial infection on the Wnt/ β -catenin pathway in bovine endometrium epithelial cells and organoids

Skus	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-0295G-FITC	FITC-labeled goat anti-rabbit secondary antibody	38402700	2.8	IF	Bovine

318. Identification of Pyroptosis-Related Genes Regulating the Progression of Chronic Rhinosinusitis with Nasal Polyps

Skus	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-5986R	anti-AIM2 antibody	38402873	2.8	WB, IHC	Human

319. Protective Effects of Myricetin and Morin on Neurological Damage in A β 1-42/A β 3+ -induced Alzheimer's Disease Model of Rats

Skus	Product Name	Pubmed id	Journal_tier	Applications	Cross Reactivity
bs-4539R	anti-Interleukins-6 (IL-6) antibody	38423257	2.8	WB	Rat

320. Mogroside V Inhibits M1 Polarization and Inflammation of Diabetic Mouse Macrophages via p38 MAPK/NF-Kb Signaling Pathway

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1035R	anti-CD86 antibody	38415803	2.8	IF	Mouse
bsm-33117M	anti-P65 antibody	38415803	2.8	WB	Mouse
bs-0982R	anti-p-p65 antibody	38415803	2.8	WB	Mouse

321. Fluorofenidone attenuates renal fibrosis by inhibiting lysosomal cathepsin-mediated NLRP3 inflammasome activation

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-6741R	anti-ASC antibody	38476910	2.7	IF, WB	rat, Mouse, Human

322. Antibiotic administration aggravates asthma by disrupting gut microbiota and the intestinal mucosal barrier in an asthma mouse model

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-10966R	rabbit anti-actin antibody	38476896	2.7	WB	Mouse

323. miR-324-3p Suppresses Hepatic Stellate Cell Activation and Hepatic Fibrosis Via Regulating SMAD4 Signaling Pathway

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bsm-33036M	anti-β-actin antibody	38407690	2.6	WB	Mouse, Human
bsm-33187M	anti-SMA antibody	38407690	2.6	WB	Mouse, Human
bsm-33170M	anti-Vimentin antibody	38407690	2.6	WB	Mouse, Human

324. Role of increased IGFBP2 in trophoblast cell proliferation and recurrent spontaneous abortion development: A pilot study

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1108R	rabbit anti-IGFBP2 antibody	38316422	2.5	IHC	Human

325. CircRNA-5335 Regulates the Differentiation and Proliferation of Sheep Preadipocyte via the miR-125a-3p/STAT3 Pathway

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1658R	anti-P-STAT3 polyclonal antibody	38393088	2.4	WB	Lamb

326. Apolipoprotein C-III amyloidosis in white lions (Panthera leo)

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-4741R	anti-human ApoC-III antibody	38345009	2.4	IHC	Lion

327. Ameliorative Effect of Macadamia Nut Protein Peptides on Acetaminophen-Induced Acute Liver Injury in Mice

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bsm-33117M	anti-NF- κ B p65 antibody	38386536	2.4	IHC	Mouse
bs-0812R	anti-IL-1 β antibody	38386536	2.4	IF	Mouse

328. In Vitro Study of Thymosin Beta 4 Promoting Transplanted Fat Survival by Regulating Adipose-Derived Stem Cells

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
Y-0266	Thymosin beta 4 (T β 4)	38409346	2.4	Other	Other

329. Short-term inhalation exposure to cigarette smoke induces oxidative stress and inflammation in lungs without systemic oxidative stress in mice

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1278R	anti-8-OH-dG antibody	10.1007/s43188-023-00223-y	2.3	IHC	Mouse

330. Prognostic biomarker CPEB3 and its associations with immune infiltration in clear cell renal cell carcinoma

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
C01-03001	normal goat serum (10%)	38476610	2.3	Other	Other

331. Inhibition of anlotinib-induced autophagy attenuates invasion and migration by regulating epithelial-mesenchymal transition and cytoskeletal rearrangement through ATG5 in human osteosarcoma cells

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0756R	anti-rabbit Vimentin antibody	38381883	2.3	WB	Human

332. PLCD3 inhibits apoptosis and promotes proliferation, invasion and migration in gastric cancer

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-12713R	anti-PLCD3 antibody	38305998	2.2	IF	Human

333. Elucidating the mechanisms of formononetin in modulating atherosclerotic plaque formation in ApoE^{-/-} mice

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-1049R	rabbit anti- α 7nAChR antibody	38388385	2.1	WB	Mouse

334. The effect of bta-miR-1296 on the proliferation and extracellular matrix synthesis of bovine mammary fibroblasts

Skus	ProductNames	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0189R	rabbit anti-alpha smooth muscle actin antibody	38409638	2.1	WB	Bovine
bs-0578R	rabbit anti-collagen I antibody	38409638	2.1	WB	Bovine
bs-6951R	rabbit anti-AKT1 + 2 + 3 antibody	38409638	2.1	WB	Bovine
bs-0876R	rabbit anti-phospho-AKT antibody	38409638	2.1	WB	Bovine

335. Adiponectin attenuates H₂O₂-induced apoptosis in chicken skeletal myoblasts through the lysosomal-mitochondrial axis

Skus	ProductNames	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-12933R	anti-chymotrypsin B antibody	38427138	2.1	WB	Chicken

336. In ovo injection of AZD6244 suppresses feather follicle development by the inhibition of ERK and Wnt/ β -catenin pathways in goose embryos (*Anser cygnoides*)

Skus	ProductNames	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-23663R	rabbit anti- β -catenin antibody	38393940	2	IHC, WB	Goose
bs-1280R	rabbit anti-TCF4 antibody	38393940	2	WB	Goose
bs-1843R	rabbit anti-LEF1 antibody	38393940	2	WB	Goose

337. Exercise before pregnancy exerts protective effect on prenatal stress-induced impairment of memory, neurogenesis, and mitochondrial function in offspring

Skus	ProductNames	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0248R	rabbit brain-derived neurotrophic factor (BDNF) antibody	38433854	1.9	WB	Mouse

338. Protective mechanisms of *Juncus effusus* and Carbonized *Juncus effusus* against d-galactosamine-induced acute liver injury in mice

Skus	ProductNames	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-2184R	anti-NQO1 antibody	38325836	1.7	WB	Mouse
bsm-33036M	anti- β -Actin antibody	38325836	1.7	WB	Mouse

339. Protective effect of exosomes derived from bone marrow mesenchymal stem cells on hypoxia reperfusion injury of cardiomyocytes

Skus	ProductNames	Pubmed id	Journal_tier	Applications	CrossReactivity
bsm-33199M	anti-cleaved caspase 3 antibody	38430039	1.6	WB	Rat

bsm-52304R	anti-Bcl-2 antibody	38430039	1.6	WB	Rat
------------	---------------------	----------	-----	----	-----

340. Human papillomavirus type 16 E7 promotes cell viability and migration in cervical cancer by regulating the miR-23a/HOXC8 axis

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-49102R	anti-HPV16 E6 antibody	38348790	1.3	WB	Human
bs-10202R	anti-HPV16 E7 antibody	38348790	1.3	WB	Human
bs-0394R	anti-HOXC8 antibody	38348790	1.3	WB	Human
bs-41373R	anti-GAPDH antibody	38348790	1.3	WB	Human
bsm-0295M-HRP	HRP labelled mouse anti-rabbit IgG antibody	38348790	1.3	WB	Human

341. Cobalt Chloride as a Hypoxia Mimicking Agent Induced HIF-1 α and mTOR Expressions of Human Umbilical Cord Mesenchymal Stem Cells

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-0737R-FITC	anti-FITC HIF-1 α antibody	10.1590/pboci.2024.020	0.6	IF	Human
bs-1992R-FITC	anti-FITC mTOR antibody	10.1590/pboci.2024.020	0.6	IF	Human

342. Catalpol attenuates EMT by inhibiting Wnt/ β -catenin and TGF- β /Smads signaling to alleviate kidney fibrosis

Sku	ProductName	Pubmed id	Journal_tier	Applications	CrossReactivity
bs-3484R	Rabbit anti-Smad3 polyclonal antibody	10.4314/tjpr.v23i1.8	0.6	WB	Rat