

bs-0296R-BF594

[Secondary Antibodies]

Rabbit Anti-Mouse IgG H&L, BF594 conjugated



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: IF (1:100-1000)
Clonality: Polyclonal		Flow-Cyt (1:100-1000)
Target: Rabbit Anti-Mouse IgG H&L		ICC/IF (1:100-1000)
Purification: affinity purified by Protein A		Excitation Spectrum: 594nm
Concentration: 2.0 mg/ml		Emission spectrum: 621nm
Storage: 10 mM TBS (pH=7.4) with 1% BSA, 0.03% Proclin300 and 50% glycerol. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		Reactivity: Mouse

— SELECTED CITATIONS —

- **[IF=8.7]** Xue Ou. et al. An MXene nanocomposite hydrogel for enhanced diabetic infected wound healing via photothermal antibacterial properties and bioactive molecule integration. Materials Today Bio. 2025 Feb 8;31:101538. IF ;Mouse. 10.1016/j.mtbio.2025.101538