

bs-0312R-FITC

[Secondary Antibodies]



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

Rabbit Anti-Avidin, FITC conjugated

— DATASHEET —

Host: Rabbit Clonality: Polyclonal Target: Rabbit Anti-Avidin Purification: affinity purified by Protein A Concentration: 2.0 mg/ml Storage: 10 mM TBS (pH=7.4) with 1% BSA, 0.03% Proclin300 and 50% glycerol. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. Background: The biological function of avidin is not known. It is a homotetrameric glycoprotein found in chicken egg white, and forms a strong non covalent specific complex with biotin (one molecule of biotin per subunit of avidin).	Isotype: IgG Applications: IF (1:100-1000) Flow-Cyt (1:100-1000) ICC/IF (1:100-1000) Excitation Spectrum: 495nm Emission spectrum: 519nm Reactivity: Avidin
--	---

— SELECTED CITATIONS —

- **[IF=2.4]** Yanru Chen. et al. Melatonin alleviates oxidative stress and inflammation of Leydig cells of Min pig through SIRT1 pathway. THERIOGENOLOGY. 2025 Feb;233:112 IF ;Pig. 39613495