

bs-10298R

[Primary Antibody]

H7N9 NA Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit Clonality: Polyclonal Target: H7N9 NA Immunogen: KLH conjugated synthetic peptide derived from H7N9 neuraminidase: 1-100/465. Purification: affinity purified by Protein A Concentration: 1mg/ml Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.	Isotype: IgG Applications: ELISA (1:5000-10000) Reactivity: (predicted: Chicken, H7N9) Predicted MW.: 51 kDa Subcellular Location: Cytoplasm ,Nucleus
--	--