

bs-41384R

[Primary Antibody]

ASFV p72 Rabbit pAb

Bioss
ANTIBODIES

www.bioss.com.cn

sales@bioss.com.cn


techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: WB (1:500-2000) ELISA (1:5000-10000)
Clonality: Polyclonal		Reactivity: ASFV
Target: ASFV p72		
Immunogen: Recombined ASFV p72 protein: 1-646/646.		
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		Predicted MW.: 73.1 kDa
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		Observed MW.: 73.1 kDa

— VALIDATION IMAGES —



Sample: Lane 1: Recombined ASFV p72 protein
Primary: Anti-ASFV p72 (bs-41384R) at 1/1000
dilution Secondary: IRDye800CW Goat Anti-
Rabbit IgG at 1/20000 dilution Predicted band
size: 73.1 kDa Observed band size: 73.1 kDa