

bs-0368G-FITC**[Secondary Antibodies]**

www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

Goat Anti-Mouse IgM, FITC conjugated**— DATASHEET —**

Host: Goat Clonality: Polyclonal Target: Goat Anti-Mouse IgM Purification: affinity purified by Protein G Concentration: 2.0 mg/ml Storage: 10 mM TBS (pH=7.4) with 1% BSA, 0.03% Proclin300 and 50% glycerol. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. Background: Immunoglobulin M (IgM) normally constitutes about 10% of serum immunoglobulins. IgM antibody is prominent in early immune responses to most antigens and is largely confined to plasma due to its large size. Monomeric IgM is expressed as a membrane bound antibody on the surface of B cells and as a pentamer when secreted by plasma cells. Due to its high valency IgM is more efficient than other isotypes in binding antigens with repeating epitopes (virus particles and red blood cells) and is more efficient than IgG in activating the complement pathway. The gene for the mu constant region contains four domains separated by short intervening sequences.	Isotype: IgG Applications: IF (1:100-1000) Flow-Cyt (1:100-1000) ICC/IF (1:100-1000) Excitation Spectrum: 495nm Emission spectrum: 519nm Reactivity: Mouse
--	---

— SELECTED CITATIONS —

- **[IF=14.7]** Cui Zhanding. et al. High-throughput screening unveils nitazoxanide as a potent PRRSV inhibitor by targeting NMRAL1. NAT COMMUN. 2024 Jun;15(1):1-12 IFA,IHC ;Pig,Monkey. 38844461
- **[IF=5.924]** Hongtao Liu. et al. Antiviral Effects of ABMA and DABMA against Influenza Virus In Vitro and In Vivo via Regulating the Endolysosomal Pathway and Autophagy. INT J MOL SCI. Int J Mol Sci. 2022 Jan;23(7):3940 IF ;Dog. 35409297
- **[IF=5.895]** Guangning Kou. et al. Sesamin Activates Skeletal Muscle FNDC5 Expression and Increases Irisin Secretion via the SIRT1 Signaling Pathway. J AGR FOOD CHEM. 2022;XXXX(XXX):XXX-XXX IF ;Mouse. 35708276
- **[IF=5.1]** Zhanding Cui. et al. Formononetin and mizoribine inhibit Porcine Reproductive and Respiratory Syndrome Virus replication in vitro. FRONT NUTR. 2025 Mar;11: IF ;. 40196742
- **[IF=5.005]** Li Y et al. Exosomes derived from Toxoplasma gondii stimulate an inflammatory response through JNK signaling pathway. Nanomedicine (Lond). 2018 May;13(10):1157-1168. IF ;Mouse. 29542367