

bsm-0294M**[Secondary Antibodies]**

www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

Mouse Anti-Goat IgG mAb**— DATASHEET —**

Host: Mouse	Isotype: IgG1	Applications: Isotype Control Blocking Assay etc. Conjugate-Dependent.
Clonality: Monoclonal	CloneNo.: 3B4	
Target: Mouse Anti-Goat IgG mAb		Reactivity: Goat
Immunogen: Native Goat IgG: full length.		
Purification: affinity purified by Protein G		
Storage: 0.01M PBS (pH7.4). Store at -20°C stable for 2 years (lyophilized powder). Avoid repeated freeze/thaw cycles.		
Background: Immunoglobulin G (IgG), is one of the most abundant proteins in serum with normal levels between 8-17 mg/mL in adult blood. IgG is important for our defence against microorganisms and the molecules are produced by B lymphocytes as a part of our adaptive immune response. The IgG molecule has two separate functions; to bind to the pathogen that elicited the response and to recruit other cells and molecules to destroy the antigen. The variability of the IgG pool is generated by somatic recombination and the number of specificities in an individual at a given time point is estimated to be 1011 variants.		

— SELECTED CITATIONS —

- **[IF=10.5]** Yuanfen Gao. et al. UV-B Responsive Flavonoid Synthesis Contributes to Tartary Buckwheat High-Altitude Adaption. PLANT BIOTECHNOL J. 2025 Aug;; WB ;Arabidopsis thaliana. 40783815