

bs-20024R**[Primary Antibody]****phospho-SMAD7 (Ser249) Rabbit pAb****Bioss**
ANTIBODIES

www.bioss.com.cn


sales@bioss.com.cn

techsupport@bioss.com.cn


400-901-9800

— DATASHEET —


Host: Rabbit Clonality: Polyclonal GeneID: 4092 Target: SMAD7 (Ser249) Immunogen: KLH conjugated synthesised phosphopeptide derived from human SMAD7 around the phosphorylation site of Ser249: SD(p-S)QL. Purification: affinity purified by Protein A Concentration: 1mg/ml Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. Background: The protein encoded by this gene is a nuclear protein that binds the E3 ubiquitin ligase SMURF2. Upon binding, this complex translocates to the cytoplasm, where it interacts with TGF-beta receptor type-1 (TGFBR1), leading to the degradation of both the encoded protein and TGFBR1. Expression of this gene is induced by TGFBR1. Variations in this gene are a cause of susceptibility to colorectal cancer type 3 (CRCS3). Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jun 2010]	Isotype: IgG SWISS: O15105	Applications: WB (1:500-2000) IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) Reactivity: Human, Mouse, Rat (predicted: Pig, Cow, GuineaPig) Predicted MW.: 46 kDa Observed MW.: 50 kDa Subcellular Location: Cytoplasm ,Nucleus
--	---	--

— VALIDATION IMAGES —

Sample: Lane 1: Mouse Cerebrium tissue lysates
Lane 2: Rat Cerebrium tissue lysates
Primary: Anti-Phospho-SMAD7 (Ser249) (bs-20024R) at 1/1000 dilution
Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution
Predicted band size: 46 kDa
Observed band size: 50 kDa



Paraformaldehyde-fixed, paraffin embedded (human brain); Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15min; Block endogenous peroxidase by 3% hydrogen peroxide for 20 minutes; Blocking buffer (normal goat serum) at 37°C for 30min; Antibody incubation with (Phospho-SMAD7 (Ser249)) Polyclonal Antibody, Unconjugated (bs-20024R) at 1:200 overnight at 4°C, followed by operating according to SP Kit(Rabbit) (sp-0023) instructions and DAB staining.



Paraformaldehyde-fixed, paraffin embedded (rat stomach); Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15min; Block endogenous peroxidase by 3% hydrogen peroxide for 20 minutes; Blocking buffer (normal goat serum) at 37°C for 30min; Antibody incubation with (Phospho-SMAD7 (Ser249)) Polyclonal Antibody, Unconjugated (bs-20024R) at 1:200 overnight at 4°C, followed by operating according to SP Kit(Rabbit) (sp-0023) instructions and DAB staining.

— SELECTED CITATIONS —

- **[IF=3.8]** Huang Lan. et al. Unveiling poly(rC)-binding protein 2 as the target protein for curcucosone C against prostate cancer: mechanism validation through click chemistry-activity based proteomics profiling approach. BMC CANCER. 2023 Dec;23(1):1-14 IHC, WB ; Mouse, Human. 37814239
- **[IF=3]** Amei Tang. et al. Tuberostemonine may alleviates proliferation of lung fibroblasts caused by pulmonary fibrosis. INT J IMMUNOPATH PH. ;(): WB ; Mouse. 39140804
- **[IF=2.2]** Yan Hong. et al. Paridis Rhizoma Saponins ameliorate CCl4-induced liver fibrosis via modulation of the Nrf2-

Important Note: This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

HO-1/NQO-1 and TGF- β 1/Smads signaling pathways. Biochemistry and Biophysics Reports. 2025 Jun;42:102040 WB,IHC
;Rat. 40458200