
Donkey Anti-Goat IgG H&L, Cy5 conjugated

Catalog Number: bs-0294D-Cy5

Target Protein: Donkey Anti-Goat IgG H&L

Concentration: 2.0 mg/ml

Form: Liquid

Host: Donkey

Clonality: Polyclonal

Isotype: IgG

Applications: Flow-Cyt (1:100-1000), IF (1:100-1000)

Excitation spectrum: 649nm

Emission spectrum: 670nm

Reactivity: Goat

Purification: affinity purified by Protein G

Storage: 10 mM TBS (pH=7.4) with 1% BSA, 0.03% Proclin300 and 50% glycerol.

Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: Immunoglobulin G (IgG), is one of the most abundant proteins in serum with normal levels between 8-17 mg/mL in adult blood. IgG is important for our defence against microorganisms and the molecules are produced by B lymphocytes as a part of our adaptive immune response. The IgG molecule has two separate functions; to bind to the pathogen that elicited the response and to recruit other cells and molecules to destroy the antigen. The variability of the IgG pool is generated by somatic recombination and the number of specificities in an individual at a given time point is estimated to be 1011 variants.

PRODUCT SPECIFIC PUBLICATIONS

[IF=5.1] Zhu Yuan-Zheng. et al. Genome-wide search links senescence-associated secretory proteins with susceptibility for coronary artery disease in mouse and human. J GERONTOL A-BIOL. 2024 Feb;; IF ; Mouse . 38416803