

Goat Anti-Rabbit IgG H&L, PE-CY5 conjugated

Catalog Number: bs-0295G-PE-CY5

Target Protein: Goat Anti-Rabbit IgG H&L

Concentration: 2.0 mg/ml

Form: Liquid

Host: Goat

Clonality: Polyclonal

Isotype: IgG

Applications: IF (1:100-1000), Flow-Cyt (1:100-1000)

Excitation spectrum: 488nm,561nm,649nm

Emission spectrum: 667nm

Reactivity: Rabbit


Purification: affinity purified by Protein G, nonspecific adsorbed

Storage: 10 mM TBS (pH=7.4) with 1% BSA, 0.03% Proclin300 and 50% glycerol.

Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: Immunoglobulin G (IgG), is one of the most abundant proteins in serum with normal levels between 8-17 mg/mL in adult blood. IgG is important for our defence against microorganisms and the molecules are produced by B lymphocytes as a part of our adaptive immune response. The IgG molecule has two separate functions; to bind to the pathogen that elicited the response and to recruit other cells and molecules to destroy the antigen. The variability of the IgG pool is generated by somatic recombination and the number of specificities in an individual at a given time point is estimated to be 10¹¹ variants.

VALIDATION IMAGES



25 ug total protein per lane of various lysates (see on figure) probed with Rabbit IgG H&L polyclonal antibody, unconjugated (bs-0295G) at 1:1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r.t. for 60 min.