

bsm-30051M**[Primary Antibody]****human CD71 Mouse mAb****BioSS**
ANTIBODIES

www.bioss.com.cn


sales@bioss.com.cn

techsupport@bioss.com.cn


400-901-9800

— DATASHEET —

Host: Mouse Clonality: Monoclonal GeneID: 7037 Target: human CD71 Purification: affinity purified by Protein G Concentration: 1mg/ml Storage: Size : 100ug 0.01M TBS (pH7.4). Size : 200ug (PBS only) 0.01M PBS Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. Background: This gene encodes a cell surface receptor necessary for cellular iron uptake by the process of receptor-mediated endocytosis. This receptor is required for erythropoiesis and neurologic development. Multiple alternatively spliced variants have been identified. [provided by RefSeq, Sep 2015]	Isotype: Mouse IgG1, k CloneNo.: 7B4 SWISS: P02786	Applications: Flow-Cyt (1ug/test) ICC/IF (1:100-500) Reactivity: Human Predicted MW.: 74/85 kDa Observed MW.: 95 kDa Subcellular Location: Secreted ,Cell membrane
--	---	--

— VALIDATION IMAGES —

Hela cell; 4% Paraformaldehyde-fixed; Triton X-100 at room temperature for 20 min; Blocking buffer (normal goat serum, C-0005) at 37°C for 20 min; Antibody incubation with (CD71) monoclonal Antibody, Unconjugated (bsm-30051M) 1:100, 90 minutes at 37°C; followed by a conjugated Goat Anti-Mouse IgG antibody at 37°C for 90 minutes, DAPI (blue, C02-04002) was used to stain the cell nuclei.



Blank control:K562. Primary Antibody (green line): Mouse Anti-human CD71 antibody (bsm-30051M) Dilution: 1ug/Test; Secondary Antibody (white blue line) : Goat anti-Mouse IgG-BF488 Dilution: 0.5ug/Test. Isotype control (orange line) : Normal Rabbit IgG Protocol
The cells were incubated in 5%BSA to block non-specific protein-protein interactions for 30 min at room temperature .Cells stained with Primary Antibody for 30 min at room temperature. The secondary antibody used for 40 min at room temperature. Acquisition of 20,000 events was performed.

— SELECTED CITATIONS —

- **[IF=1.18]** Ma Caiyun. et al. In Vitro Culture and Multipotency Evaluation of Broiler Umbilical Cord Mesenchymal Stem Cells. BRAZ ARCH BIOL TECHN. 2023 Feb;66:e23220031 IF,FCM ;Human. 10.1590/1678-4324-2023220031