

DATASHEET

Mouse Anti-Bovine IgG H&L / PE-Cy7 **PE-Cy7标记的小鼠抗牛IgG H&L**

Catalog Number: bs-0326M-PE-Cy7

Target Protein: Mouse Anti-Bovine IgG H&L / PE-Cy7

Concentration: 2.0 mg/ml

Storage: 10 mM TBS (pH=7.4) with 1% BSA, 0.03% Proclin300 and 50% glycerol.

Background:

VALIDATION IMAGES

PRODUCT SPECIFIC PUBLICATIONS

Important Note: This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.