

bs-0887R

[Primary Antibody]

MeV H Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit Clonality: Polyclonal Target: MeV H Immunogen: KLH conjugated synthetic peptide derived from Measles virus hemagglutinin protein: 101-200/617. Purification: affinity purified by Protein A Concentration: 1mg/ml Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. Background: No data available.	Isotype: IgG Applications: ELISA (1:5000-10000) Reactivity: (predicted: Measles virus) Predicted MW.: 68 kDa Subcellular Location: Nucleus
--	---