

bs-0336R-HRP

[Secondary Antibodies]

Rabbit Anti-Monkey IgM, HRP conjugated

Bioss
ANTIBODIES

www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: WB (1:1000-10000)
Clonality: Polyclonal		IHC-P (1:100-500)
Target: Rabbit Anti-Monkey IgM		IHC-F (1:100-1000)
Purification: affinity purified by Protein A		ELISA (1:1000-10000)
Concentration: 2.0 mg/ml		Reactivity: Monkey
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C. Avoid repeated freeze/thaw cycles.		
Background: IgM normally constitutes about 10% of serum immunoglobulins. IgM antibody is prominent in early immune responses to most antigens and is largely confined to plasma due to its large size. Monomeric IgM is expressed as a membrane bound antibody on the surface of B cells and as a pentamer when secreted by plasma cells. Due to its high valency IgM is more efficient than other isotypes in binding antigens with repeating epitopes (virus particles and red blood cells) and is more efficient than IgG in activating the complement pathway. The gene for the mu constant region contains four domains separated by short intervening sequences.		

— SELECTED CITATIONS —

- **[IF=16.806]** Nan Gao. et al. Development of Neutralization Breadth against Diverse HIV-1 by Increasing Ab-Ag Interface on V2. 2022 Mar 23 ELISA ;Gp120 Monomer,Bg505 Trimer. 35319830