

bs-0326R-HRP**[Secondary Antibodies]****Rabbit Anti-Bovine IgG H&L, HRP conjugated****BioSS**
ANTIBODIES

www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: WB (1:1000-10000)
Clonality: Polyclonal		IHC-P (1:100-500)
Target: Rabbit Anti-Bovine IgG H&L		IHC-F (1:100-1000)
Purification: affinity purified by Protein A		ELISA (1:1000-10000)
Concentration: 2.0 mg/ml		Reactivity: Bovine
Storage: 10 mM TBS (pH=7.4) with 1% BSA, 0.03% Proclin300 and 50% glycerol. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		
Background: Immunoglobulin G (IgG), is one of the most abundant proteins in serum with normal levels between 8-17 mg/mL in adult blood. IgG is important for our defence against microorganisms and the molecules are produced by B lymphocytes as a part of our adaptive immune response. The IgG molecule has two separate functions; to bind to the pathogen that elicited the response and to recruit other cells and molecules to destroy the antigen. The variability of the IgG pool is generated by somatic recombination and the number of specificities in an individual at a given time point is estimated to be 1011 variants.		

— SELECTED CITATIONS —

- **[IF=3.471]** Tongsheng Qi. et al. Seroepidemiology of Neosporosis in Various Animals in the Qinghai-Tibetan Plateau. FRONT VET SCI. 2022 Jul 19;9:953380 ELISA ;Neospora caninum. 35928116
- **[IF=3.231]** Jinchao Zhang. et al. Serological Analysis of IgG and IgM Antibodies against Anaplasma spp. in Various Animal Species of the Qinghai-Tibetan Plateau. ANIMALS. 2022 Jan;12(19):2723 ELISA ;Bovine. 36230463
- **[IF=2.289]** Zhen Han. et al. Identification of a novel variant erythrocyte surface antigen-1 (VESA1) in Babesia orientalis. 2021 Jul 05 WB ;Mouse. 34219188