

bsm-54451R**[Primary Antibody]****BioSS**
ANTIBODIES

www.bioss.com.cn


sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

Fbx32 Recombinant Rabbit mAb**— DATASHEET —**

Host: Rabbit Clonality: Recombinant GeneID: 114907 Target: Fbx32 Immunogen: A synthesized peptide derived from human FBXO32: 20-84. Purification: affinity purified by Protein A Concentration: 1mg/ml Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4°C. Background: Fbx32 is an E3 ubiquitin ligase that initiates ATP dependent ubiquitin-mediated proteolysis and promotes muscle atrophy. It is highly expressed during muscle atrophy, whereas mice deficient in this gene were found to be resistant to atrophy. It is also thought to recognize and bind to some phosphorylated proteins and promote their ubiquitination and degradation during skeletal muscle atrophy. Fbx32 interacts with MyoD by ubiquitination via a sequence found in transcriptional coactivators and therefore may play an important role in the course of muscle differentiation by determining the abundance of MyoD.	Isotype: IgG CloneNo.: 2A3 SWISS: Q969P5	Applications: WB (1:1000-2000) ICC/IF (1:50-200) Reactivity: Mouse, Rat (predicted: Human) Predicted MW.: 42 kDa Observed MW.: 42 kDa Subcellular Location: Cytoplasm ,Nucleus
---	---	---

— VALIDATION IMAGES —

25 ug total protein per lane of various lysates (see on figure) probed with Fbx32 monoclonal antibody, unconjugated (bsm-54451R) at 1:10000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r.t. for 60 min.

— SELECTED CITATIONS —

- **[IF=8.713]** Zhao-Bo Luo. et al. Fecal transplant from myostatin deletion pigs positively impacts the gut-muscle axis. ELIFE. 2023; 12: e81858 WB ;Mouse. 37039469
- **[IF=6.228]** Xuzhu Lin. et al. Undercarboxylated osteocalcin and ibandronate combination ameliorates hindlimb immobilization-induced muscle wasting. J PHYSIOL-LONDON. 2023 Apr;; WB ;Human. 36999349
- **[IF=5.085]** Shuang-Yan Chang. et al. miR-320 regulates myogenesis by targeting growth factor receptor-bound protein-2 and ameliorates myotubes atrophy. INT J BIOCHEM CELL B. 2022 Jun;147:106212 WB ;Mouse. 10.1016/j.biocel.2022.106212
- **[IF=4.755]** Zhang Lijiao. et al. Dysregulated myokines and signaling pathways in skeletal muscle dysfunction in a

Important Note: This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

cigarette smoke-induced model of chronic obstructive pulmonary disease. FRONT PHYSIOL. 2022 Aug;0:1600 WB ;Mouse. 36091368

- **[IF=4.7]** Danyang Li. et al. IL-36R deletion mitigates cigarette smoke-induced airway inflammation and skeletal muscle dysfunction. INT IMMUNOPHARMACOL. 2025 Oct;164:115317 WB ;Mouse. 40773893