

bs-2001R**[Primary Antibody]****H1N1 Hemagglutinin 1 Rabbit pAb****Bioss**
ANTIBODIES

www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: ELISA (1:5000-10000)
Clonality: Polyclonal		Reactivity: (predicted: Influenza A virus H1N1)
Target: H1N1 Hemagglutinin 1		Predicted MW.: 63 kDa
Immunogen: KLH conjugated synthetic peptide derived from Influenza A Virus Hemagglutinin: 1-100/566.		Subcellular Location: Cell membrane ,Cytoplasm
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		
Background: Influenza A virus is a major public health threat. Novel influenza virus strains caused by genetic drift and viral recombination emerge periodically to which humans have little or no immunity, resulting in devastating pandemics. Influenza A can exist in a variety of animals; however it is in birds that all subtypes can be found. These subtypes are classified based on the combination of the virus coat glycoproteins hemagglutinin (HA) and neuraminidase (NA) subtypes. During 1997, an H5N1 avian influenza virus was determined to be the cause of death in 6 of 18 infected patients in Hong Kong. There was some evidence of human to human spread of this virus, but it is thought that the transmission efficiency was fairly low. HA interacts with cell surface proteins containing oligosaccharides with terminal sialyl residues. Virus isolated from a human infected with the H5N1 strain in 1997 could bind to oligosaccharides from human as well as avian sources, indicating its species jumping ability.		

— SELECTED CITATIONS —

- **[IF=5.44]** Li, Yanxia, et al. "Highly sensitive fluorescent immunosensor for detection of influenza virus based on Ag autocatalysis." Biosensors and Bioelectronics (2013). ELISA ;Virus. 24292140
- **[IF=2.961]** Bao Qiwen. et al. Electrochemical biosensor based on antibody-modified Au nanoparticles for rapid and sensitive analysis of influenza A virus. IONICS. 2023 Mar;;1-9 Other ;. 10.1007/s11581-023-04944-w