

bs-0365R-HRP

[Secondary Antibodies]

human IgG Fc, HRP conjugated



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit

Isotype: IgG

Clonality: Polyclonal

Target: human IgG Fc

Concentration: 1mg/ml

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Background: The Fc region is the constant region on an immunoglobulin molecule, the area that is exactly the same on all antibodies. The Fc region is found on the heavy chains and is not involved in binding antigens.

Applications: **WB** (1:1000-3000)

IHC-P (1:100-500)

IHC-F (1:500-2000)

ELISA (1:1000-6000)

Reactivity: Human

Subcellular Location: Secreted