

bs-4620R

[Primary Antibody]



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

H1N1+H5N1+H3N2 nucleocapsid protein Rabbit pAb

DATASHEET

Host: Rabbit

Isotype: IgG

Clonality: Polyclonal

Target: H1N1+H5N1+H3N2 nucleocapsid protein

Immunogen: KLH conjugated synthetic peptide derived from H1N1/H5N1/H3N2 nucleocapsid protein: 151-250/498.

Purification: affinity purified by Protein A

Concentration: 1mg/ml

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Applications: ELISA (1:5000-10000)

Reactivity: (predicted: H1N1, H5N1, H3N2)

Predicted MW.: 56 kDa