

bs-23966R**[Primary Antibody]**

Smad4 Rabbit pAb

BioSS
ANTIBODIES

www.bioss.com.cn

sales@bioss.com.cn


techsupport@bioss.com.cn

400-901-9800

DATASHEET

Host: Rabbit	Isotype: IgG	Applications: WB (1:500-2000)
Clonality: Polyclonal		Reactivity: Human, Mouse (predicted: Rat, Rabbit, Pig, Sheep, Cow, Zebrafish, Chicken, Dog, Horse)
GeneID: 4089	SWISS: Q13485	Predicted MW.: 60 kDa
Target: Smad4		Observed MW.: 62 kDa
Immunogen: KLH conjugated synthetic peptide derived from human Smad4 : 451-552/552.		Subcellular Location: Cytoplasm ,Nucleus
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		
Background: Smad 4 is a member of the Mothers Against Dpp (MAD)-related family of proteins. So far, eight Smads have been identified and can be divided in 3 subgroups based on their structure and functions; pathway-restricted, common mediator and inhibitory Smad. Smad 4 is the common Smad (co-Smad). Previously identified as the tumor suppressor DPC4 (deleted in pancreatic carcinoma, locus 4), Smad 4 is functionally distinct among the Smad family, and is required for the assembly and transcriptional activation of diverse, Smad-DNA complexes. In contrast to the R-Smads, Smad 4 is not regulated by phosphorylation, but acts as a common mediator of TGF-Beta, activin, and bone morphogenetic protein signaling responses. Smad 4 is frequently inactivated in pancreatic, biliary and colorectal tumors.		

VALIDATION IMAGES



25 ug total protein per lane of various lysates (see on figure) probed with Smad4 polyclonal antibody, unconjugated (bs-23966R) at 1:1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r.t. for 60 min.

SELECTED CITATIONS

- **[IF=8.5]** Genghua Chen. et al. Bulk and single-cell alternative splicing analyses reveal roles of TRA2B in myogenic differentiation. CELL PROLIFERAT. 2023 Sep;;e13545 WB ;Chicken. 37705195
- **[IF=8.3]** Cai Bolin. et al. MYH1G-AS is a chromatin-associated lncRNA that regulates skeletal muscle development in chicken. CELL MOL BIOL LETT. 2024 Dec;29(1):1-25 WB ;Chicken. 38177995
- **[IF=5.959]** Wang Y et al. SPARCL1 promotes C2C12 cell differentiation via BMP7-mediated BMP/TGF-β cell signaling pathway. Cell Death Dis. 2019 Nov 7;10(11):852. WB ;Mouse. 31699966
- **[IF=4.175]** Huajun Wang. et al. LncRNA NEAT1 promotes proliferation, migration, invasion and epithelial-mesenchymal

Important Note: This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

transition process in TGF- β 2-stimulated lens epithelial cells through regulating the miR-486-5p/SMAD4 axis. *Cancer Cell Int.* 2020 Dec;20(1):1-12 WB ;Human. 33292220

- **[IF=2.795]** Chunyu Zhang et al. WISP1 promotes bovine MDSC differentiation via recruitment of ANXA1 for the regulation of the TGF- β signalling pathway. *Mol Cell Biochem.* 2020 Jul;470(1-2):215-227. WB ;Bovine. 32458119