
Recombinant Toxoplasma gondii SAG2

Catalog Number: bs-41560P

Concentration: 0.56mg/ml

Species: Toxoplasma gondii

Predicted MW: 19

Tags: C-His

Purity: >95% as determined by SDS-PAGE

Form: Lyophilized or Liquid

Storage: Phosphate buffered saline.

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Background: Toxoplasma gondii is an obligate intracellular protozoan parasite, it has been known associating with severe congenital defects and life threatening infection in immunocompromised hosts. In human, T. gondii infection in pregnant women can result in congenital infection, such as late-onset eye disease, it is also a frequent cause of encephalitis in AIDS patients: in animal, the infection can cause a variety diseases and symptoms.