

bs-41401R

[Primary Antibody]

PCV2 capsid protein Rabbit pAb

BioSS
ANTIBODIES

www.bioss.com.cn

sales@bioss.com.cn


techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: WB (1:500-2000) ELISA (1:5000-10000)
Clonality: Polyclonal		Reactivity: PCV2
Target: PCV2 capsid protein		
Immunogen: Recombined PCV2 Cap protein: full length protein .		
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		Predicted MW.: 28 kDa
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		Observed MW.: 41.5 kDa

— VALIDATION IMAGES —



Sample: Lane 1: Recombinant PCV2 Cap protein,
Trx & His Primary: Anti-PCV2 capsid protein
(bs-41401R) at 1/1000 dilution Secondary:
IRDye800CW Goat Anti-Rabbit IgG at 1/20000
dilution Predicted band size: 28 kDa Observed
band size: 41.5 kDa