

DATASHEET

Goat anti Human IgE/PE

Catalog Number: bs-2226G-PE

Target Protein: Goat anti Human IgE/PE

Concentration: 1.0 mg/ml

Host: Goat

Clonality: Polyclonal

Cross Reactive (predicted: Human)

Species:

Isotype: IgG

Swiss Prot: P0DOX4

Purification: affinity purified by Protein A

Storage: 10 mM TBS (pH=7.4) with 1% BSA, 0.03% Proclin300 and 50% glycerol. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: