

bs-0293G-BF594

[Secondary Antibodies]

Goat Anti-Rat IgG H&L, BF594 conjugated

BioSS
ANTIBODIES

www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Goat	Isotype: IgG	Applications: IF (1:200-1000)
Clonality: Polyclonal		Flow-Cyt (1:50-200)
Target: Goat Anti-Rat IgG H&L		ICC/IF (1:100-1000)
Purification: affinity purified by Protein G		Excitation Spectrum: 594nm
Concentration: 1.0 mg/ml		Emission spectrum: 621nm
Storage: 10 mM TBS (pH=7.4) with 1% BSA, 0.02% Proclin300 and 50% glycerol. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		Reactivity: Rat

— SELECTED CITATIONS —

- **[IF=14.7]** Yixuan Liu. et al. A protease-cleavable liposome for co-delivery of anti-PD-L1 and doxorubicin for colon cancer therapy in mice. Nature Communications. 2025 Mar 24;16(1):2854. Fluorescence resonance energy transfer assay ;Mouse. 40128211
- **[IF=8.3]** Gaiyue Yue. et al. Salvia miltiorrhiza Bunge Promotes Osteogenesis and Angiogenesis by Activating HIF-1 α /VEGF and Slit3 Pathways. PHYTOMEDICINE. 2026 Jan;;157842 IF ;Mouse. 41621348