
human IgG-Fc/AP

Catalog Number: bsm-0365M-AP

Target Protein: human IgG-Fc/AP

Concentration: 1mg/ml

Host: Mouse

Clonality: Monoclonal

Applications: WB

Reactivity: Human (predicted:)

Predicted MW: 13.8 kDa

Isotype: IgG1

Entrez Gene: N/A

Swiss Prot: N/A

Purification: affinity purified by Protein A

Storage: 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Background: The Fc region is the constant region on an immunoglobulin molecule, the area that is exactly the same on all antibodies. The Fc region is found on the heavy chains and is not involved in binding antigens.