

bsm-34316R**[Primary Antibody]**

Bioss
ANTIBODIES

www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

AIF Recombinant Rabbit mAb**— DATASHEET —****Host:** Rabbit**Clonality:** Recombinant**GeneID:** 9131**Target:** AIF**Immunogen:** A synthesized peptide derived from human AIF: 510-565.**Purification:** affinity purified by Protein A**Concentration:** 1mg/ml

Storage: pH7.4, 150mM NaCl, 0.02% Proclin300 and 50% Glycerol.
Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.


Background: This gene encodes a flavoprotein essential for nuclear disassembly in apoptotic cells, and it is found in the mitochondrial intermembrane space in healthy cells. Induction of apoptosis results in the translocation of this protein to the nucleus where it affects chromosome condensation and fragmentation. In addition, this gene product induces mitochondria to release the apoptogenic proteins cytochrome c and caspase-9. Mutations in this gene cause combined oxidative phosphorylation deficiency 6, which results in a severe mitochondrial encephalomyopathy. Alternative splicing results in multiple transcript variants. A related pseudogene has been identified on chromosome 10. [provided by RefSeq, May 2010].

Isotype: IgG**CloneNo.:** 20H17**SWISS:** O95831**Applications:** WB (1:1000-5000)**IHC-P** (1:200-800)**IHC-F** (1:200-800)**IF** (1:200-800)**Flow-Cyt** (1:50-100)**ICC/IF** (1:50-200)**Reactivity:** Human, Mouse, Rat


Predicted
MW.: 67 kDa

Observed
MW.: 67 kDa


Subcellular
Location: Cytoplasm ,Nucleus

— VALIDATION IMAGES —


25 ug total protein per lane of various lysates (see on figure) probed with AIF monoclonal antibody, unconjugated (bsm-34316R) at 1:1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r.t. for 60 min.




Paraformaldehyde-fixed, paraffin embedded Human Endometrial Cancer; Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15 min; The section was incubated with AIF Monoclonal Antibody, Unconjugated (bsm-34316R) at 1:200 overnight at 4°C, followed by conjugation to the bs-0295G-HRP and DAB (C-0010) staining.




Paraformaldehyde-fixed, paraffin embedded Human Kidney; Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15 min; Antibody incubation with AIF Monoclonal Antibody, Unconjugated (bsm-34316R) at 1:200 overnight at 4°C, followed by conjugation to the SP Kit(Rabbit, SP-0023) and DAB (C-0010) staining.



4% Paraformaldehyde-fixed HeLa (H) cell; Triton X-100 at r.t. for 20 min; Antibody incubation with (AIF) monoclonal Antibody, unconjugated (bsm-34316R) 1:100, 90 min at 37°C; followed by BF488 conjugated Goat Anti-Rabbit IgG antibody (green, bs-60295G-BF488) at 37°C for 90 min, DAPI (blue, C02-04002) was used to stain the cell



Paraformaldehyde-fixed, paraffin embedded Human Colon cancer; Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15 min; Antibody incubation with AIF Monoclonal Antibody, Unconjugated (bsm-34316R) at 1:200 overnight at 4°C, followed by a BF488 conjugated Goat Anti-Rabbit IgG (bs-0295G-




Paraformaldehyde-fixed, paraffin embedded Human Kidney; Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15 min; Antibody incubation with AIF Monoclonal Antibody, Unconjugated (bsm-34316R) at 1:200 overnight at 4°C, followed by a BF488 conjugated Goat Anti-Rabbit IgG (bs-0295G-

Important Note: This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

nuclei. PBS instead of the primary antibody was used as the blank control.

BF488) antibody at 37°C for 90 minutes, DAPI (blue, C02-04002) was used to stain the cell nuclei.

BF488) antibody at 37°C for 90 minutes, DAPI (blue, C02-04002) was used to stain the cell nuclei.



The HeLa (H) cells were fixed with 4% PFA (10 min at r.t.) and then permeabilized with 90% ice-cold methanol for 20 min at -20°C, the cells then were incubated in 5%BSA to block non-specific protein-protein interactions (30 min at r.t.), followed by secondary antibody incubation for 40 min at room temperature. Primary Antibody (green):Rabbit Anti-AIF antibody (bsm-34316R,1:100); Isotype Control (orange): Rabbit IgG (bs-0295P). Blank control (black): PBS. Acquisition of 20,000 events was performed.