


bsm-43589R**[Primary Antibody]****CD86 Recombinant Rabbit mAb****Bioss**
ANTIBODIES

www.bioss.com.cn


sales@bioss.com.cn

techsupport@bioss.com.cn


400-901-9800

— DATASHEET —**Host:** Rabbit**Isotype:** IgG**Clonality:** Recombinant**CloneNo.:** 8F6**GeneID:** 942**SWISS:** P42081**Target:** CD86**Immunogen:** Recombinant human CD86 protein: 24-239/329. < Extracellular >**Purification:** affinity purified by Protein A**Concentration:** 1mg/ml**Storage:** Size : 25ul/50ul/100ul/200ul
0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Size : 200ug (PBS only)
0.01M PBS
Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.**Background:** This gene encodes a type I membrane protein that is a member of the immunoglobulin superfamily. This protein is expressed by antigen-presenting cells, and it is the ligand for two proteins at the cell surface of T cells, CD28 antigen and cytotoxic T-lymphocyte-associated protein 4. Binding of this protein with CD28 antigen is a costimulatory signal for activation of the T-cell. Binding of this protein with cytotoxic T-lymphocyte-associated protein 4 negatively regulates T-cell activation and diminishes the immune response. Alternative splicing results in several transcript variants encoding different isoforms.[provided by RefSeq, May 2011].**Applications:** **Flow-Cyt** (1:50-100)
ICC/IF (1:50-200)**Reactivity:** Human**Predicted**
MW.: 37 kDa**Subcellular**
Location: Cell membrane**— VALIDATION IMAGES —**


4% Paraformaldehyde-fixed Daudi (H) cell;
Antibody incubation with (CD86) monoclonal Antibody, unconjugated (bsm-43589R) 1:100, 90 min at 37°C; followed by conjugated Goat Anti-Rabbit IgG antibody (green, bs-60295G-BF488) at 37°C for 90 min, DAPI (blue, C02-04002) was used to stain the cell nuclei. PBS instead of the primary antibody was used as the blank control.



4% Paraformaldehyde-fixed MJ (H) cell;
Antibody incubation with (CD86) monoclonal Antibody, unconjugated (bsm-43589R) 1:100, 90 min at 37°C; followed by BF488 conjugated Goat Anti-Rabbit IgG antibody (green, bs-60295G-BF488) at 37°C for 90 min, DAPI (blue, C02-04002) was used to stain the cell nuclei. PBS instead of the primary antibody was used as the blank control.



The MJ (H) cells were incubated in 5%BSA to block non-specific protein-protein interactions (30 min at r.t.), followed by secondary antibody incubation for 40 min at room temperature. Primary Antibody (green): Rabbit Anti-CD86 antibody (bsm-43589R,1:100); Isotype Control (orange): Rabbit IgG (bs-0295P). Blank control (black): PBS. Acquisition of 20,000 events was performed.



The Daudi (H) cells were incubated in 5%BSA to block non-specific protein-protein interactions (30 min at r.t.), followed by secondary antibody incubation for 40 min at room temperature.

Important Note: This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

Primary Antibody (green): Rabbit Anti-CD86 antibody (bsm-43589R,1:100); Isotype Control (orange): Rabbit IgG (bs-0295P). Blank control (black): PBS. Acquisition of 20,000 events was performed.