

bs-0886R

[Primary Antibody]

MeV F Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: ELISA (1:5000-10000)
Clonality: Polyclonal		Reactivity: (predicted: Measles virus)
Target: MeV F		
Immunogen: KLH conjugated synthetic peptide derived from Measles virus fusion protein: 451-550/550.		Predicted MW.: 57 kDa
Purification: affinity purified by Protein A		Subcellular Location: Nucleus
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		
Background: No data available.		

— SELECTED CITATIONS —

- **[IF=3.616]** Katelyn D. Miller. et al. Murine BST2/tetherin promotes measles virus infection of neurons. Virology. 2021 Nov;563:38 IF ;Mouse. 34416448