

bs-10381R**[Primary Antibody]****phospho-Smad1 (Ser463) Rabbit pAb****Bioss**
ANTIBODIES

www.bioss.com.cn


sales@bioss.com.cn

techsupport@bioss.com.cn


400-901-9800

— DATASHEET —


Host: Rabbit	Isotype: IgG	Applications: WB (1:500-2000)
Clonality: Polyclonal		IHC-P (1:100-500)
GeneID: 4086	SWISS: Q15797	IHC-F (1:100-500)
Target: Smad1 (Ser463)		IF (1:100-500)
		Flow-Cyt (2ug/Test)
Immunogen: KLH conjugated Synthesised phosphopeptide derived from human Smad1 around the phosphorylation site of Ser463: IS(p-S)VS.		Reactivity: Human, Mouse, Rat (predicted: Pig)
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		Predicted MW.: 51 kDa
		Observed MW.: 51 kDa
Background: The protein encoded by this gene belongs to the SMAD, a family of proteins similar to the gene products of the Drosophila gene 'mothers against decapentaplegic' (Mad) and the C. elegans gene Sma. SMAD proteins are signal transducers and transcriptional modulators that mediate multiple signaling pathways. This protein mediates the signals of the bone morphogenetic proteins (BMPs), which are involved in a range of biological activities including cell growth, apoptosis, morphogenesis, development and immune responses. In response to BMP ligands, this protein can be phosphorylated and activated by the BMP receptor kinase. The phosphorylated form of this protein forms a complex with SMAD4, which is important for its function in the transcription regulation. This protein is a target for SMAD-specific E3 ubiquitin ligases, such as SMURF1 and SMURF2, and undergoes ubiquitination and proteasome-mediated degradation. Alternatively spliced transcript variants encoding the same protein have been observed. [provided by RefSeq].		Subcellular Location: Nucleus

— VALIDATION IMAGES —


Sample: Heart (Mouse) Lysate at 40 ug Primary:
Anti- Phospho-Smad1 (Ser463) (bs-10381R) at
1/300 dilution Secondary: IRDye800CW Goat
Anti-Rabbit IgG at 1/20000 dilution Predicted
band size: 51 kD Observed band size: 51 kD




Sample: Muscle (Mouse) Lysate at 40 ug Primary:
Anti-Phospho-Smad1 (Ser463) (bs-10381R) at
1/500 dilution Secondary: IRDye800CW Goat
Anti-Rabbit IgG at 1/20000 dilution Predicted
band size: 51 kD Observed band size: 51 kD



Paraformaldehyde-fixed, paraffin embedded
(Rat heart); Antigen retrieval by boiling in
sodium citrate buffer (pH6.0) for 15min; Block
endogenous peroxidase by 3% hydrogen
peroxide for 20 minutes; Blocking buffer (normal
goat serum) at 37°C for 30min; Antibody
incubation with (Phospho-Smad1 (Ser463))
Polyclonal Antibody, Unconjugated (bs-10381R)
at 1:400 overnight at 4°C, followed by operating
according to SP Kit(Rabbit) (sp-0023)
instructions and DAB staining.



Paraformaldehyde-fixed, paraffin embedded (human brain glioma); Antigen retrieval by microwave in sodium citrate buffer (pH6.0) ; Block endogenous peroxidase by 3% hydrogen peroxide for 30 minutes; Blocking buffer (3%BSA) at RT for 30min; Antibody incubation with (Smad1 (Ser463)) Polyclonal/Monoclonal Antibody, Unconjugated (bs-10381R) at 1:400 overnight at 4°C, followed by conjugation to the secondary antibody (labeled with HRP) and DAB staining.



Blank control (black line) : HeLa. Primary Antibody (green line): Rabbit Anti-Phospho-Smad1 (Ser463) antibody (bs-10381R) Dilution: 2ug/Test; Secondary Antibody (white/blue line) : Goat anti-rabbit IgG-FITC Dilution: 0.5ug/Test. Isotype control (orange line) : Normal Rabbit IgG Protocol The cells were fixed with 4% PFA (10min at room temperature) and then permeabilized with 90% ice-cold methanol for 20 min at -20°C. The cells were then incubated in 5% BSA to block non-specific protein-protein interactions for 30 min at room temperature. Cells stained with Primary Antibody for 30 min at room temperature. The secondary antibody used for 40 min at room temperature. Acquisition of 20,000 events was performed.