

bsm-56001R**[Primary Antibody]****Bioss**
ANTIBODIES

www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

C-Myc Recombinant Rabbit mAb**— DATASHEET —****Host:** Rabbit**Isotype:** IgG1**Clonality:** Recombinant**CloneNo.:** 18H19**Target:** C-Myc**Immunogen:** A synthesized peptide derived from human c Myc: 1-50.**Purification:** affinity purified by Protein A**Concentration:** 1mg/ml**Storage:** 0.1M Phosphate Buffered Saline, pH7.4.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.


Background: The protein encoded by this gene is a multifunctional, nuclear phosphoprotein that plays a role in cell cycle progression, apoptosis and cellular transformation. It functions as a transcription factor that regulates transcription of specific target genes. Mutations, overexpression, rearrangement and translocation of this gene have been associated with a variety of hematopoietic tumors, leukemias and lymphomas, including Burkitt lymphoma. There is evidence to show that alternative translation initiations from an upstream, in-frame non-AUG (CUG) and a downstream AUG start site result in the production of two isoforms with distinct N-termini. The synthesis of non-AUG initiated protein is suppressed in Burkitt's lymphomas, suggesting its importance in the normal function of this gene. [provided by RefSeq, Jul 2008].

Applications: **WB** (1:2000-10000)
Flow-Cyt (1µg/Test)
ICC/IF (1:50-200)

Reactivity: Human

Predicted
MW.: 51 kDa


Subcellular
Location: Nucleus

— VALIDATION IMAGES —

25 µg total protein per lane of various lysates (see on figure) probed with c-Myc monoclonal antibody, unconjugated (bsm-56001R) at 1:2000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r.t. for 60 min.



4% Paraformaldehyde-fixed HeLa (H) cell; Triton X-100 at r.t. for 20 min; Antibody incubation with (c-Myc) monoclonal Antibody, unconjugated (bsm-56001R) 1:100, 90 min at 37°C; followed by conjugated Goat Anti-Rabbit IgG antibody (green, bs-60295G-BF488) at 37°C for 90 min, DAPI (blue, C02-04002) was used to stain the cell nuclei. PBS instead of the primary antibody was used as the blank control.



The HeLa (H) cells were fixed with 4% PFA (10 min at r.t.) and then permeabilized with 90% ice-cold methanol for 20 min at -20°C, the cells then were incubated in 5% BSA to block non-specific protein-protein interactions (30 min at r.t.), followed by secondary antibody incubation for 40 min at room temperature. Primary Antibody (green): Rabbit Anti-c-Myc antibody (bsm-56001R): 1 µg/10⁶ cells; Isotype Control (orange): Rabbit IgG (bs-0295P). Blank control (black): PBS. Acquisition of 20,000 events was performed.