

bsm-41573M

[Primary Antibody]

Rotavirus Mouse mAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Mouse

Isotype: IgG1

Clonality: Monoclonal

SWISS: P04509

Target: Rotavirus

Purification: affinity purified by Protein A

Storage: PBS pH 7.0.
Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Applications: ELISA (1:5000-20000)

Reactivity: (predicted: Rotavirus)

Predicted MW.: 45 kDa