

bs-10299R

[Primary Antibody]

H7N9 Matrix Protein 1 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: WB (1:500-2000)
Clonality: Polyclonal		IHC-P (1:100-500)
Target: H7N9 Matrix Protein 1		IHC-F (1:100-500)
Immunogen: KLH conjugated synthetic peptide derived from H7N9 Matrix Protein 1: 161-252/252.		IF (1:100-500)
Purification: affinity purified by Protein A		ELISA (1:5000-10000)
Concentration: 1mg/ml		Reactivity: (predicted: Human)
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		Predicted MW.: 28 kDa
		Subcellular Location: Cytoplasm ,Nucleus

— SELECTED CITATIONS —

- **[IF=3.293]** Jiao Hu. et al. Baculovirus-derived influenza virus-like particle confers complete protection against lethal H7N9 avian influenza virus challenge in chickens and mice. Vet Microbiol. 2022 Jan;264:109306 IF ;Spodoptera frugiperda. 34923247