

**bsm-33346R****[ Primary Antibody ]****BioSS**  
ANTIBODIES

www.bioss.com.cn


sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

**Flag Tag Recombinant Rabbit mAb****— DATASHEET —**

|                                                                                                                                                                      |                       |                                                                                                                                                                                                                                                                                                            |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Host:</b> Rabbit                                                                                                                                                  | <b>Isotype:</b> IgG   | <b>Applications:</b> <b>WB</b> (1:5000-50000)<br><b>IHC-P</b> (1:200-1000)<br><b>IHC-F</b> (1:200-1000)<br><b>IF</b> (1:200-1000)<br><b>Flow-Cyt</b> (1 $\mu$ g/Test)<br><b>ICC/IF</b> (1:100-500)<br><b>IP</b> (1:20-50)<br><br><b>Reactivity:</b> Species independent<br><br><b>Observed MW.:</b> 48 kDa |
| <b>Clonality:</b> Recombinant                                                                                                                                        | <b>CloneNo.:</b> 4E16 |                                                                                                                                                                                                                                                                                                            |
| <b>Target:</b> Flag Tag                                                                                                                                              |                       |                                                                                                                                                                                                                                                                                                            |
| <b>Purification:</b> affinity purified by Protein A                                                                                                                  |                       |                                                                                                                                                                                                                                                                                                            |
| <b>Concentration:</b> 1mg/ml                                                                                                                                         |                       |                                                                                                                                                                                                                                                                                                            |
| <b>Storage:</b> 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.<br>Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |                       |                                                                                                                                                                                                                                                                                                            |
| <b>Background:</b> This is a useful tool for the localisation and characterisation of DDDDK tagged proteins.                                                         |                       |                                                                                                                                                                                                                                                                                                            |

**— VALIDATION IMAGES —**

5 ng rhTNFSF13B-flag protein (bs-43132P) per lane probed with Flag Tag monoclonal antibody respectively, unconjugated (bsm-33346R) at 1:30000 dilution and 4°C overnight incubation. Followed by corresponding conjugated secondary antibody incubation at r.t. for 60 min.

**— SELECTED CITATIONS —**

- **[IF=2.4]** Jingwen Ma. et al. Enhancement of protective efficacy of recombinant attenuated Salmonella typhimurium delivering H9N2 avian influenza virus hemagglutinins(HA) antigen vaccine candidate strains by C-C motif chemokine ligand 5 in chickens(chCCL5). VET MICROBIOL. 2024 Nov;298:110264 WB ;Human. 39395372