

bsm-33004R**[Primary Antibody]****Bioss**
ANTIBODIES

www.bioss.com.cn


sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

His tag Recombinant Rabbit mAb**— DATASHEET —**

Host: Rabbit	Isotype: IgG	Applications: WB (1:1000-10000) IHC-P (1:200-1000) IHC-F (1:200-1000) IF (1:200-1000) Flow-Cyt (1 μ g/Test) ICC/IF (1:100-500) IP (1:20-50) Reactivity: Species independent Observed MW.: 30 kDa
Clonality: Recombinant	CloneNo.: 16A3	
Target: His tag		
Purification: affinity purified by Protein G		
Storage: PBS (pH=7.4) Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		
Background: The H-H-H-H-H motif is used as a tag on many recombinant proteins to facilitate purification. The antibody recognizes the His-tag fused to the amino- or carboxy- termini of targeted proteins in transfected or transformed cells.		

— VALIDATION IMAGES —

20 ng rGFP-His protein (bs-33009P) per lane probed with His tag monoclonal antibody respectively, unconjugated (bsm-33004R) at 1:1000 dilution and 4°C overnight incubation. Followed by corresponding conjugated secondary antibody incubation at r.t. for 60 min.