

[Primary Antibody]

ATG5 Recombinant Rabbit mAb



BioSS
ANTIBODIES

www.bioss.com.cn

sales@bioass.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET

Host: Rabbit

Clonality: Recombinant

Isotype: IgG

CloneNo.: 4G5

Target: ATG5

Immunogen: A synthesized peptide derived from human ATG5: 1-64.

Purification: affinity purified by Protein A

Concentration: 1mg/ml

Storage: 0.01M TBS(pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Background: In yeast, autophagy is an essential process for survival during nutrient starvation and cell differentiation. The process of autophagy is characterized as a non-selective degradation of cytoplasmic proteins into membrane structures called autophagosomes, and it is dependent on several proteins, including the autophagy proteins APG5 and APG7. Yeast Apg7 and the human homolog, APG7, share similarities with the ubiquitin-activating enzyme E1 in *Saccharomyces cerevisiae* and are likewise responsible for enzymatically activating the autophagy conjugation system. Apg5 and the human homolog, APG5 (also designated apoptosis-specific protein or APS), function as substrates for the autophagy protein Apg12. These proteins are covalently bonded together to form Apg12/APG5 conjugates, which are required for the progression of autophagy.

Applications: WB (1:500-2000)

IHC-P (1:50-200)

IHC-F (1:50-200)

IF (1:50-200)


IP (1:20-100)

Reactivity: Human, Mouse, Rat

Predicted
MW.: 32 kDa

Subcellular Location: Cytoplasm


— VALIDATION IMAGES



25 ug total protein per lane of various lysates (see on figure) probed with ATG5 monoclonal antibody, unconjugated (bsm-52596R) at 1:1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r.t. for 60 min.



Paraformaldehyde-fixed, paraffin embedded Human Breast Cancer; Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15 min; Antibody incubation with ATG5 Monoclonal Antibody, Unconjugated(bsm-52596R) at 1:200 overnight at 4°C, followed by conjugation to the bs-0295G-HRP and DAB (C-0010) staining.



Paraformaldehyde-fixed, paraffin embedded
Human placenta; Antigen retrieval by boiling in
sodium citrate buffer (pH6.0) for 15 min;
Antibody incubation with ATG5 Monoclonal
Antibody, Unconjugated (bsm-52596R) at 1:200
overnight at 4°C, followed by conjugation to the
bs-0295G-HRP and DAB (C-0010) staining.

— SELECTED CITATIONS

- **[IF=1.8]** Huang Pan. et al. Imperatorin promotes melanin degradation in keratinocytes through facilitating autophagy via the PI3K/Akt signaling pathway. ARCH DERMATOL RES. 2025 Dec;317(1):1-11 **WB** ;Human. 39636461

Important Note: This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.