

bsm-61071R**[Primary Antibody]**

BioSS
ANTIBODIES

www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn


400-901-9800

Caspase-3 Recombinant Rabbit mAb


— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: WB (1:1000-5000) IHC-P (1:200-800) IHC-F (1:200-800) IF (1:200-800)
Clonality: Recombinant	CloneNo.: 4D2	
GeneID: 836	SWISS: P42574	
Target: Caspase-3		
Immunogen: A synthesized peptide derived from human Caspase 3: 176-277/277.		
Purification: affinity purified by Protein A		Reactivity: Human, Mouse, Rat
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		Predicted MW.: 32 kDa Observed MW.: 34 kDa
Background: The caspase family of cysteine proteases play a key role in apoptosis. Caspase 3 is the most extensively studied apoptotic protein among caspase family members. Caspase 3 is synthesized as inactive pro enzyme that is processed in cells undergoing apoptosis by self proteolysis and/or cleavage by other upstream proteases (e.g. Caspases 8, 9 and 10). The processed form of Caspase 3 consists of large (17kDa) and small (12kDa) subunits which associate to form an active enzyme. Caspase 3 is cleaved at Asp28 Ser29 and Asp175 Ser176. The active Caspase 3 proteolytically cleaves and activates other caspases (e.g. Caspases 6, 7 and 9), as well as relevant targets in the cells (e.g. PARP and DFF). Alternative splicing of this gene results in two transcript variants which encode the same protein. In immunohistochemical studies Caspase 3 expression has been shown to be widespread but not present in all cell types (e.g. commonly reported in epithelial cells of skin, renal proximal tubules and collecting ducts). Differences in the level of Caspase 3 have been reported in cells of short lived nature (eg germinal centre B cells) and those that are long lived (eg mantle zone B cells). Caspase 3 is the predominant caspase involved in the cleavage of amyloid beta 4A precursor protein, which is associated with neuronal death in Alzheimer's disease.		Subcellular Location: Cytoplasm


— VALIDATION IMAGES —




25 ug total protein per lane of various lysates (see on figure) probed with Caspase-3 monoclonal antibody, unconjugated (bsm-61071R) at 1:1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r.t. for 60 min.




25 ug total protein per lane of various lysates (see on figure) probed with Caspase-3 monoclonal antibody, unconjugated (bsm-61071R) at 1:2000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r.t. for 60 min.



Paraformaldehyde-fixed, paraffin embedded Human Breast Cancer; Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15 min; The section was incubated with Caspase-3 Monoclonal Antibody, Unconjugated (bsm-61071R) at 1:200 overnight at 4°C, followed by conjugation to the bs-0295G-HRP and DAB (C-0010) staining.



Paraformaldehyde-fixed, paraffin embedded Human Tonsil; Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15 min; Antibody incubation with Caspase-3 Monoclonal Antibody, Unconjugated (bsm-61071R) at 1:200 overnight at 4°C, followed by conjugation to the SP Kit (Rabbit, SP-0023) and DAB (C-0010) staining.



Paraformaldehyde-fixed, paraffin embedded Human Tonsil; Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15 min; The section was incubated with Caspase-3 Antibody, Unconjugated (bsm-61071R) at 1:200 overnight at 4°C. Followed by conjugated Goat Anti-Rabbit IgG antibody (Green, bs-0295G-BF488), DAPI (blue, C02-04002) was used to stain the cell nuclei.

— SELECTED CITATIONS —

- **[IF=19]** Zhongming Huang. et al. Dual-Functional Sonodynamic Therapeutic Platform for Complete Osteosarcoma Elimination With Synchronized Bone Defect Regeneration. ADV FUNCT MATER. 2026 Feb;;e31043 IF ;Human. 10.1002/adfm.202531043
- **[IF=17.1]** Lei Liu. et al. Myricetin Oligomer Triggers Multi-Receptor Mediated Penetration and Autophagic Restoration of Blood-Brain Barrier for Ischemic Stroke Treatment. ACS NANO. 2024;XXXX(XXX):XXX-XXX WB,IF ;Mouse,Rat. 38533773
- **[IF=9.7]** Jie Meng. et al. pH-triggered Schottky heterojunctions for NIR-II-activated and tumor-specific pyroelectrodynamics and photothermal therapy. J COLLOID INTERF SCI. 2026 Jan;702:138919 IHC ;Mouse. 40915294
- **[IF=8.1]** Liu Anhong. et al. Zinc ions trigger PANoptosis-like cell death in magnetic hyperthermia therapy of magnesium based implant for hepatocellular carcinoma. APOPTOSIS. 2025 Jun;;1-13 WB ;Mouse. 40537603
- **[IF=4.9]** Yuting Wang. et al. PPAR-γ Inhibits Chronic Apical Periodontitis by Facilitating Macrophage Efferocytosis. INT J MOL SCI. 2025 Jan;26(20):10157 IF ;Rat. 41155449