

CD97 Recombinant Rabbit mAb

Catalog Number: bsm-61103R

Target Protein: CD97

Form: Liquid

Host: Rabbit

Clonality: Recombinant

Isotype: IgG

Applications: WB (1:1000-1:2000), IHC-P (1:100-1:200), IHC-F (1:100-1:200), IF (1:50-1:200), Flow-Cyt (1:20-1:100), ICC/IF (1:50-1:200), IP (1:20-1:50)

Reactivity: Human, Mouse

Predicted MW: 92

Detected MW: 92

Entrez Gene: 976

Swiss Prot: P48960

Source: A synthesized peptide derived from human ADGRE5: 450-500.


Purification: affinity purified by Protein A

Storage: 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.

Background: Receptor potentially involved in both adhesion and signaling processes early after leukocyte activation. Plays an essential role in leukocyte migration.

VALIDATION IMAGES



Western blot analysis of U937 cell lysate. Using CD97 (bsm-61103R) monoclonal antibody at 1:1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r.t. for 60 min.