

YB1 Recombinant Rabbit mAb

Catalog Number: bsm-61192R

Target Protein: YB1

Form: Liquid

Host: Rabbit

Clonality: Recombinant

Clone No.: 15E6

Isotype: IgG

Applications: WB (1:500-1:2000), IHC-P (1:100-1:200), IHC-F (1:100-1:200), IF (1:50-1:200), Flow-Cyt (1:20-1:100), ICC/IF (1:50-1:200), IP (1:20-1:50)

Reactivity: Human, Mouse, Rat

Predicted MW: 36

Detected MW: 50

Entrez Gene: 4904

Swiss Prot: P67809


Purification: affinity purified by Protein A

Storage: 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.

Background: DNA- and RNA-binding protein involved in various processes, such as translational repression, RNA stabilization, mRNA splicing, DNA repair and transcription regulation.

VALIDATION IMAGES



25 ug total protein per lane of various lysates (see on figure) probed with YB1 monoclonal antibody, unconjugated (bsm-61192R) at 1:2000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r.t. for 60 min.